

**OLYMPUS**<sup>®</sup>

Your Vision, Our Future

System Microscope

**BX43/BX46/BX53**

**BX3 Series**

Comfort and Performance  
in Efficient Harmony



**UIS2**  
World-leading optics



Comfort &



**BX46**  
Clinical Microscope



**BX43**  
System Microscope



**BX53**  
System Microscope

# A Revolutionary New Standard in Operating Comfort and Imaging Efficiency

The Olympus BX43/BX46/BX53 series defines a revolutionary new standard in laboratory and clinical microscopy. Ergonomic design can provide enhanced operator comfort during long hours of use, with an intuitive control layout for fast, efficient observation and imaging.

Efficiency





## Comfortable Posture for the Operator

### Binocular Tube Lifts to Meet You Halfway

The U-TTLBI binocular tube lifts, extends and tilts with simple adjustments for optimum posture. This new design allows long hours of observation to be more comfortable, regardless of physique.



### Tilting Binocular tube Lineup Caters To Your Individual Needs


A diverse lineup of tilting tubes is now available to meet various demands, including cost performance with U-TBI-3 and erect image observation with U-ETBI. Eyepoint adjusters/U-EPA2 and U-EPAL-2 accommodate height differences, resulting in a height-appropriate eye-point fit.



## Diminishes Eye Fatigue and Repetitive Movements

### Light Intensity Manager Controls Brightness


The BX43/BX46 assigns objectives their own brightness, freeing the operator from light intensity adjustments with changes in magnification. Uniform brightness is maintained with objective changeovers, eliminating light intensity control and reducing eye fatigue. Brightness can also be set to suit operator preferences.




### Olympus LED\*, Optimal Lighting for Pathology and Cytology

The BX43/BX46 offers a solution for detecting purples, cyans, and reds that are often difficult to view with white LEDs. Olympus white LED with high color rendering index has wavelength characteristics close to those of halogen, realizing halogen lamp-equivalent color reproductions. LED lighting retains color temperature even preventing images from turning reddish or bluish in color.

\*LED was jointly developed by Olympus Corporation and CCS Inc.



Spectral Characteristics\*



\*This graph shows the spectral characteristics of each light source regularized with the luminosity curve. It does not compare the strength of light for each light source.



## Accessible Comfort in Conference and Imaging

### Tilting Trinocular Tube for Optimal Posture

A camera can be attached to the tilting trinocular tube/U-TTR-2 which accommodates operator physique and posture. The optical path switch can be attached to either side of the tube.



### Hand Switch for Image Capture

The U-HSEXP hand switch for exposure attaches easily on the microscope, and allows the operator to capture an image without having to take his/her eyes away from the specimen. The exposure switch is located close to the focus handle creating minimal movement during operation.




### DP21 Digital Camera Takes Comfort Further

The Olympus DP21, a stand-alone digital camera that is also computer-friendly, provides accurate color reproductions and smooth, high-definition live image displays that are ideal for conference use. The convenient handset ensures simple operation during image acquisition and measurement.



BX53+Digital Camera DP21 (Stand-alone) Configuration




Stomach (HE Stain)

### Capturing Images Made Easy

The Olympus digital camera, using the cellSens imaging software platform, facilitates image capture. cellSens is simple to operate, allowing the operator to customize his/her workflow.



BX53+Digital Camera DP73 Configuration



Lung (HE Stain)

### Accurate Post-imaging Measurements

Integration with coded revolving nosepiece allows sharing and recording of objective magnification. The coded revolving nosepiece eliminates errors that occur when the wrong magnification is manually recorded by the operator.




BX43



# BX43



## Comfort and Efficiency in Operation



Cervical Cytology (Papanicolaou Stain)

### Light Intensity Manager Controls Brightness


This feature functions with the BX43-5RES, a coded 5-position nosepiece for BX43, eliminating light intensity adjustments when switching magnifications. The long-life LED ensures uniform brightness with easy maintenance.




Hematology (Giemsa Stain)

### Low Magnification Condenser U-LC

Changing objective magnification from 2x to 100x (dry) is now possible without having to change the condenser or move the top lens.



Kidney (Fibrin, PTAH Stain)



### Various Units to Accommodate Observation Style


A diverse lineup of module units including ergonomic observation tubes and stages are available to suit individual applications.



BX46

# BX46


## Improved User Comfort and Efficient Operation



Stomach (HE Stain)

### Ergonomic Tube Provides Comfort in Posture

The binocular tube/U-TTLBI tilts, extends and lifts to accommodate operator physique and posture. This feature enhances comfort during long hours of observation.




Breast (Anti HER2)

### Revolving Nosepiece with Light Intensity Manager Functionality

The light intensity manager sets brightness settings according to the specific objective. Uniform brightness is achieved from low to high magnifications, eliminating light intensity adjustments, therefore reducing eye fatigue. The long-life LED ensures uniform brightness with easy maintenance.

### Light Action, Low Torque Stage — Reduced Operator Force to Move the Specimen

Quick yet thorough screenings with simple finger taps alleviate operator fatigue. Low-position handles require only minimal movements, with the lower arm resting comfortably on the desk.



Cervical Cytology (Papanicolaou Stain)

### Low, Fixed Stage for Ease in Changing Specimens


Working comfort increases with our low-position fixed stage and nosepiece focus to suit the individual user. Specimen is easily set with minimum hand operations and visibility is improved to reduce operator fatigue.




BX53

# BX53


## Advanced Modularity to Suit Observation Style



Liver (HBV, Victoria Blue Stain)



Stomach (HE Stain)



Breast (HER2, FISH)

### Ideal Scalability Answers Various Demands

In addition to integrations with modular components for brightfield observation, various other modular units such as fluorescence, phase contrast, polarization and darkfield are available to enhance scalability. The 100 W halogen-transmitted light source can meet various imaging and multi observation modes.

### Further Comfort in Fluorescence Observation

Olympus introduces special high performance filters for all fluorescence mirror units. The fluorescence illuminators can carry eight mirror units simultaneously to seamlessly image multi-stained specimens. Furthermore, no tools are required for mirror unit replacements.



### Customizable Control Layout

The BX53 is designed for maximum flexibility in placement, with a centrally positioned light intensity dial and a fluorescence illuminator shutter that can be operated with either hand. In addition, the fine focus handles can be attached on either the right or left side of the microscope according to operator preference.



### Energy-saving Switch Turns off Automatically

The motion sensor detects when an operator leaves and automatically turns off the transmitted light lamp after around 30 minutes. The energy-saving switch conserves energy and lamp lifetime.



# A Microscope Designed to Meet Your Specifications and Needs

The Olympus infinity-corrected optical system UIS2 satisfies future scalability needs. Inserting an optical element into the infinity space causes no additional image distortion or deterioration in image quality.

**UIS2**  
World-leading optics

## ● UPLSAPO Series

Thanks to the application of the Olympus original UW multi-coatings, these Super Apochromat objectives compensate for both spherical and chromatic aberrations from the UV to the near infrared region. Their sensitivity to fluorescence emissions ensures the acquisition of sharp, clear images, without color shift, even in brightfield observation. For quality and performance, they offer solutions for digital imaging needs.



## ● PLAPON Series

Designed for unsurpassed resolution and contrast, these Plan Apochromat objectives reduce chromatic aberration to an absolute minimum. 1.25x and 2x objectives are available.



## ● UPLFLN (UPLFLN-PH) Series

These plan objectives also provide flat images with high transmission up to the near infrared region of the spectrum. With their high signal to noise ratio, excellent resolution and high contrast imaging, they are especially effective in brightfield observation. The UPLFLN-PH series is optimized for phase contrast observation.



## ● PLN (PLN-PH) Series

Ideal for a range of clinical and research applications, these high quality objectives feature excellent flatness up to FN 22 in transmitted brightfield (phase contrast) observation. The PLN-PH series is specifically designed for phase contrast observation.



## ● No Cover Objectives

The Olympus coverglass-free objectives are designed for glass slides without a cover slip. This feature is largely designed for observation of blood smear specimens.




## Brightfield

### Brighter Images, with Superb Resolution/Flatness at All Magnifications

Olympus offers a diverse line of condensers including: the U-SC3 swing-out condenser, suitable for observations from 1.25x to 100x; the U-LC, for consecutive observations from 2x to 100x (dry); the U-AAC, whose Aplanat-Achromat design reduces chromatic aberration; and the U-ULC-2 special condenser for ultra low magnifications.

\*Select the U-ULC2 condenser for optimal digital imaging with the 1.25x objective.



Pancreas (HE)




①U-SC3 ②U-ULC-2 ③U-AC2 ④U-AAC  
⑤U-LC

## Fluorescence

### Olympus Takes Fluorescence Observation to Another Plane

Two types of reflected illuminators are available, the universal reflected illuminator/BX3-URA and the coded fluorescence illuminator/BX3-RFAS. A total of eight fluorescence mirror units can be attached for comfortable multi-color fluorescence observations. High-performance filters provide efficient, bright and even fluorescence observations.



Kidney (Anti IgG)




BX3-URA

BX3-RFAS

## Phase Contrast

### High-contrast, High-resolution Imaging

High contrast phase imaging allows close observation of the cell interior and of live bacteria. Using UPLFLN-PH or PLN-PH series objectives, phase contrast observation from 10x up to 100x is available. With the U-PCD2 phase/darkfield condenser, users can view specimens in brightfield or darkfield. Simultaneous observation with reflected light fluorescence microscopy is also possible.



Asbestos




U-PCD2

## Polarized Light

### High-resolution View of Double Refraction Structure in Cells

Tooth, bone, muscle tissue, nerve tissue, actomyosin fiber and mitotic spindle can all be observed, without staining. There are intermediate attachments (U-OPA/U-CPA) for orthoscopic and orthoscopic/conoscopic viewing. Various compensators make it possible to observe a wide range of retardation. Also available is a condenser exclusively for polarization observation, revolving nosepiece, rotating stage, objectives, simple polarizing attachment and analyzer to detect uric acid crystal.



Urate Crystals




①U-POC-2 ②U-CPA ③U-OPA ④U-AN360P-2  
⑤U-P4RE ⑥U-GAN ⑦BX45-PO

## Darkfield

### Excellent Darkfield Effect from Low to High Magnifications

Choose from the 10x to 100x dry darkfield condenser/U-DCD or the 20x to 100x oil immersion darkfield condenser/U-DCW.

\*Please consult your nearest Olympus representative for applicable objectives.



Spirogyra




U-DCD

U-DCW

# Group Observation Systems

Olympus discussion systems are invaluable for lab training and education. There is a wide choice, designed for 2 to 10 participants. The pointer is powered by LED, so there is no concern for sudden lamp failure.



- Lead-and-arsenic-free Eco-glass for optics, such as lenses and prisms
- Exclusion of hexavalent chrome from metal materials and surface metal treatments
- Exclusion of lead from electric components and solders
- Exclusion of chlorofluorocarbons in production
- Compliance with laws and regulations : RoHS (EU), Chinese RoHS (China), WEEE (EU)

## Transport

- Adoption of cardboard packaging

## Development/ Production

## Enhanced Performance Even for the Environment

Olympus Group will develop products, services and production technologies with a careful and conscientious regard for safety and environmental protection, so that the society and environment we pass on will be sound and full of health. BX43/BX46/BX53 are manufactured under Olympus' own green designing standards, established in reference to Type II environmental label indication stipulations in the international standard ISO 14021.

## Use


- LED illumination White LED as light source, BX43 and BX46 reduce power consumption by 10%
- ECO (energy-saving) modes The BX53 sensor turns off the transmitted light lamp after around 30 minutes.

## Waste Disposal

- Waste sorting
- Improved recycling rate




**BX46 system diagram**



\*1 Vignetting may occur in combination with an additional intermediate attachment.  
\*3 Attached to BX46F.

\*2 Only U-EPA-2 and U-EPAL-2 are able to use as an additional intermediate attachment.



\*1 Slight vignetting may occur in combination with an additional intermediate attachment or observation method. \*2 Require an additional intermediate attachment or fluorescence illuminator. \*3 Cannot be used with U-TTLBI. \*4 Compatible with FN 22. \*5 An auxiliary lens is equipped.

**CAMERAS**

- U-TV0.25XC**  
0.25X  
C-Mount  
Adapter
- U-TV0.35XC-2**  
0.35X  
C-Mount  
Adapter
- U-TV0.5XC-3**  
0.5X  
C-Mount  
Adapter
- U-TV0.63XC-3**  
0.63X  
C-Mount  
Adapter
- U-TV1XC**  
C-Mount  
adapter 1X  
(XY adjustment)

- U-ECA**  
Magnification changer 2x
- U-CA**  
Magnification changer
- U-DP\*1\*3**  
Dual port
- U-DP1XC**  
Dual port 1x
- U-APT**  
Arrow pointer

- U-TBI-3-CLI\*1**  
Tilting binocular tube
- U-TTBI**  
Ergonomic  
binocular tube
- U-ETBI**  
Ergonomic erect  
image binocular  
tube
- U-TTLBI\*2**  
Tilting,  
telescopic,  
lifting binocular tube

- U-TRU\*1\*3**  
Trinocular intermediate unit
- U-TRUS\*3**  
Trinocular intermediate unit

- U-5RE-2**  
Quintuple  
revolving  
nosepiece
- BX43-5RES**  
Coded  
5-position  
nosepiece  
for BX43

- U-LLGAD**  
Liquid light  
guide  
adapter
- U-LLG150/  
U-LLG300**  
Liquid light  
guide  
(1.5 m/3 m)
- U-HGLGPS**  
Light source

- U-ANT**  
Analyzer for  
transmitted light
- U-KPA**  
Intermediate attachment for  
simple polarizing observation
- U-EPA2**  
Eyepoint adjuster
- U-EPAL-2**  
Eyepoint adjuster

- BX3-URA**  
Universal reflected illuminator
- U-PO**  
Polarizer
- BX3-6ND6**  
ND filter
- BX3-25ND25**  
ND filter
- U-AN-2**  
Analyzer slider
- BX3-RFAS\*4**  
Coded  
fluorescence  
illuminator
- U-DULHA**  
Double  
lamp house  
adapter

- U-LH75XEAP0**  
75 W xenon  
apo lamp  
housing
- U-LH100HGAP0**  
100 W mercury  
apo lamp  
housing
- U-LH100HG**  
100 W mercury  
lamp housing
- U-RX-T**  
Power supply unit  
for xenon lamp
- U-RFL-T**  
Power supply unit  
for mercury lamp

- U-CPA**  
Intermediate attachment for  
conoscopic and  
orthoscopic observation
- U-AN360P-2**  
Rotatable  
analyzer
- U-OPA**  
Intermediate attachment for  
orthoscopic observation
- U-DAL10X**  
Drawing attachment 10x
- U-DA**  
Drawing attachment
- U-DO3**  
Dual observation attachment
- U-SDO3**  
Side by side observation  
attachment

- U-LC\*5**  
Low magnification  
condenser
- U-AC2**  
Abbe  
condenser
- U-SC3**  
Swing-out  
condenser
- U-AAC**  
Achromatic/  
Aplanatic  
condenser
- U-ULC-2**  
Ultra low  
condenser
- U-PCD2**  
Phase/darkfield  
condenser
- U-POC-2**  
Polarizing  
condenser
- U-DCD**  
Darkfield  
condenser,  
dry
- U-DCW**  
Darkfield  
condenser,  
oil

- U-TLO**  
Oil top lens
- U-TLD**  
Dry top lens
- U-UCD8-2**  
8-position  
universal  
condenser
- U-CO1.25X**  
Low magnification  
conversion  
lens for UCD
- U-IFRES**  
Interface for  
coded nosepiece
- U-CBS**  
Control box for  
coded function
- U-HSEXP**  
Hand switch for  
exposure
- DP21-SAL**  
Standalone Connection Kit
- PC (Software)**

**BX53 system diagram**

◆ **WHN10XX-H, CROSS WHN10X**  
Eyepieces  
**U-CT30-2**  
Centering telescope

\* **WHN10X, WHN10X-H, CROSS WHN10X**  
Eyepieces  
**U-CT30-2**  
Centering telescope

▲ **SWH10X-H, CROSS SWH10X, MICRO SWH10X**  
Eyepieces  
**U-CT30-2**  
Centering telescope

**TR-Adapter**

**U-FMT**  
F-mount Adapter  
**U-TMAD**  
T-mount Adapter

**U-SMAD**  
Sony-mount Adapter

**U-BMAD**  
Bayonet-mount Adapter

**U-CMAD3**  
C-Mount Adapter

**U-TVZ**  
Zoom Adapter

**U-TV1X-2**  
TV Adapter

**U-DPCAD**  
Dual port tube with C mounts

**U-CMDPTS**  
C mount adapter for U-DPTS

**U-DPTS**  
Multi double port tube

**U-CMDPTS**  
C mount adapter for U-DPTS

**U-TV0.63XB**  
B4-Mount Adapter

**U-SWETTR-5**  
Super widefield erect image tilting trinocular observation tube

**U-SWTR-3**  
Super widefield trinocular tube

**U-ETR-4\*1**  
Erect image trinocular tube

**U-TR30NIR**  
Trinocular tube

**U-TR30-2**  
Trinocular tube

**U-TTR-2**  
Tilting trinocular tube

**U-BI30-2**  
Binocular tube

**U-TBI-3\*1**  
Tilting binocular tube

**U-GAN**  
Gout analyzer

**U-ANT**  
Analyzer for transmitted light

**U-DICT**  
DIC slider for transmitted light

**U-DICTS**  
Shift DIC slider for transmitted light

**U-DICTHR**  
High resolution DIC slider for transmitted light

**U-DICTHC**  
High contrast DIC slider for transmitted light

**U-ANT**  
Analyzer for transmitted light

**U-TAD**  
Plate adapter

**COMPENSATORS**

**U-FMP**  
Mechanical stage

**U-CST**  
Centering target

**U-SHG**  
Rubber grip  
**U-SHGT**  
Rubber grip

**U-HLD-4, U-HLDT-4**  
Specimen holder

**U-HRD-4, U-HRDT-4**  
Specimen holder

**U-HLS-4, U-HLST-4**  
Specimen holder

**U-SRG2**  
Rotatable graduated stage

**U-SRP**  
Precision rotatable stage


**U-SP**  
Plain stage

**U-SVRO**  
Oil rectangular stage with right-hand control

**U-SVLO**  
Oil rectangular stage with left-hand control

**U-SVRB-4**  
Mechanical stages with right-hand control

**U-SVLB-4**  
Mechanical stages with left-hand control



\*1 Slight vignetting may occur in combination with an additional intermediate attachment or observation method. \*2 Require an additional intermediate attachment or fluorescence illuminator. \*3 Cannot be used with U-TTLBI. \*4 Compatible with FN 22. \*5 Cannot be used with BX3-URA. \*6 Stand is a standard equipment of the U-MDOSV and U-MDO10R3. \*7 An auxiliary lens is equipped.

**CAMERAS**

- U-TV0.25XC**  
0.25X  
C-Mount  
Adapter
- U-TV0.35XC-2**  
0.35X  
C-Mount  
Adapter
- U-TV0.5XC-3**  
0.5X  
C-Mount  
Adapter
- U-TV0.63XC-3**  
0.63X  
C-Mount  
Adapter
- U-TV1XC**  
C-Mount  
adapter 1X  
(XY adjustment)

- U-TBI-3-CLI\*1**  
Tilting binocular tube
- U-TTBI**  
Ergonomic binocular tube
- U-ETBI**  
Ergonomic erect image binocular tube
- U-TTLBI\*2**  
Tilting, telescopic, lifting binocular tube

- U-P6RE**  
Centerable sextuple revolving nosepiece
- U-P4RE**  
Centerable revolving nosepiece
- U-5RE-2**  
Quintuple revolving nosepiece

- U-LLGAD**  
Liquid light guide adapter
- U-LLG150/U-LLG300**  
Liquid light guide (1.5 m/3 m)
- U-HGLGPS**  
Light source

- U-LHEAD\*5**  
Extension adapter for lamp housing
- U-DULHA**  
Double lamp house adapter
- U-AW**  
Motorized attenuator wheel

- U-LH75XEAP0**  
75 W xenon apo lamp housing
- U-LH100HGAP0**  
100 W mercury apo lamp housing
- U-LH100HG**  
100 W mercury lamp housing
- U-RX-T**  
Power supply unit for xenon lamp
- U-RFL-T**  
Power supply unit for mercury lamp

- U-LC\*7**  
Low magnification condenser
- U-AC2**  
Abbe condenser
- U-SC3**  
Swing-out condenser
- U-AAC**  
Achromatic/Aplanatic condenser
- U-ULC-2**  
Ultra low condenser
- U-PCD2**  
Phase/darkfield condenser
- U-POC-2**  
Polarizing condenser

- U-CO1.25X**  
Low magnification conversion lens for UCD
- BX3-UCD8A**  
Motorized universal condenser
- U-UCD8-2**  
8-position universal condenser
- U-TLO**  
Oil top lens
- U-TLD**  
Dry top lens
- U-DCD**  
Darkfield condenser, dry
- U-DCW**  
Darkfield condenser, oil

- U-ECA**  
Magnification changer 2x
- U-ANT**  
Analyzer for transmitted light
- U-KPA**  
Intermediate attachment for simple polarizing observation

- U-CA**  
Magnification changer
- U-DP\*1+3**  
Dual port
- U-DP1XC**  
Dual port 1x

- U-APT**  
Arrow pointer
- U-EPA2**  
Eyepoint adjuster
- U-EPAL-2**  
Eyepoint adjuster

- U-TRU\*1+3**  
Trinocular intermediate unit

- U-TRUS\*1**  
Trinocular intermediate unit

- U-CPA**  
Intermediate attachment for conoscopic and orthoscopic observation
- U-AN360P-2**  
Rotatable analyzer

- U-OPA**  
Intermediate attachment for orthoscopic observation

- U-DAL10X**  
Drawing attachment 10x
- U-DA**  
Drawing attachment

- U-DO3**  
Dual observation attachment

- U-SDO3**  
Side by side observation attachment

- U-MDOV\*6**  
Multi observation side viewer

- U-MDOB3**  
Multi observation body

- U-MDO10B3**  
Multi observation body for 10 persons

- U-MDO10R3\*6**  
Multi observation body for 10 persons

- U-HSEXP**  
Hand switch for exposure

- U-CBS**  
Control box for coded function

- U-HSCBM**  
Hand switch for CBM

- U-IFCBL200**  
Interface cable, 200 cm

- U-CBM**  
Control box

- DP21-SAL**  
Standalone Connection Kit

PC (Software)

### BX43 specifications

|                     |                           |   |
|---------------------|---------------------------|---|
| Microscope frame    | Optical system            | UIS2 optical system   |
|                     | Focus                     | Vertical stage movement: 25 mm stage stroke with coarse adjustment limit stopper, Torque adjustment for coarse adjustment knobs, Stage mounting position variable, High sensitivity fine focusing knob (minimum adjustment gradations: 1 µm)  |
|                     | Illuminator               | Built-in Koehler illumination for transmitted light, light intensity manager switch<br>High color reproductivity LED light source, 6 V 30 W halogen bulb (pre-centered)   |
| Revolving nosepiece |                           | Interchangeable reversed quintuple/coded quintuple/sextuple/septuple/coded sextuple/coded septuple nosepiece  |
| Observation tube    | Widefield (FN 22)         | <ul style="list-style-type: none"> <li>• Widefield tilting, telescopic and lifting binocular</li> <li>• Widefield tilting trinocular • Widefield trinocular • Widefield erect image trinocular</li> <li>• Widefield tilting binocular • Widefield ergo binocular • Widefield binocular</li> </ul>   |
|                     | Super widefield (FN 26.5) | <ul style="list-style-type: none"> <li>• Super widefield trinocular • Super widefield erect image tilting trinocular</li> </ul>   |
| Stage               |                           | Ceramic-coated coaxial stage with left or right hand low drive control: with rotating mechanism and torque adjustment mechanism, optional rubber grips available (Non stick grooved coaxial, plain, rotatable stages are also available)  |
| Condenser           |                           | <ul style="list-style-type: none"> <li>• Abbe (NA 1.1), for 4x–100x</li> <li>• Swing out Achromatic (NA 0.9), for 1.25x–100x (swing-out: 1.25x–4x)</li> <li>• Achromatic Aplanatic (NA 1.4), for 10x–100x</li> <li>• Phase contrast, darkfield (NA 1.1), [phase contrast: for 10x–100x, darkfield: for 10x–100x (up to NA 0.80)]</li> <li>• Universal (NA 0.9), for 1.25x–100x [swing-out: 1.25x–4x, with oil top lens:(NA 1.4)]</li> <li>• Low (NA 0.75), for 2x–100x (Dry)</li> <li>• Ultra low (NA 0.16), for 1.25x–4x</li> <li>• Darkfield dry (NA 0.8–0.92), for 10x–100x</li> <li>• Darkfield oil (NA 1.20–1.40), for 10x–100x</li> </ul> |

### BX46 specifications

|                     |                   |  |
|---------------------|-------------------|--|
| Microscope frame    | Optical system    | UIS2 optical system  |
|                     | Focus             | Fixed low stage nosepiece focus<br>15 mm focus stroke with coarse adjustment limit stop<br>Torque adjustment for coarse adjustment knobs<br>High sensitivity fine focusing knob (adjustment gradations: 1 µm)                                      |
|                     | Illuminator       | Built-in Koehler illumination for transmitted light, light intensity manager switch<br>High color reproductivity LED light source  |
| Revolving nosepiece |                   | Fixed reversed coded quintuple nosepiece   |
| Observation tube    | Widefield (FN 22) | <ul style="list-style-type: none"> <li>• Widefield tilting trinocular • Widefield trinocular</li> <li>• Widefield tilting binocular • Widefield tilting, Telescopic, Lifting binocular • Widefield ergo binocular • Widefield binocular</li> </ul> |
| Stage               |                   | Ceramic-coated coaxial stage with left or right hand low drive control, rotating mechanism and torque adjustment mechanism (Low torque, Plain, Rotating stages are also available)   |
| Condenser           |                   | Built-in condenser (NA 0.9) 1.25x–100x (swing out: 1.25x–2x)   |

### BX53 specifications

|                           |                           |   |
|---------------------------|---------------------------|---|
| Microscope frame          | Optical system            | UIS2 optical system   |
|                           | Focus                     | Vertical stage movement: 25 mm stage stroke with coarse adjustment limit stopper, Torque adjustment for coarse adjustment knobs, Stage mounting position variable, High sensitivity fine focusing knob (minimum adjustment gradations: 1 µm)  |
|                           | Illuminator               | Built-in Koehler illumination for transmitted light, Light preset switch, Light intensity LED indicator,<br>Built-in filters (LBD-IF, ND6, ND25, optional)<br>12 V 100 W halogen bulb (pre-centered)  |
| Revolving nosepiece       |                           | Interchangeable reversed quintuple/coded quintuple/sextuple/septuple/coded sextuple/coded septuple nosepiece  |
| Observation tube          | Widefield (FN 22)         | <ul style="list-style-type: none"> <li>• Widefield tilting trinocular • Widefield trinocular</li> <li>• Widefield tilting binocular • Widefield tilting, Telescopic, Lifting binocular • Widefield ergo binocular • Widefield binocular</li> </ul>  |
|                           | Super widefield (FN 26.5) | <ul style="list-style-type: none"> <li>• Super widefield trinocular • Super widefield erect image tilting trinocular</li> </ul>   |
| Stage                     |                           | Ceramic-coated coaxial stage with left or right hand low drive control: with rotating mechanism and torque adjustment mechanism, optional rubber grips available (Non stick grooved coaxial, plain, rotatable stages are also available)  |
| Condenser                 |                           | <ul style="list-style-type: none"> <li>• Abbe (NA 1.1), for 4x–100x</li> <li>• Swing out Achromatic (NA 0.9), for 1.25x–100x (swing-out: 1.25x–4x)</li> <li>• Achromatic Aplanatic (NA 1.4), for 10x–100x</li> <li>• Phase contrast, darkfield (NA 1.1), [phase contrast: for 10x–100x, darkfield: for 10x–100x (up to NA 0.80)]</li> <li>• Universal (NA 0.9), for 1.25x–100x [swing-out: 1.25x–4x, with oil top lens:(NA 1.4)]</li> <li>• Low (NA 0.75), for 2x–100x (Dry)</li> <li>• Ultra low (NA 0.16), for 1.25x–4x</li> <li>• Darkfield dry (NA 0.8–0.92), for 10x–100x</li> <li>• Darkfield oil (NA 1.20–1.40), for 10x–100x</li> </ul> |
| Fluorescence illuminator  |                           | <ul style="list-style-type: none"> <li>• Multi-purpose coded type (FN 22, 8-position mirror unit turret, 4-position ND slider)</li> <li>• Economical type (FN 26.5, 8-position mirror unit turret)</li> </ul>   |
| Fluorescence light source |                           | 100 W Hg apo lamp housing and transformer, 100 W Hg lamp housing and transformer, 75 W Xe lamp housing and transformer or 130 W Hg light guide illumination   |


The U-CBM is designed for the BX3 use in industrial environments for the EMC performance (IEC61326-1 Class A device). Using it in a residential environment may affect other equipment in the environment.

### BX43/BX46/BX53 common specifications


|                       |  |
|-----------------------|--|
| Operating environment | <ul style="list-style-type: none"> <li>• Indoor use</li> <li>• Ambient temperature : 5 ° to 40 °C (41 ° to 104 °F)</li> <li>• Maximum relative humidity : 80 % for temperatures up to 31° C (88 °F), decreasing linearly through 70 % at 34 °C (93 °F), 60 % at 37 °C (99 °F), to 50 % relative humidity at 40 °C (104 °F)</li> <li>• Supply voltage fluctuations : Not to exceed ±10 % of the normal voltage</li> </ul> |
|-----------------------|--|

**BX43 dimensions**


(unit: mm)

**BX46 dimensions**


(unit: mm)

**BX53 dimensions**


(unit: mm)

**BX53 FL dimensions**


(unit: mm)

**BX53+U-DO dimensions**

(unit: mm)

**BX53+U-MD010 dimensions**

(unit: mm)



- **OLYMPUS CORPORATION is ISO14001 certified.**
- **OLYMPUS CORPORATION is FM553994/ISO9001 certified.**
- **OLYMPUS CORPORATION is MD540624/ISO13485 certified.**
- **Illumination devices for microscope have suggested lifetimes. Periodic inspections are required. Please visit our web site for details.**
- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Images on the PC monitors are simulated.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

---

**OLYMPUS®**

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo, Japan  
**OLYMPUS AMERICA INC.**  
3500 Corporate Parkway, Center Valley, Pennsylvania 18034-0610, U.S.A.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS AUSTRALIA PTY. LTD.**  
31 Gilby Road, Mt. Waverley, VIC 3149, Melbourne, Australia.

**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.  
**OLYMPUS (CHINA) CO., LTD.**  
A8F Ping An International Financial Center, No. 1-3, Xinyuan South Road,  
Chaoyang District, Beijing, China, 100027  
**OLYMPUS KOREA CO., LTD.**  
Olympus-Tower, 114-9 Samseong-Dong, Gangnam-Gu, Seoul, Korea

[www.olympus.com](http://www.olympus.com)

Printed in Japan M1696E-032013



# OLYMPUS

Your Vision, Our Future


LED Illumination Light Sources

FluoLED Range

Transmitted Fluorescence

Fluorescence for All





### FluoLED

LED cassette



### FluoLED MultiFluo

Illuminator with three LED cassettes



## LEDS: THE FUTURE OF ILLUMINATION

LEDs are small solid-state electronic devices which can be made to emit light in specific wavelength bands without the use of colour filters. They have already found widespread acceptance for diverse lighting requirements such as digital clocks, traffic lights and household appliances, and were first used in microscopy for white light illumination in ring lights and illumination bases for stereo microscopes. Recent advances now mean they are powerful enough to be used as sources of fluorescent illumination for microscopy. They have many advantages over arc burner-based systems, such as lower cost of purchase and ownership. Furthermore, they reach operating conditions almost instantaneously and can be switched on and off very quickly. There is no need for neutral-density filters since light intensity can be changed by varying the voltage, without any colour effects as with tungsten bulbs or arc burners.

LEDs offer excellent illumination stability, and are also extremely safe to use and dispose of since they do not explode or have any harmful components. They produce light more efficiently than other light sources and this combined with their long lifetimes (2-year warranties offered on FluoLED units) ensure they offer cost and energy savings.

### FluoLED

The FluoLED system has been specifically designed for use with transparent specimens, such as bacteria and thin sections, which require transmitted fluorescence illumination. The system provides very good signal-to-noise ratios (S/N) with high intensities and is very easy to retrofit to the Olympus CX upright microscope range. With the unit fitted, normal brightfield light can still be used by just sliding out a mirror.

Importantly, the light source is simple and intuitive to use as well as being alignment-free, meaning it can be installed and used instantly. Very few other optical items are required since no excitation or dichroic filters are necessary except for in the MultiFluo system, which uses dichroic filters due to the multiple input wavelengths. As a result, signal separation is very easy to achieve.

All these features make it the perfect entry-level fluorescence system for educational, clinical and routine use – day in, day out.

### Multiple versions


The FluoLED system is available in three versions, all of which attach directly to Olympus CX microscopes to provide transmitted fluorescence whilst still maintaining the microscopes' capacity to perform normal transmitted brightfield illumination. The fixed-wavelength FluoLED Easy Blue system provides a low-cost, entry-level system with consistent intensity. The package includes a 480 nm LED cassette and 510 nm emission filter, making it perfect for blue excitable fluorescent dyes such as FITC. The FluoLED single-channel system provides on/off and intensity control for one of the seven interchangeable LED cassettes, which provides excellent flexibility for fluorescence dye detection. The FluoLED MultiFluo (available for Olympus CX31/CX41 only) provides variable intensity control of two or three of the seven available LED cassettes for multicolour observations.

## Fluorescence for all

The ease of use, low running costs and the adaptation to small laboratory microscopes of the FluoLED systems make fluorescence microscopy available to a new range of user groups. For example, using fluorescent markers in the detection of infections such as tuberculosis could deliver results with higher fidelity and efficiency than standard staining procedures. Furthermore, the on-site investigation at crime scenes could benefit from access to a fluorescence microscope for immediate analysis and would ensure that decisions are made more quickly and therefore more effectively. For the most remote locations, the FluoLED EasyBlue and Single systems can be powered via batteries or even solar power, opening up a new world of possibilities.


## Your fluorescence microscope

All Olympus CX2 microscopes feature peerless UIS2 optics which ensure the highest optical standard. The most cost-efficient fluorescent system available from Olympus consists of the CX21 microscope, which provides the best image flatness in its class, in combination with FluoLED EasyBlue. For additional flexibility, the CX21 can also be combined with the FluoLED Single system. The FluoLED Single and multicolour FluoLED MultiFluo can also be fitted onto the CX31 with its excellent field number (FN20) and integrated 30 W adjustable halogen lamp. For the ultimate FluoLED system, the Single and MultiFluo can be fitted to the highly flexible, modular CX41 system microscope for all inspection and education applications in life science and medicine.



### FluoLED Single

Illuminator with LED cassette



## Specifications

|                               | FluoLED EasyBlue<br>CX21     | FluoLED EasyBlue<br>CX31/CX41 | FluoLED Single<br>CX21      | FluoLED Single<br>CX31/41   | FluoLED MultiFluo<br>CX31/CX41 |
|-------------------------------|------------------------------|-------------------------------|-----------------------------|-----------------------------|--------------------------------|
| <b>LED cassettes</b>          | 480 nm incl.                 | 480 nm incl.                  | 1                           | 1                           | Up to 3                        |
| <b>Intensity control</b>      | Fixed driver                 | Fixed driver                  | 1-channel electronic driver | 1-channel electronic driver | 3-channel electronic driver    |
| <b>Emission filter slider</b> | None, 510 nm LP filter incl. | None, 510 nm LP filter incl.  | 3 positions                 | 3 positions                 | 6 positions                    |
| <b>UV blocking filter</b>     | Not necessary                | Not necessary                 | Included                    | Included                    | Included                       |
| <b>Battery pack</b>           | Optional                     | Optional                      | Optional                    | Optional                    | Not possible                   |
| <b>Solar battery charger</b>  | Optional                     | Optional                      | Optional                    | Optional                    | Not possible                   |

Note: For FluoLED Single and MultiFluo: LED cassettes and emission filters to be ordered separately. Additional dichroic filters must be ordered for FluoLED MultiFluo.  
 Note: FluoLED® is a registered trademark of FRAEN Corporation Srl.

The manufacturer reserves the right to make technical changes without prior notice.


[www.olympus-europa.com](http://www.olympus-europa.com)

**OLYMPUS**

OLYMPUS LIFE SCIENCE EUROPA GMBH

Postfach 10 49 08, 20034 Hamburg, Germany  
 Wendenstraße 14–18, 20097 Hamburg, Germany  
 Phone: +49 (0) 40 23 77 30, Fax: +49 (0) 40 2 37 73 46 47  
 E-mail: [microscopy@olympus-europa.com](mailto:microscopy@olympus-europa.com)

2006-09






New avenues of research are opening in biological and medical fields. As research demands become more specialized and diversified, biological microscopes must offer the capabilities to meet these needs.

Olympus microscopes and their accessories are developed to meet the ever-changing needs of research applications.

Our accomplishments in microscope development date back more than three-quarters of a century. Olympus has accumulated a broad range of advanced optical and precision technologies and we are renowned for our innovative approach to microscopy.

An outstanding example of Olympus ingenuity is the superior UIS2 infinity-corrected optical system employed on the BX2 and IX2 microscopes. Olympus has also won acclaim for its system versatility and broad range of advanced accessories.

Our microscopes are evolving with enhanced performance and operational ease. Olympus continues to answer research demands in the biological and medical field of today and pave the way for future advances with increasingly sophisticated research equipment.



# Contents



Fluoview FV1000 Confocal Laser Scanning Biological Microscope ... 3

Fluoview FV300 Confocal Laser Scanning Biological Microscope ... 3

BX61 Motorized System Microscope ... 4

BX51 System Microscope ... 4

BX41 Laboratory Microscope ... 5

BX45 Ergonomic Laboratory Microscope/  
BX45A Automated Ergonomic Laboratory Microscope ... 5

Accessories for BX2 ... 6-7



BX51WI Fixed Stage Upright Microscope/  
BX61WI Fixed Stage Upright Microscope with Motorized Focusing ... 7

CX41 System Microscope/CX31 Biological Microscope ... 8

CX21 Biological Microscope ... 8

IX81 Motorized Inverted System Microscope ... 9

IX71 Research Inverted System Microscope/  
IX51 Inverted System Microscope ... 9

Accessories for IX2 ... 10

CKX41/CKX31 Inverted Microscopes ... 11

Accessories for CKX41/CKX31 ... 11

ON3 Series Micromanipulators ... 12

BX51-P Polarizing Microscope ... 13

CX31-P Polarizing Microscope ... 13

SZX16/SZX10 Research Stereomicroscope System ... 14

SZX7 Stereomicroscope System ... 15

SZ61/SZ51 Stereomicroscope System ... 15

Accessories for SZX/SZ ... 16

Macro View MVX10 Research macro zoom system microscope ... 17

DP30BW High Sensitivity cooled CCD camera ... 18

DP71 Digital Camera ... 18

DP20 Digital Camera ... 18

DP12 Digital Camera ... 18

UIS2 Series Objectives ... 19-21

UIS Series Objectives ... 20, 22



# FLUOVIEW FV1000

## CONFOCAL LASER SCANNING BIOLOGICAL MICROSCOPE

**UIS2**  
World-leading optics



FV1000+IX81

The FV1000 has an original spectral detection system which uses a high speed diffraction grating combined with a variable slit to deliver superior linear spectral distribution. This enables high-precision, high-resolution, high-speed spectroscopy in observations ranging from milliseconds to hours. SIM (SIMultaneous) Scanner System synchronizes laser light stimulation and confocal imaging. The FV1000 incorporates 2 laser scanners for simultaneous observation and laser light stimulation.

The FV1000 is the most suitable choice of microscope for FRAP, FLIP and photo activation.

|                        |                       | Spectral type fluorescence detector   | Filter type fluorescence detector   |
|------------------------|-----------------------|---|---|
| Laser light            | Visible light laser   | Multi-line Ar laser, HeNe(G) laser, HeNe (R) laser<br>HeNe(G) laser (543nm, 1mW), HeNe (R) laser(633nm, 10mW)<br>Visible light laser platform with implemented AOTF system  |   |
|                        | Violet laser (option) | Laser Diode 405nm, 440nm  |   |
|                        |                       | Laser port: Standard 3 laser ports, VIS - UV - IR   |   |
| Scanning and Detection | Detector module       | Standard 3 confocal Channels.<br>Spectral detection :<br>CH1 and CH2 equipped with independent grating.<br>CH3 with 6 position barrier filter turret  | Standard 3 confocal Channels.<br>CH1 to CH3 each with 6 position barrier filter turret. |
|                        | Scanning method       | Dual galvano mirror scanner (X, Y)  |   |
|                        |                       | Scanning modes: pixel size: 64 x 64 — 4096 x 4096<br>Pixel Dwell time: 2 to 200 microsec with unidirectional<br>0.5 or 1 microsec with fast bidirectional scanning<br>X, Y, T, Z, λ. (any combination)<br>Line scanning: Straight line with free orientation, free line |   |
|                        |                       | Field Number (N.A.) 18  |   |
|                        |                       | Optical Zoom 1X — 50X in 0.5X increment   |   |
|                        |                       | Z-drive Motorized focus module of the microscope, minimum increment 10 nm   |   |
|                        |                       | Transmitted light detector unit External transmitted photomultiplier detector   |   |
| Microscope             | Motorized microscope  | Inverted IX81, Upright BX61, Upright focussing nosepiece & fixed stage BX61WI   |   |
| Optional unit          | SIM Scanner           | 2 Galvano scanning mirrors, pupil projection lens, built-in laser shutter, 1 laser port<br>Fiber introduction of near UV laser diode or visible light laser, Optional: 2nd AOTF laser combiner  |   |

\*Please refer to FV1000 catalog for further details

3

# FLUOVIEW FV300

## CONFOCAL LASER SCANNING BIOLOGICAL MICROSCOPE

**UIS2**  
World-leading optics



FV300+BX61WI

The FV300 gives both individual and group users the right solution to match their research needs and budget.

The system is compatible with the Olympus research range of microscopes offering high resolution confocal sectioning with the ability to conduct time-lapse experiments.

The FV300 offers a wider number of options with the ability to upgrade the systems for the future.

- Highest image quality (12 bit, 2048x2048 pixel resolution) with economical cost.
- Simultaneous capturing of 2 fluorescence and 1 transmitted light detector images.
- Simple, straight optical systems for easy system construction.

|  |  |
|--|--|
| Up to 3 channel *2-channel fluorescence (+1-channel transmitted light)   |  |
| Visible light laser, *Ar laser, *Multi laser, *Green HeNe laser, *Red HeNe laser, Helium Cadmium laser                               |  |
| 1 laser port for visible light laser   |  |
| Manual operating scanning unit   |  |
| 5-position single pinhole turret   |  |
| Scanning unit, *Galvanometer mirror scanner (both X and Y), *Photo detector: PMT   |  |
| Automatic laser control/laser combiner   |  |
| *Each laser light path is equipped with continuously variable neutral density filters or AOTF and a shutter                          |  |
| All laser lines are combined and introduced into via a single fiber optic  |  |
| Intuitive user friendly software   |  |
| Scanning mode:   | 1-dimension: Point scanning<br>2-dimension (space): X-Y, rect, X-Z, linear line-Z and free line-Z<br>2-dimension (time): X-t, linear line-t and free line-t<br>3-dimension (space): X-Y-Z and rect-Z<br>3-dimension (time): X-Y-t, X-Z-t and rect-t<br>4-dimension: X-Y-Z-t and rect-Z-t |
| Image size and scanning speed:   | 256X256(0.45s) - 2048X2048(10.835s)  |
| *High speed mode:  | 512X512(0.25s)   |
| Image acquisition:   | 12bit  |
| Image analysis:  | 2D measurement, image filtering, various image display mode, 3D animation display etc.   |
| Microscope<br>Upright microscope (BX50, BX51, BX61), Inverted microscope (IX70, IX71, IX81), Fixed stage upright microscope (BX61WI) |  |

\*Please refer to FV300 catalog for further details



# BX61

## MOTORIZED SYSTEM MICROSCOPE

**UIS2**  
World-leading optics



Motorized system microscope BX61 in the BX2 series features outstanding reliability in highly advanced microscopy such as three dimensional imaging capture with high-end fluorescence imaging and confocal microscope. Standard features include motorized focusing and light adjustment and a stage escape mechanism. Several key macro microscope operations are available by software-controlled setting, and executed by hand switch, personal computer.

- Uses flexible software which can easily be customized for specific purposes.
- Many operating procedures (including switching between observation methods) are stored as individual macros and allotted to buttons on the microscope, hand switches and keys on the monitor of the personal computer.
- A wide variety of separate modules make expansion very easy.

|                  |  |
|------------------|--|
| Illumination     | Transmitted light 12V100W halogen Koehler illumination   |
| Focusing         | Motorized focus<br>Full stroke: 14mm, minimum fine adjustment: 0.01µm  |
| Observation tube | Widefield binocular (F.N.22), widefield tilting binocular (F.N.22), widefield trinocular (F.N.22), widefield tilting/telescoping binocular (F.N.22), super widefield trinocular (F.N.26.5)                   |
| Nosepiece        | Interchangeable reversed quintuple/sextuple/septuple, motorized sextuple with slider slot for DIC, septuple for DIC/simple POL   |
| Stage            | Ceramic-coated coaxial with left or right hand low drive control, non-stick grooved coaxial, plain, rotatable  |
| Condenser        | Abbe (N.A.1.1), swing out Achromatic (N.A.0.9), Achromatic Aplanatic (N.A.1.4), Universal (N.A.1.4/0.9)  |
| Other features   | Coarse/fine changeover button, stage shunting button, stage up/down button, built-in filters (LBD-IF, ND6, ND25, option)   |
| Accessories      | Motorized fluorescence illuminator, motorized universal condenser, motorized transmitted filter wheel, motorized reflected filter wheel, motorized observation filter wheel, hand switch, control unit, etc. |

\*Please refer to BX51/BX61 catalog for further details

# BX51

## SYSTEM MICROSCOPE

**UIS2**  
World-leading optics



This leading model from the BX2 series offers improved ergonomic and system performance, and is widely used in both routine work and specialized research. It is equipped with Olympus' original UIS2 optic system, and a high-rigidity Y-shape frame with newly refined ergonomics. Excellent features provide the flexibility needed for compliance with a variety of applications. In addition, many kinds of filter sliders and accessories are all designed with multiple openings on the bodies and illuminators.

- Can be combined with a swing-out condenser and a septuple revolving nosepiece with DIC slider to enable continuous observations from 1.25X to 100X.
- A newly developed 8-position universal condenser gives even greater freedom to combine observation methods.

|                  |  |
|------------------|--|
| Illumination     | Transmitted light 12V100W halogen Koehler illumination   |
| Focusing         | Coarse & fine coaxial handle:<br>full stroke: 25mm, minimum fine adjustment: 1µm   |
| Observation tube | Widefield binocular (F.N.22), widefield tilting binocular (F.N.22), widefield trinocular (F.N.22), widefield tilting/telescoping binocular (F.N.22), super widefield trinocular (F.N.26.5) |
| Nosepiece        | Interchangeable reversed quintuple/sextuple/septuple   |
| Stage            | Ceramic-coated coaxial with left or right hand low drive control, non-stick grooved coaxial, plain, rotatable  |
| Condenser        | Abbe (N.A.1.1), swing out Achromatic (N.A.0.9), Achromatic Aplanatic (N.A.1.4), Universal (N.A.1.4/0.9)  |
| Other features   | Coarse/fine changeover button, stage shunting button, stage up/down button, built-in filters (LBD-IF, ND6, ND25, option)   |
| Accessories      | Reflected light fluorescence attachment, DIC attachment, phase-contrast attachment, multi-viewing attachment, polarizing attachment, etc.  |

\*Please refer to BX51/BX61 catalog for further details

# BX41

## LABORATORY MICROSCOPE

**UIS2**  
World-leading optics



This key model in the BX2 series offers excellent cost-performance and meets a wide variety of needs in the routine inspection of laboratories and clinics. The extensive range of features benefits from consistently simple operation, enhanced by Olympus' original and ergonomically-advanced Y-shape frame design and UIS2 optical system. This model allows easy system expansion: by combining multiple units together, it can meet the majority of general working requirements.

- Newly developed 30W halogen light source provides as much brightness as a conventional 50W halogen lamp.
- Equipped with an easy-to-operate rackless stage with no protrusion of the X-axis guide.
- The front-located power switch and light intensity volume control make it possible to operate the microscope with the arms resting on the desk.
- A 3-filter cassette enables quick, easy exchange.
- The detachable revolving nosepiece allows quick exchange of objectives to suit different observation methods.

|                  |  |
|------------------|--|
| Illumination     | Transmitted light 6V30W halogen Koehler illumination   |
| Focusing         | Vertical stage movement: 25mm; Stage stroke with coarse adjustment limit stop; Torque adjustment for coarse adjustment knobs; Stage mounting position variable; High sensitivity fine focusing knob (adjustment graduations 1µm) |
| Observation tube | Widefield binocular (F.N.22), widefield tilting binocular (F.N.22), widefield trinocular (F.N.22), widefield ergo binocular (F.N.22) super widefield trinocular (F.N.26.5)   |
| Nosepiece        | Interchangeable reversed sextuple/quintuple  |
| Stage            | Ceramic coated coaxial with right or left hand low drive control and rotating mechanism, non stick grooved coaxial, plain, rotatable   |
| Condenser        | Abbe(N.A.1.1), swing-out Achromat (N.A.0.9), Achromatic Aplanat (N.A.1.4), phase-contrast/darkfield (N.A.1.1), darkfield dry (N.A.0.8-0.92), darkfield oil (N.A.1.20-1.40), ultra low (N.A.0.16)                                 |
| Other feature    | Light pre-set switch   |
| Accessories      | Reflected light fluorescence attachment, DIC attachment, phase-contrast attachment, multi-viewing attachment, etc.   |

\*Please refer to BX41/BX45/BX45A/BX51 catalog for further details

5

# BX45/BX45A

## ERGONOMIC LABORATORY MICROSCOPE AUTOMATED ERGONOMIC LABORATORY MICROSCOPE

**UIS2**  
World-leading optics



These high-performance laboratory microscopes are equipped with UIS2 optics that provide excellent image quality, and the same ergonomically advanced Y-shape frame as other models in the BX2 series.

- Use no-cover objectives which are suitable for observing smear specimens (e.g. blood) with no cover glass.
- 2x, 4x, 10x and 20x objectives equipped with ND filter enable the same level of brightness even if the magnification is changed from 2x to 40x. No brightness adjustment (e.g. of light intensity) is required.
- Employ a rackless stage providing precise movement and smooth stopping with no protrusion of the X-axis guide.
- A convenient marking unit is provided.

### BX45

Equipped with a 3-position universal condenser which complies with brightfield, phase-contrast and darkfield observations, and a gout inspection analyzer that uses simple polarizing observation. Enables optimum illumination for observations from 1.25X to 40X.

### BX45A

Equipped with a low stage position fixed at 128mm from the desktop, a revolving nosepiece with variable up/down movement, and a motorized 2-position revolving nosepiece which can make quick changes between objectives.

|                  |  |
|------------------|--|
| Illumination     | Transmitted light 6V30W halogen Koehler illumination   |
| Focusing         | Fixed low stage nosepiece focus  |
| Observation tube | Widefield binocular (F.N.22), widefield tilting binocular (F.N.22), widefield trinocular (F.N.22), widefield ergo binocular (F.N.22) |
| Nosepiece        | Fixed motorized 2-position (BX45A)<br>Fixed reversed quintuple (BX45)  |
| Stage            | Ceramic-coated coaxial with right or left hand low drive control, rotating and torque adjustment mechanism, plain, rotation          |
| Condenser        | Brightfield N.A.0.9 (BX45A)<br>Fixed 3-position universal condenser N.A.0.9 (BX45)   |
| Other features   | One-touch marking unit, Light preset switch, Torque adjustment   |

\*Please refer to BX41/BX45/BX45A/BX51 catalog for further details

## Accessories for BX2



### **BX-RFAA** Motorized fluorescence illuminator

Up to 6 fluorescence mirror units can be attached simultaneously. Mirror unit can be exchanged automatically with corresponding shutter adjustment.



### **U-AFP1** Auto focus unit

Maintains steady auto focusing with 1.25X to 100X objectives. All observation methods are applicable except phase-contrast. Combined with BX61TRF and BX62TRF. Personal computer and adapter (BX2-UCB) are necessary to use this unit.



### **BX-RFA, BX-URA2** BX fluorescence illuminator BX reflected light illuminator

Fully integrated into the microscope arm, two illuminators add to the stability of the imaging platform. BX-RFA for research needs and BX-URA2 for routine observation.



### **U-RSL6, U-RSL6EM** 6-position filter slider

Equipped with three single, two dual and one triple band exciter, the 6-position filter slider provides all excitation modes typically wanted with just one filter set.



### **U-FWR, U-FWO, U-FWT** Filter wheels

Motorized exchange of 6 filters. Three kinds of filters can be attached simultaneously: U-FWR ( $\phi 32$ , 25) for excitation, U-FWO ( $\phi 32$ , 25) for absorption and U-FWT ( $\phi 32$ ) for transmitted light.



### **U-UCD8-2** 8-position universal condenser

The universal condenser simultaneously accepts up to 6 DIC prisms and 2 other optical components at maximum. The condenser numerical aperture of 0.9 or 1.4 (oil) can be selected through interchangeable top elements.



### **U-RFSS** Rectangular field stop

Designed for use with CCD cameras, prevents photobleaching of the specimen outside of the imaging area.



### **U-EXBAG, U-EXBAUB, U-EXBAUG** Excitation balancer

Used singly or in tandem, the excitation balancer curtails the individual excitation bandwidths of the fluorochromes under observation.



### **U-UCD8A-2** Motorized universal condenser

Different combinations of designated optical components allow correspondence with various kinds of transmitted light observations. Automatic control of optical component exchange, top lens swing in/out and aperture iris diaphragm.



### **U-PCD2** Phase/darkfield condenser

High contrast phase imaging allows close observation of the cell interior and of live bacteria. Standard brightfield and dry darkfield as well as simultaneous reflected light fluorescence observations are possible.



### **U-KPO** Simple polarizing attachment

Simple polarizing observation can be accomplished with the combination of U-KPA intermediate attachment for simple polarizing observation, U-ANT analyzer for transmitted light and U-POT polarizer.

## Accessories for BX2

### U-SDO3, U-MDO10B3

#### Multi observation body

A single image can be viewed with the same orientation and brightness by up to ten persons simultaneously to facilitate discussion between researchers, and for training and education. No additional power unit is required since the power unit of the arrow pointer is incorporated in the main body.

•Also available for five persons (U-MDO) and two persons in face-to-face (U-BDO).



BX51+U-SDO3


BX51+U-MDO10B3

7

## BX51WI/BX61WI

### FIXED STAGE UPRIGHT MICROSCOPE/ FIXED STAGE UPRIGHT MICROSCOPE WITH MOTORIZED FOCUSING

UIS  
UNIVERSAL  
INFINITY SYSTEM



Designed for neuroscience and cell biology applications, the BX51WI offers front focus operation and a complete absence of vibration, even when switching the objectives on the nosepiece and filters on the turret and in the optical path of the intermediate tubes. The availability of two-wavelength IR-DIC (775nm and 900nm) has huge advantages for researchers performing electrophysiological work using brain slice samples. Our unique two-position nosepiece can easily switch between objectives, and prevent air bubbles forming. To avoid having to change objectives, Olympus offers a single objective (XLUMPLFL20X/W) for low to high magnifications, combined with an intermediate magnification changer. The BX51WI provides dramatically bright fluorescence images, equivalent to our BX2 models.

|                  |   |
|------------------|---|
| Illumination     | Transmitted light 12V100W halogen Koehler illumination (BX51WI)<br>Transmitted light 12V100W halogen Koehler illumination;<br>Light adjustment: less than DC2V~12V (continuous adjustment)<br>Brightness adjustment, light preset switch (BX61WI)   |
| Focusing         | Nosepiece focus by roller guide (rack & pinion);<br>Stroke per rotation: fine: 0.1mm, coarse:15mm;<br>Maximum stroke: 25mm; Coarse lower limit stopper mechanism,<br>Torque adjustment mechanism for coarse focus (BX51WI)<br>Motorized focusing using stepping motor and ball screw<br>Nosepiece focus by cross roller guide:<br>Minimum graduation: fine: 1µm (Sensitivity 1µm)<br>Resolution 0.01µm: Maximum stage movement speed 3mm/s<br>Stroke: 25mm, stage escape mechanism (BX61WI) |
| Observation tube | Trinocular (F.N.22), erect image trinocular (F.N.22),<br>double port magnification change unit (F.N.22)   |
| Nosepiece        | Swing, slide, single, position, swing-slide   |
| Stage            | Mechanical, bridge  |
| Condenser        | 8-position universal, long working distance oblique,<br>long working distance DIC, swing-out  |

\*Please refer to BX51WI/BX61WI catalog for further details

# CX41/CX31

## SYSTEM MICROSCOPE/ BIOLOGICAL MICROSCOPE

**UIS2**  
World-leading optics



The CX41/CX31 offer extended capabilities to match a wide range of applications from routine clinical work to educational use. They not only feature an ergonomically designed frame for maximum operating comfort and enhanced rigidity, but also offer the convenience of extra-bright illumination.

### CX41

Featuring powerful 6V30W halogen Koehler illumination and outstanding flat images in this class of microscope, the CX41 is applicable for a wide range of observation methods and photomicroscopy.

### CX31

An ergonomically designed frame and bright 6V30W halogen illumination make the CX31 ideal for routine clinical work and educational applications.

|                  |   |
|------------------|---|
| Illumination     | Built-in transmitted Koehler illuminator<br>6V30W halogen bulb<br>100-120V/220-240V~ 0.85/0.45A 50/60Hz   |
| Focusing         | Stage height movement by roller guide (rack & pinion)<br>Stroke per rotation: 36.8mm<br>Full stroke range: 25mm<br>Tension adjustment on coarse focus adjustment knob<br>Upper limit stopper (CX41)<br>Upper limit stopped by simplified pre-focusing dial (CX31) |
| Observation tube | Binocular/Tilting binocular/Trinocular  |
| Nosepiece        | Fixed quadruple nosepiece with inward tilt  |
| Stage            | Size: 188(W) X 134 (Y) mm<br>Traveling range: 76mm(X) x 50mm(Y)<br>Specimen holder: Double slide holder   |
| Condenser        | Abbe condenser, with built-in daylight filter(CX31 only)  |
| Accessories      | Dual-observation attachment, phase-contrast attachment, drawing attachment, simple polarizing attachment, digital camera adapter etc.   |

\*Please refer to CX41, CX31 catalogs for further details

# CX21

## BIOLOGICAL MICROSCOPE

**UIS2**  
World-leading optics



The CX21 demonstrates the ideal combination of advanced performance and operability for multiple inspection and educational purposes in the medical field. Incorporating the UIS2 optical system and employing Plan objective lenses as standard, it delivers class-leading standards of image clarity and flatness in a wide range of observation methods. Other characteristics include excellent durability, and ergonomic design features to reduce fatigue during long observations. To maintain performance in any working environment, an effective anti-fungal treatment is applied to the objectives, eyepieces and microscope tube.

|                      |  |
|----------------------|--|
| Optical system       | UIS2 (Universal Infinity System) optical system  |
| Illumination System  | Built-in transmitted illumination system<br>6V20W halogen bulb<br>100-240V 50/60Hz universal voltage   |
| Focusing             | Stage height movement (coarse movement stroke 20mm)<br>Fine focus graduation: 2.5µm  |
| Revolving nosepiece  | Fixed quadruple nosepiece  |
| Stage                | Wire movement mechanical fixed stage: 120 X 132mm<br>Traveling range: 76mm(X) X 30mm(Y)<br>Single specimen holder  |
| Observation tube     | 30° inclined binocular tube<br>Interpupillary distance adjustment range 48-75mm  |
| Condenser            | Abbe type with aperture iris diaphragm N.A.: 1.25  |
| Objective lens       | Plan Achromatic objectives (anti-fungus)<br>4X N.A.: 0.10 W.D.: 18.5mm<br>10X N.A.: 0.25 W.D.: 10.6mm<br>40X N.A.: 0.65 W.D.: 0.6mm<br>100X N.A.: 1.25 W.D.: 0.13mm (option) |
| Eyepiece (10X)       | Field Number (F.N.): 18 (anti-fungus)  |
| Optional accessories | Mirror unit, 15X eyepiece (F.N. 12, anti-fungus), cord rest, wooden storage box, filar micrometer, wire pointer, filter holder, darkfield stop                               |

\*Please refer to CX21 catalog for further details

# IX81

## MOTORIZED INVERTED SYSTEM MICROSCOPE

**UIS2**  
World-leading optics



This model allows researchers to customize the motorized system according to their own specific purpose, with operating control handled from the front. By using special software via a personal computer, it is also possible to exercise accurate control of multi-dimensional analyses, ranging from 2D to 6D. The full range of IX81 performance functions, including observations, measurements and manipulation, can be monitored via the numerous input/output ports, which allow connection of various kind of light sources and motorized modules.

- Purpose-selectable motorized units and easy operation right by the operator's hand.
- Multi-dimensional analysis by PC control.
- Maximum installation of experimental equipment and minimum layout limitations.
- Sharp, fade-free fluorescent images and faster observations.
- Optimized resolution and contrast in Nomarski DIC observation, for both thick- and thin-cell specimens.
- Combining different light sources and video systems to obtain even clearer images.
- Prolonged active-cell observation with highly reliable data.
- Special microscope body for FV1000/300 is available.

|                  |  |
|------------------|--|
| Illumination     | Transmitted Koehler light 12V100W halogen  |
| Focusing         | Motorized focus: Stroke: 9mm Resolution: 0.1µm   |
| Observation tube | Tilting binocular (F.N. 22), trinocular (F.N. 22)  |
| Nosepiece        | Motorized sextuple with simple waterproof mechanism  |
| Stage            | Cross with flexible right handle, plain, mechanical  |
| Condenser        | Motorized long working distance universal, long working distance universal, DIC, mid long working distance, 8-position universal, ultra long working distance universal      |
| Other features   | Video port for primary image, integrated magnification change lens   |
| Accessories      | Motorized/manual reflected light fluorescence attachment, DIC attachment, external power supply unit, side-viewing attachment, incubator, heat stage, micromanipulator, etc. |

\*Please refer to IX71/IX81 catalog for further details

# IX71/IX51

## RESEARCH INVERTED SYSTEM MICROSCOPE/ INVERTED SYSTEM MICROSCOPE

**UIS2**  
World-leading optics

**eco**  
ECO-PRODUCTS  
IX51



Designed to provide the high performance and versatility needed by researchers involved in live cell experiments, the IX71 offers highly precise temperature control and resistance to heat and vibration, enabling work on live cells with much less risk of damage and reducing the incidence of failure in prolonged experiments.

A 30W illumination pillar type IX51, which has no intermediate magnification changer, is also available.

- More free space and a better working environment, with flexible use of several cameras and light sources.
- Easy front operation allows auxiliary equipment to be placed near the microscope.
- Flexible system expansion allows many different fluorescence applications without major remodeling.
- Obtaining high-quality, purpose-specific images with different cells and different types of container.
- Rigid construction and comprehensive system features to analyze time lapse changes in highly active cell conditions.
- Special microscope body for FV1000/300 is available.

|                  |  |
|------------------|--|
| Illumination     | Transmitted light 12V100WHAL for IX71, 6V30WHAL for IX51   |
| Focusing         | Vertical movement of nosepiece (stage fixed); coarse & fine coaxial handle; full stroke: 9mm; minimum fine adjustment: 1.2µm   |
| Observation tube | Tilting binocular (F.N. 22), binocular (F.N. 22), trinocular (F.N. 22)   |
| Nosepiece        | Sextuple, simple waterproof mechanism  |
| Stage            | Cross with flexible right handle, cross with short left handle, plain, mechanical  |
| Condenser        | Motorized long working distance universal, long working distance universal, DIC, mid long working distance, 8-position universal ultra long working distance universal |
| Other features   | Video port for primary image, integrated magnification change lens   |
| Accessories      | Reflected light fluorescence attachment, DIC attachment, external power supply unit, side-viewing attachment, incubator, heat stage, micromanipulator, etc.            |

\*Please refer to IX71/IX81 catalog and IX51 catalog for further details

## Accessories for IX2

### TIRFM Total internal reflected fluorescence microscopy

An exclusive high N.A.(1.65, 1.45) objective and reflected light illuminator are provided, allowing exchange between evanescent wave and normal fluorescence observation.



IX71 TIRFM version



### U-DPCAD Dual port with C-mount adapter

This double port tube allows the attachment of two cameras (both primary images).



### IX2-RFA Fluorescence illuminator

Can be mounted with six different mirror units. An original mirror unit can be tailor-made from generic mirror units. UV cut filter integrated.



### IX2-RFAL L-shaped fluorescence illuminator

Allows easy centering and AS/FS operation from the front and also permits attachment of a large format camera to the back port.



### IX2-GS Gliding stage

To follow the quick movement of *caenorhabditis elegans*, this stage is designed to move smoothly and freely throughout the plane.



### IX2-LWUCD Long working distance universal condenser

The new DIC system is especially effective in obtaining high-contrast, high-resolution images in 20X and 40X DIC observations.



### IX2-MLWCD Mid long working distance condenser

The relief contrast condenser is designed to produce contrast and shading effects, similar to DIC, yet within the confines of plastic sample vessels.



### IX2-DICD+IX2-TLW DIC condenser + water top lens for DIC

By combining the IX2-TLW top lens, DIC prism and DIC slider, this model provides excellent operability in injection and patch clamping operations.



### IX2-TVRAc Motorized bottom port unit with C-mount

Entirely aberration-free primary images from UIS2 objectives are directed to C-mount CCD camera.



### IX2-RFACA Motorized fluorescent cube turret

Accepts up to 6 fluorescence filter cubes, making it easy to switch between them during fluorescence observation of multi stained specimens. (Manual shutter included.)



### IX2-LWUCDA2 Motorized long working distance universal condenser

Simultaneously accepts up to 6 optical components at maximum. Motorized exchange through PC possible.

# CKX41/CKX31

## INVERTED MICROSCOPES

**UIS2**  
World-leading optics



CKX41/CKX31 are designed to check the viability of cultured cells more quickly and efficiently. Its unique, centering-free phase annulus (common for 10X, 20X and 40X) provides faster phase contrast observations with no need for adjustment. The solid frame has a compact, space-saving design which is ideal for standard workbench surfaces.

### CKX41

Observation tube is exchangeable, a trinocular tube is also mountable. A glass stage insert plate provides quick recognition of objectives.

### CKX31

The CKX31 is a standard type with a fixed binocular tube and a powerful 6V30W halogen illumination and ideal for routine cultured cell observation.

|                  |  |
|------------------|--|
| illumination     | Transmitted light 6V30W halogen illumination   |
| Focusing         | Vertical movement of nosepiece (stage fixed); coarse & fine coaxial handle; full stroke: 9mm; minimum fine adjustment: 1.2µm                 |
| Observation tube | Binocular (CKX31 frame/CKX41 frame with tilting binocular tube)<br>Trinocular (CKX41 frame with trinocular tube)                             |
| Nosepiece        | Fixed quintuple  |
| Stage            | Plain stage (160mm X 240mm), attachable mechanical stage   |
| Condenser        | Detachable long working distance condenser (N.A. 0.3, W.D. 72mm)   |
| Accessories      | Glass stage insert plate, eyepoint adjuster, relief contrast system, Terasaki holder, 35mm dia., petri dish holder, slide glass holder, etc. |

\*Please refer to CKX41/CKX31 catalog for further details

11

## Accessories for CKX41/CKX31



### IX2-SLP

#### Phase contrast slider (pre-centered)

Centering-free type phase slider. A common phase annulus for 10x, 20x and 40x enables fast and easy operation for routine use.



### IX2-SL

#### Phase contrast slider (centerable)

Centerable type phase slider. The centering unit for phase annuli is available for precise adjustment.



### CKX-RCD

#### Relief Contrast condenser

The reduced halo of the relief contrast improves information on the the cell's interior. Maintains the same shadow direction even if the magnification is changed.



### CKX-RFA

#### Fluorescence illuminator

Provides fluorescence observation (B and G excitation) for CKX41.



### CK40-CPG30

#### Glass stage insert plate

Easy recognition of objectives. An insert with ø30mm opening is available (CKX41 only).



### IX2-BCTP

#### Hemacytometer holder

The mechanical stage offers excellent inspection performance with hemacytometer holder or other micro plates.



# ON3 Series

## MICROMANIPULATORS

As a joint development with Narishige Scientific lab, the ON3 micromanipulators offer high-precision and easy operation for IVF, injection and physiological experiments.

### ON3-99D (1:1)

#### Oil-hydraulic micromanipulator system

The ON3-99D consists of a pair of three-axis motorized positioners, drop handle joystick micromanipulators and UT-D universal joints. The ON3-99D is ideal for ICSI applications.


Two types of injectors, the IM-9B microinjector (for sperm injection) and the IM-9C pneumatic injector (for oocyte holding) are available separately.

#### Drop handle joystick

- Fine movement range: 10mm
- Full rotation of knob: 250 $\mu$ m
- Minimum graduation: 2.5 $\mu$ m
- Joystick (for X and Y movement): 400 $\mu$ m max. (movement ratio and lever tension adjustable)

#### Motorized positioner

- Coarse movement range: 23mm (movement speed adjustable)



### ONO-301D

#### Drop handle joystick micromanipulator (1:1)

Thanks to its symmetrical design, this micromanipulator can be attached on either the right or the left side. In combination with the optional universal joint (UT-D) and return mechanism (UT-R), it also provides a pipette return function.

- Accessory: IP plate
- Fine movement range: X, Y and Z axes 10mm
- Full rotation of knob: 250 $\mu$ m
- Minimum graduation: 2.5 $\mu$ m

Photo: ONM-2D, ONO-301D, UT-D and UT-R configuration

### MHW-3

#### Three-axis water hydraulic drum type micromanipulator

The MHW-3 is useful for physiological experiments as well as fine injection or suction. The solid design including mechanical coarse positioner provides reliable stability and precise movement. The optional MHW-4 single axis micromanipulator is mountable.

- Fine movement range: 2mm
- Coarse movement range: 30mm
- Full rotation of knob: 50 $\mu$ m
- Minimum graduation: 0.2 $\mu$ m

# BX51-P

## POLARIZING MICROSCOPES

**UIS2**  
World-leading optics



This series employs UIS2 optics to achieve unsurpassed performance in polarized light observation. These units deliver optimum compensation for optical aberrations to achieve images of unprecedented sharpness. Six compensators are available to allow observation and measurement at various retardation levels.

- Conoscopic/orthoscopic version and orthoscopic version available.
- 6 different kinds of compensators are available for BX51-P.
- Accessories and video/camera system of BX2 Series mountable on BX51-P.

|                              |  |
|------------------------------|--|
| Illumination                 | Transmitted light 12V100W halogen Koehler illumination   |
| Polarizing intermediate tube | Swing-out focusable Bertrand lens with slot for 360° rotatable analyzer for conoscopic & orthoscopic observation (U-CPA)                               |
| Test plate                   | 1 wavelength (1λ), 1/4 wavelength (1/4λ)   |
| Compensators                 | Berek, Senarmont, Brace-Koehler, quartz wedge, etc. (6 types available)  |
| Focusing                     | Coarse & fine coaxial handle; full stroke: 25mm; minimum fine adjustment: 1μm  |
| Observation tube             | Trinocular (F.N. 22)   |
| Nosepiece                    | Detachable quadruple nosepiece with centering adjustment function  |
| Stage                        | Circular rotatable stage with centering adjustment function and attachable mechanical stage. 360° graduated in 1° increments, lockable in any position |
| Condenser                    | Achromat strain-free condenser with built-in 360° rotatable polarizer (N.A. 0.18-0.9)  |

\*Please refer to BX51-P catalog for further details

13

# CX31-P

## POLARIZING MICROSCOPE

**UIS2**  
World-leading optics



The CX31-P is a high-quality polarizing microscope that's ideal for training, with the wide-ranging functions and superior durability required in every field of research.

Its excellent optical performance is matched with the versatility to meet the demands of many different kinds of applications, from double-refraction examination of the structure and characteristics of transparent specimens to complex analyses of rocks, fibers, macromolecules and new materials.

|                                     |   |   |
|-------------------------------------|---|---|
| Optics                              | Objective   | ACHN-P series, UPLFL-P series   |
|                                     | Eyepiece  | WHN10X, WHN10X-H, CROSSWHN10X, WHB10X3, WHB10X2-H   |
| Observation tube                    | Binocular   | U-BI30P, U-CBI30-2  |
|                                     | Trinocular  | U-TR30-2, U-CTR30-2   |
| Conoscopic intermediate tube (U-PA) | Bertrand lens   | Incorporated, detachable, focusable   |
|                                     | Changeover between orthoscopic/conoscopic observation | Engage or disengage of Bertrand lens<br>Position: ● IN<br>Position: ○ OUT   |
|                                     | Analyzer  | Incorporated, detachable, 180° rotatable, lockable in any position, 2° increments, minimum retardation resolution 6', using vernier scale |
|                                     | Slot for compensators                                 | Tint plate (U-TP530), 1/4 wavelength retardation plate (U-TP137) and various compensators attachable                                      |

\*Please refer to BX31-P catalog for further details

# SZX16/SZX10

## RESEARCH STEREO MICROSCOPE SYSTEM



SZX16 with fluorescence system

SZX10

The system modularity allows users to create the application dedicated configurations they need. Offering optical and mechanical excellence and stability, and a wide range of modularity components, the SZX2 Series is today's first choice in research stereo microscopy.

### SZX16

Offering a zoom ratio of 1:16.4, the SZX16 is ideal for the most demanding applications. New SDF objective lenses provide the highest NA with 900lp/mm resolution. Optimum specimen viewing from large field overview to microstructure, along with instant zoom function to select observation points, is assured.

### SZX10

A zoom ratio of 1:10 is suitable for operations like specimen selection or dissection. SZX10 provides wide viewing and assures fewer oversights while relieving fatigue. Choose from a wide range of accessories to suit your sample needs.

|                    | SZX16   | SZX10  |
|--------------------|---|--|
| Optical system     | Telescope type system   |  |
| Zoom range         | 0.7x-11.5x (zoom ratio 1: 16.4),<br>Click stop equipped (releasable)  | 0.63x-6.3x (zoom ratio 1: 10.0)  |
| Aperture diaphragm | Built-in  |  |
| Total mag. range   | 2.1x-690x   | 3.15x-378x   |
| Working distance   | 141 (with SDFPLFLO.3x) -<br>20mm (with SDFPLAPO2xPFC)   | 171 (with DFPL 0.5x) -<br>33.5mm (with DFPL2x-4)   |
| Observation tube   | SZX2-TTR/SZX2-TTRPT: tilting trinocular, 5 - 45° variable inclination<br>SZX2-TR30/SZX2-TR30PT: 30 degree trinocular, 30° inclination   | SZX-BI30: 30° binocular,<br>30° inclination<br>SZX-BI45: 45° binocular,<br>5 - 45° variable inclination<br>SZX-TBI: tilting binocular,<br>5 - 45° variable inclination |
| Objective          | SDFPLFLO.3x, SDFPLAPO0.5xPF,<br>SDFPLAPO0.8x, SDFPLAPO1xPF,<br>SDFPLAPO1.6xPF, SDFPLAPO2xPFC  | DFPLO.5x-4, DFPL0.75x-4<br>DFPLAPO1x-4, SZX-ACH1x,<br>DFPLAPO1.25x, SZX-ACH1.25x-2<br>DFPL1.5x-4, DFPL2x-4   |
| Eyepiece           | WHN10x-H (FN 22)<br>WHSZ15x-H (FN 16), WHSZ20x-H (FN 12.5), WHSZ30x-H (FN 7)  | WHSZ10x-H (FN 22)  |
| Focusing           | SZX2-F0: Focusing unit, coarse handle stroke 80mm<br>SZX2-FOF: Fine focusing unit, coarse handle stroke 80mm,<br>fine handle stroke 80mm<br>SZX2-FOFH: Fine focusing unit for heavy loading, stroke 80mm<br>fine handle stroke 80mm<br>SZX-FOA2: Motorized focus unit, focusing stroke 75mm |  |
| Accessories        | Fluorescence illuminator, coaxial illuminator, light beam splitter, revolving nosepiece, large stage plate, stage adapter, etc.   | Eyepoint adjuster, arrow pointer,<br>drawing attachment, side by side<br>discussion tube, etc.   |

\*Please refer to SZX16/SZX10 catalog for further details

### SZX2-ILLB

#### High-level transmitted light illumination base

Provides effective contrast from oblique illumination and easily selected "High" and "low" contrast settings. Light volume and color temperature are adjusted by means of built-in filters (LBD/ND).



### SZX2-ILLD

#### Brightfield/darkfield transmitted light illumination base

Enables darkfield observation under illumination twice as bright as conventional models. Flat and thin specimens like brain tissue slices are vividly displayed on a black background.



### SZX2-ILLT

#### Slim LED transmitted light illumination base

With a slim design of 41mm, this transmitted light illumination base has a lower height to enable a low eyepoint and easy access to base-mounted samples during observation and operation. The LED 4-position turret enables contrast adjustment between brightfield, oblique, and darkfield illumination with a simple turn.



# SZX7

## STEREOMICROSCOPE SYSTEM



The SZX7 minimizes strain and fatigue while fulfilling the key mission of Olympus microscope designers — to provide the optimal image for any specimen. The clear, accurate performance of the high-level Galilean optical system is complemented by a distortion free objective lens series with maximum numerical aperture. The SZX7 microscope body is manufactured using newly developed lead-free optics, demonstrating Olympus' commitment to protect the environment.

|                      |  |
|----------------------|--|
| Optical system       | Galilean type optical system   |
| Zoom microscope body | Zoom range 0.8x-5.6x (zoom ratio 1:7), Lead-free   |
| Observation tube     | 1. SZX-BI45: Binocular, 45° inclination<br>2. SZX-TBI: Binocular, 5°-45° variable inclination<br>3. SZX2-TR30: Trinocular, 30° inclination<br>4. SZX2-TR30PT: Trinocular, 30° inclination<br>All observation tubes: Lead-free<br>Interpupillary distance adjustable range: 50 to 76 mm |
| Objective            | DFPL0.5x-4, DFPL0.75x-4, DFPLAPO1x-4, SZX-ACH1x, DFPLAPO1.25x, SZX-ACH1.25x-2, DFPL1.5x-4, DFPL2x-4<br>All objectives: lead-free   |
| Eyepieces            | "ComfortView" WHSZ series<br>All eyepieces: Lead-free  |

\*Please refer to SZX7 catalog for further details

# SZ61/SZ51

## ZOOM STEREOMICROSCOPE



Incorporating new improvements to the highly-regarded Greenough optical system, the SZ61 and SZ51 successfully meet the demand for a variety of observation and documentation options in a genuinely compact microscope design. Clear, sharp image reproduction is matched by new ergonomic design elements which maximize comfort and ease of use. The SZ61 and SZ51 microscope bodies are manufactured using newly developed lead-free optics, demonstrating Olympus' commitment to protect the environment.



SZ61/SZ51

SZ61TR

SZ61-60/SZ51-60

|                           | SZ61  | SZ61-60 | SZ61TR                  | SZ51 | SZ51-60 |
|---------------------------|---|---------|-------------------------|------|---------|
| Optical system            | Greenough type optical system   |         |                         |      |         |
| Zoom ratio                | 1:6.7   |         |                         | 1:5  |         |
| Working distance          | 110mm   |         |                         |      |         |
| Tube inclination angle    | 45°   | 60°     | 45°                     | 60°  |         |
| Video camera adaptability | —   |         | C-mount (0.5x built in) | —    |         |
| Optical component         | Lead-free   |         |                         |      |         |
| Auxiliary objective       | Mounting by screwing into the thread at the bottom of frame (M48 threadx0.75) |         |                         |      |         |
| Eyepieces                 | "ComfortView" WHSZ series<br>All eyepieces: Lead-free                         |         |                         |      |         |

\*Please refer to SZ61/SZ51 catalog for further details

## Accessories for SZX/SZ

### SZ2-ILST

#### LED illuminator stand

The world's first LED stand features a thin design to keep sample positions low and to optimize operability. Simultaneous transmitted and reflected light are available on this stand. LED light offers both long lifetime and constant color temperature at any intensity.



### SZ2-ILA

#### Transmitted illumination attachment

Used with the SZ2-ST stand, this cost-effective illumination stand provides bright, uniform illumination from low to high magnifications. Tilttable mirror provides direct and oblique illumination for low contrast specimen. Available 22W and 100W lamphouses provide necessary power for a variety of illumination needs.



### SZ2-LGDI

#### Interlock dual light guide

Standard oblique semi-rigid fiber optic light guide. The light source position on the rear side of the stand saves desk space.



### SZ2-LGSF

#### Flexible light guide

A single fiber optic guide is fixed at the back of the objective so as not to disturb microscope operation.



### SZ2-STU3

#### Table clamp stand

The versatile SZ2-STU3 can be fixed to the side of a desk to expand the working area. Also, this stand exhibits operational versatility in anatomical and electrophysiological experimentation.

- For use with SZ2-STB1 bonder arm



### SZ2-LGR

#### Ring light guide

Used with the GB illumination system, this ring light guide provides bright and uniform images.



### SZ2-STU2

#### Universal stand type 2

This versatile stand allows smooth adjustment of both arm angle and length. This stand is perfect for photomicrography and video mounting of large specimens.

- For use with SZ2-ST5 arm



# Macro View MVX10

## RESEARCH MACRO ZOOM SYSTEM MICROSCOPE



Developed as a system specially designed for macro fluorescence imaging, the MVX10 employs a single-zoom optical system, and has world-class features including a high resolution of 1500 lines/mm, a zoom ratio of 10, and NA of 0.5 (with 2x objective). This enables seamless observation of bright fluorescence images, from macro to micro, and provides extraordinarily high resolution.

|   |                                     |   |              |              |
|---|-------------------------------------|---|--------------|--------------|
| Zoom microscope body<br>MVX-ZB10                | Optical system                      | Mono-zoom variable magnification system   |              |              |
|   | Zoom range                          | 0.63x-6.3x (zoom ratio 1:10)  |              |              |
|   | Aperture iris diaphragm             | Built-in  |              |              |
| Observation head<br>MVX-TTRS                    | Features                            | Tilting binocular head that allows switching between standard and stereo observation  |              |              |
|   | Field number (FN)                   | 22  |              |              |
|   | Tilting angle                       | 0 — 23° continuously variable system  |              |              |
|   | Light path selection                | 2-step binocular 100%/photo 100%  |              |              |
| Reflected light<br>fluorescence unit<br>MVX-RFA | Illumination mode                   | Coaxial reflected light   |              |              |
|   | Filter selection                    | Turret 3 filter + BF  |              |              |
|   | Fluorescence mirror unit            | For CFP, GFP, YFP, RFP separation high quality mirror unit  |              |              |
|   | Light source                        | 100W mercury apo lamp housing and power source, 100W mercury lamp housing and power source, or 75W xenon apo lamp housing and power source          |              |              |
| Magnification changer<br>MVX-CA2X               | Magnification                       | 1x, 2x selection  |              |              |
| Objectives (when used with eyepiece WHN10X)     |                                     | MVPLAPO 0.63X   | MVPLAPO 1X   | MVPLAPO 2XC  |
|   | Total magnification                 | 4.0 — 40x   | 6.3 — 63x    | 12.5 — 125x  |
|   | Working distance (WD)               | 87mm  | 65mm         | 20mm         |
|   | Numerical aperture (NA)             | 0.15  | 0.25         | 0.5          |
|   | Field of view                       | 55 — 5.5mm  | 34.9 — 3.5mm | 17.6 — 1.7mm |
| Stands,<br>transmitted illuminators             | Stands,<br>transmitted illuminators | High-level transmitted light illumination base SZX2-ILLB, Brightfield/darkfield transmitted light illumination base SZX2-ILLD, Large stand SZX2-STL |              |              |
|   | Focusing unit                       | Fine focusing unit SZX2-FOFH, motorized focusing unit SZX-FOA2  |              |              |
|   | Stage                               | Large stage plate, thermoplate, CO <sub>2</sub> incubator   |              |              |

\*Please refer to MVX10 catalog for further details

## DP30BW

HIGH SENSITIVITY  
COOLED CCD CAMERA

Using its Peltier-cooled system, the DP30BW offers quiet, vibration-free operation. Combined with the built-in shutter and new background subtract function (noise at long exposures is reduced by using this), these features enable high-quality recording of even weak fluorescence images.

\* Please refer to DP30BW catalog for further details




DP30BW+BX51

## DP71

DIGITAL CAMERA

By combining Olympus digital technologies together with high-speed processing hardware, even an image of 12.5 million pixels can be captured at high speed, around 3 seconds, while fully maintaining image quality, accuracy and color fidelity. High sensitivity and low noise ensure that even images derived from relatively faint fluorescence can be acquired clearly.

\* Please refer to DP71 catalog for further details



DP71+BX51+PC

## DP20

DIGITAL CAMERA

Image display, storage, and simple measurement can be done from a compact, palm-size handset control unit. The outstanding operability of the functional key layout allows starting up quickly and continuously shoot images at 1 second intervals. The UXGA (1600 x 1200 pixel) compatible monitor enables real-time 15 fps display, and the system is also ideal for high-resolution monitoring applications without a PC.

\* Please refer to DP20 catalog for further details



DP20+BX51

## DP12

DIGITAL CAMERA

The overall design is compact, with a palm-size multi function control unit integrating a 3.5" LCD monitor with 200,000-pixel display, and a small footprint that makes it easy to install and lay out any necessary auxiliary equipment. The 3.34 million-pixel and 1/1.8 inch progressive scanning CCD system ensures highly precise digital images which can be stored at a maximum resolution of 2048x 1536.

\* Please refer to DP12 catalog for further details



DP12+BX51+PC



**UIS2 optics inherit high expandability**


As heir to Olympus' infinity-corrected optical system, in which the tube lens is built into the observation tube, UIS2 optics display no image deterioration even when many different optical components or equipment are inserted in the parallel light path. This inherent expandability gives users ample freedom to construct the system in a way that meets their specific requirements.

**UW (Ultra wideband) multi-coatings reduce autofluorescence and improves S/N ratio**

By using carefully selected raw materials for glass, and applying advanced UW multi-coatings technology, Olympus has reduced objective autofluorescence and significantly improved the S/N ratio.


**Flat, high transmission over wide wavelength range from UV to IR**

UW multi-coatings also yield flat, high transmission over a wide wavelength range, ensuring high performance in research tasks using different types of fluorochromes.



**Complete chromatic aberration compensation up to near infrared region**

UPLSAPO objectives completely eliminate chromatic aberration up to the near infrared region, matching the ability of Super Apochromat objectives to provide clear images without overlapping colors or color shift. As a result, a single objective can perform imaging from UV to IR wavelengths.





## UIS2/UIS Series Objectives

### Universal objectives



#### UPLSAPO series

The top-performance universal Plan Super Apochromat objectives offer an unbeatable solution to every kind of digital imaging need.



#### UPLFLN series

These affordable Semi-Apochromat universal objectives deliver superb resolution, contrast and flatness for any microscopic technique.



#### UPLFL-P series

These strain-free Semi-Apochromat universal objectives reduce internal strain to an absolute minimum and are best suited for polarizing and Nomarski DIC microscopies.

### Brightfield objectives



#### PLAPON series

Designed for unsurpassed resolution and contrast, these Plan Apochromat objectives keep chromatic aberration down to an absolute minimum.



#### PLN series

These cost-effective Achromat objectives ensure field flatness up to F.N. 22 and are widely used in research, educational and routine work applications.

### Objectives for special purpose



#### UPLFLN-PH series

The newly designed phase annuli reduce flare and halo to a minimum and ensure high resolution and contrast for unstained specimens, e.g. living cells and microorganisms.



#### No cover objectives

These no cover objectives are specially designed for microscopy without a cover slip such as for blood smear specimens.



#### UAPO/340 series

These objectives feature a highest transmission of 340nm wavelength light, ensuring maximum performance in fluorescence microscopy through UV excitation including  $CA^{2+}$  photometry.

### Objectives for BX51WI/BX61WI



#### UMPLFLW, LUMPLFLW series

These objectives address the need for high transmission from the near UV to visible light. For ratio imaging (fura-2, 340nm transmission requirement) fluorescence and DIC observation.



#### LUMPLFLW/IR series

These objectives are specially designed for visible band near IR spectral regions. Near IR-DIC imaging deep within thick brain sections can be observed.



#### XLFLUOR/340 series, XLUMPLFL20XW

XLFLUOR/340 series objectives are designed for low magnification fluorescence observation. High N.A. long W.D. XLUMPLFL20XW objective allows the measurement of cell membrane electric potential.

### Objectives for inverted microscopes



#### LUCPLFLN-PH series

These objectives are exclusively designed for culture specimens. An excellent phase-contrast image is assured regardless of the thickness and material of the vessel.



#### LUCPLFLN series

These Semi-Apochromat objectives are dedicated for tissue culture and offer excellent contrast and resolution in brightfield, Nomarski DIC and fluorescence observations.



#### LCACHN series

These Achromat phase-contrast objectives are designed for cell culture observations and are best suited for various clinical examinations and cell testing.

UIS2 objectives \*

\* All UIS2 objectives and WHN eyepieces: lead-free eco-glass

| Description     |                  | N.A.      | W.D. (mm) | F.N. | Cover glass | Immersion | Spring | Correction ring | Iris diaphragm | Water proof & oil proof function | For upright microscope | For inverted microscope |
|-----------------|------------------|-----------|-----------|------|-------------|-----------|--------|-----------------|----------------|----------------------------------|------------------------|-------------------------|
| UPLSAPO         | UPLSAPO 4X       | 0.16      | 13        | 26.5 | —           |           |        |                 |                |                                  | ○                      | ○                       |
|                 | UPLSAPO 10X      | 0.40      | 3.1       | 26.5 | 0.17        |           |        |                 |                |                                  | ○                      | ○                       |
|                 | UPLSAPO 20X      | 0.75      | 0.6       | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      | ○                       |
|                 | UPLSAPO 20XO     | 0.75      | 0.6       | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      | ○                       |
|                 | UPLSAPO 40X      | 0.90      | 0.18      | 26.5 | 0.11-0.23   |           | ○      | ○               |                |                                  | ○                      | ○                       |
|                 | UPLSAPO 60XW     | 1.20      | 0.28      | 26.5 | 0.15-0.2    | Water     | ○      | ○               |                | ○                                | ○                      | ○                       |
|                 | UPLSAPO 60XO     | 1.35      | 0.15      | 26.5 | 0.17        | Oil       | ○      |                 |                | ○                                | ○                      | ○                       |
| UPLSAPO 100XO   | 1.40             | 0.13      | 26.5      | 0.17 | Oil         | ○         |        |                 | ○              | ○                                | ○                      |                         |
| PLAPON          | PLAPON 1.25X     | 0.04      | 5         | 26.5 | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLAPON 2X        | 0.08      | 6.2       | 26.5 | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLAPON 60XO      | 1.42      | 0.15      | 26.5 | 0.17        | Oil       | ○      |                 |                | ○                                | ○                      | ○                       |
|                 | PLAPON 60XOTIRFM | 1.45      | 0.1       | 26.5 | 0.13-0.19   | Oil       | ○      | ○               |                | ○                                | ○                      | ○                       |
| UPLFLN          | UPLFLN 4X        | 0.13      | 17        | 26.5 | —           |           |        |                 |                |                                  | ○                      | ○                       |
|                 | UPLFLN 10X       | 0.30      | 10        | 26.5 | —           |           |        |                 |                |                                  | ○                      | ○                       |
|                 | UPLFLN 20X       | 0.50      | 2.1       | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      | ○                       |
|                 | UPLFLN 40X       | 0.75      | 0.51      | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      | ○                       |
|                 | UPLFLN 40XO      | 1.30      | 0.2       | 26.5 | 0.17        | Oil       | ○      |                 |                | ○                                | ○                      | ○                       |
|                 | UPLFLN 60X       | 0.90      | 0.2       | 26.5 | 0.11-0.23   |           | ○      | ○               |                |                                  | ○                      | ○                       |
|                 | UPLFLN 60XOI     | 1.25-0.65 | 0.12      | 26.5 | 0.17        | Oil       | ○      |                 | ○              | ○                                | ○                      | ○                       |
|                 | UPLFLN 100XO2    | 1.30      | 0.2       | 26.5 | 0.17        | Oil       | ○      |                 |                | ○                                | ○                      | ○                       |
|                 | UPLFLN 100XOI2   | 1.3-0.6   | 0.2       | 26.5 | 0.17        | Oil       | ○      |                 | ○              | ○                                | ○                      | ○                       |
|                 | UPLFLN 4XPH      | 0.13      | 0.17      | 26.5 | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | UPLFLN 10XPH     | 0.30      | 10        | 26.5 | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | UPLFLN 20XPH     | 0.50      | 2.1       | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | UPLFLN 40XPH     | 0.75      | 0.51      | 26.5 | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | UPLFLN 60XOIPH   | 1.25-0.65 | 0.2       | 26.5 | 0.17        | Oil       | ○      |                 | ○              |                                  | ○                      |                         |
| UPLFLN 100XO2PH | 1.30             | 0.2       | 26.5      | 0.17 | Oil         | ○         |        |                 |                | ○                                |                        |                         |
| PLN             | PLN 2X           | 0.06      | 5.8       | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN 4X           | 0.10      | 18.5      | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN 10X          | 0.25      | 10.6      | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN 20X          | 0.40      | 1.2       | 22   | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | PLN 40X          | 0.65      | 0.6       | 22   | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | PLN 50XOI        | 0.90-0.50 | 0.2       | 22   | —           | Oil       | ○      |                 | ○              |                                  | ○                      |                         |
|                 | PLN 100XO        | 1.25      | 0.15      | 22   | —           | Oil       | ○      |                 |                |                                  | ○                      |                         |
| PLN-PH          | PLN 10XPH        | 0.25      | 10.6      | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN 20XPH        | 0.40      | 1.2       | 22   | 0.17        |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN 40XPH        | 0.65      | 0.6       | 22   | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | PLN 100XOPH      | 1.25      | 0.15      | 22   | —           | Oil       | ○      |                 |                |                                  | ○                      |                         |
| PLN & ACHN-P    | PLN 4XP          | 0.10      | 18.5      | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | ACHN 10XP        | 0.25      | 6         | 22   | —           |           |        |                 |                |                                  | ○                      |                         |
|                 | ACHN 20XP        | 0.40      | 3         | 22   | 0.17        |           |        |                 |                |                                  | ○                      |                         |
|                 | ACHN 40XP        | 0.65      | 0.45      | 22   | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
|                 | ACHN 100XOP      | 1.25      | 0.13      | 22   | —           | Oil       | ○      |                 |                |                                  | ○                      |                         |
| PLFLN-CY        | PLFLN10XCY       | 0.3       | 10        | 26.5 |             |           |        |                 |                |                                  | ○                      |                         |
| PLN-CY          | PLN2XCY          | 0.06      | 5.8       | 22   |             |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN4XCY          | 0.1       | 18.5      | 22   |             |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN10XCY         | 0.25      | 10.6      | 22   |             |           |        |                 |                |                                  | ○                      |                         |
|                 | PLN20XCY         | 0.4       | 1.2       | 22   | 0.17        |           | ○      |                 |                |                                  | ○                      |                         |
| LUCPLFLN        | LUCPLFLN 20X     | 0.45      | 6.6-7.8   | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 40X     | 0.60      | 2.7-4     | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 60X     | 0.70      | 1.5-2.2   | 22   | 0.1-1.3     |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 20XPH   | 0.45      | 6.6-7.8   | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 20XRC   | 0.45      | 6.6-7.8   | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 40XPH   | 0.60      | 3.0-4.2   | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 40XRC   | 0.60      | 3.0-4.2   | 22   | 0-2         |           |        | ○               |                |                                  |                        | ○                       |
|                 | LUCPLFLN 60XPH   | 0.70      | 1.5-2.2   | 22   | 0.1-1.3     |           |        | ○               |                |                                  |                        | ○                       |
| UPLFLN-PH       | UPLFLN 4XPH      | 0.13      | 17        | 26.5 | —           |           |        |                 |                |                                  |                        | ○                       |
|                 | UPLFLN 10XPH     | 0.30      | 10        | 26.5 | —           |           |        |                 |                |                                  |                        | ○                       |
| UPLFLN-PHP      | UPLFLN 4XPHP     | 0.13      | 16.4      | 22   | —           |           |        |                 |                |                                  |                        | ○                       |
| CPLFLN          | CPLFLN 10XPH     | 0.30      | 9.5       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | CPLFLN 10XRC     | 0.30      | 9         | 22   | 1.5         |           |        |                 |                |                                  |                        | ○                       |
| LCACHN          | LCACHN 20X PH    | 0.40      | 3.2       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | LCACHN 20X PHP   | 0.40      | 3.2       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | LCACHN 20XRC     | 0.40      | 2.8       | 22   | 1.5         |           |        |                 |                |                                  |                        | ○                       |
|                 | LCACHN 40XPH     | 0.55      | 2.2       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | LCACHN 40XPHP    | 0.55      | 2.2       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | LCACHN 40XPHP    | 0.55      | 1.9       | 22   | 1.5         |           |        |                 |                |                                  |                        | ○                       |
| CACHN & CPLN    | CACHN 10XPHP     | 0.25      | 8.8       | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | CPLN 10XPH       | 0.25      | 10        | 22   | 1           |           |        |                 |                |                                  |                        | ○                       |
|                 | CPLN 10XRC       | 0.25      | 9.7       | 22   | 1.5         |           |        |                 |                |                                  |                        | ○                       |

## UIS objectives

| Description                       |                  | N. A.     | W.D.<br>(mm) | F. N. | Cover<br>glass | Immersion | Spring | Correction<br>ring | Iris<br>diaphragm | Water proof &<br>oil proof cap | For upright<br>microscope | For inverted<br>microscope |
|-----------------------------------|------------------|-----------|--------------|-------|----------------|-----------|--------|--------------------|-------------------|--------------------------------|---------------------------|----------------------------|
| UPLAPO                            | UPLAPO 10XO3     | 0.40      | 0.24         | 26.5  | 0.17           | Oil       | ○      |                    |                   | ○                              | ○                         | ○                          |
|                                   | UPLAPO 10XW3     | 0.40      | 0.43         | 26.5  | 0.17           | Water     | ○      |                    |                   | ○                              | ○                         | ○                          |
|                                   | UPLAPO 40XOI3    | 1.00-0.50 | 0.12         | 26.5  | —              | Oil       | ○      |                    | ○                 | (○)                            | ○                         | ○                          |
| PLAPO                             | PLAPO 40X        | 0.95      | 0.13         | 26.5  | 0.11-0.23      |           | ○      | ○                  |                   |                                | ○                         | ○                          |
| UPLFL-P                           | UPLFL 4XP        | 0.13      | 13           | 26.5  | —              |           |        |                    |                   |                                | ○                         |                            |
|                                   | UPLFL 10XP       | 0.30      | 3.1          | 26.5  | —              |           |        |                    |                   |                                | ○                         |                            |
|                                   | UPLFL 20XP       | 0.50      | 1.6          | 26.5  | 0.17           |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | UPLFL 40XP       | 0.75      | 0.51         | 26.5  | 0.17           |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | UPLFL 100XO3P    | 1.30      | 0.1          | 26.5  | 0.17           | Oil       | ○      |                    |                   |                                | ○                         |                            |
| PLFL                              | PLFL 100X        | 0.95      | 0.2          | 26.5  | 0.14-0.2       |           | ○      | ○                  |                   |                                | ○                         |                            |
| UAPO                              | UAPO 10X/340     |           |              |       |                |           |        |                    |                   |                                |                           |                            |
|                                   | UAPO 20X3/340    | 0.75      | 0.55         | 22    | 0.17           |           | ○      |                    |                   | ○                              |                           | ○                          |
|                                   | UAPO 40X3/340    | 0.90      | 0.2          | 22    | 0.11-0.23      |           | ○      | ○                  |                   | ○                              |                           | ○                          |
|                                   | UAPO 40XOI3/340  | 1.35-0.65 | 0.1          | 22    | 0.17           | Oil       | ○      |                    | ○                 | ○                              |                           | ○                          |
|                                   | UAPO 20XW3/340   | 0.70      | 0.4          | 22    | 0.17           | Water     | ○      |                    |                   | ○                              |                           | ○                          |
|                                   | UAPO 40XW3/340   | 1.15      | 0.25         | 22    | 0.13-0.25      | Water     | ○      | ○                  |                   | ○                              |                           | ○                          |
| APO                               | APO 100XOHR      | 1.65      | 0.1          | 22    | 0.15           | Oil       | ○      |                    |                   | ○                              |                           | ○                          |
| Low magnification<br>fluorescence | XLFLUOR 2X/340   | 0.14      | 21**         | 22    | 0-5 (Water)    |           |        |                    |                   |                                |                           |                            |
|                                   | XLFLUOR 4X/340   | 0.28      | 29.5**       | 22    | 0-5 (Water)    |           |        |                    |                   |                                |                           |                            |
| Super high N.A.                   | XLUMPLFL 20XW    | 0.95      | 2            | 22    |                | Water     |        |                    |                   |                                |                           |                            |
| No cover<br>objective             | MPLAPO 50X       | 0.95      | 0.3          |       | 0              |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | MPLAPO 60X       | 0.90      | 0.4          |       | 0              |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | MPLAPO 100XO     | 1.40      | 0.1          |       | 0              | Oil       | ○      |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 40X       | 0.75      | 0.63         |       | 0              |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 50X       | 0.80      | 0.66         |       | 0              |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 100X      | 0.95      | 0.31         |       | 0              |           | ○      |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 10XW      | 0.30      | 3.3          | 26.5  | —              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 20XW      | 0.50      | 3.3          | 26.5  | —              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 40XW      | 0.80      | 3.3          | 26.5  | 0              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | UMPLFL 60XW      | 0.90      | 2            | 26.5  | 0              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | LUMPLFL 40XW/IR2 | 0.80      | 3.3          | 26.5  | 0              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | LUMPLFL 60W/IR2  | 0.90      | 2            | 26.5  | 0              | Water     |        |                    |                   |                                | ○                         |                            |
|                                   | LUMPLFL 100XW    | 1.00      | 1.5          | 26.5  | 0              | Water     |        |                    |                   |                                | ○                         |                            |
| LSM objective                     | PLAPO 40XWLSM    | 0.90      | 0.16         | 22    | 0.17           | Water     | ○      |                    |                   | ○                              | ○                         | ○                          |
|                                   | PLAPO 60XWLSM    | 1.00      | 0.15         | 22    | 0.17           | Water     | ○      |                    |                   | ○                              | ○                         | ○                          |
|                                   | PLAPO 40XOLSM    | 1.10      | 0.13         | 22    | 0.17           | Oil       | ○      |                    |                   | ○                              | ○                         | ○                          |

\*\* Include 5mm water (○): oil proof cap applicable

## Macro Imaging



## Time Lapse



## Live-cell Confocal Microscopy

## Patch Clamping

Image data courtesy of:

Pr. Kentaro Arikawa, Ms. Mituyo Kinoshita,  
Laboratory of Neuroethology, Graduate School of Integrated Science, Yokohama City University (cover page left top)

Atsushi Miyawaki M.D., Ph.D, Ms. Ryoko Ando,  
RIKEN Brain Science Institute Laboratory for Cell Function Dynamics (cover page above middle)

Mark Ellisman Ph.D, Hiroyuki Hakozaki M.S., Natalie Maclean M.S.,  
University of California, San Diego, NCMIR (back page left top)

Tohru Murakami M.D., Ph.D,  
Department of Neuromuscular and Developmental Anatomy, Gunma University Graduate School of Medicine (back page left bottom)



Products with Olympus Eco-products mark in this catalog are environmentally conscious products made according to OLYMPUS standard.  
This does not apply to some accessories.  
Please visit our website for more information:

**ECO-PRODUCTS** <http://www.olympus.co.jp/en/eco-products>



ISO9001 Certification  
OLYMPUS CORPORATION Micro-Imaging System Division



ISO14001 Certification  
OLYMPUS CORPORATION Ina-Plant



ISO9001/ISO14001 Certification  
Olympus Optical Technology Philippines Inc.

**OLYMPUS CORPORATION has obtained ISO9001/ISO14001**

Specifications are subject to change without any obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo, Japan  
**OLYMPUS LIFE AND MATERIAL SCIENCE EUROPA GMBH**  
Postfach 10 49 08, 20034, Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
3500 Corporate Parkway, Center Valley, Pennsylvania 18034-0610, U.S.A.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS AUSTRALIA PTY. LTD.**  
31 Gilby Road, Mt. Waverley, VIC 3149, Melbourne, Australia.

**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.  
**OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.**  
12-13F, NCI Tower, A12 Jianguomenwai Avenue, Chaoyang District,  
Beijing, 100022, China

[www.olympus.com](http://www.olympus.com)



This catalog is printed by environmentally-friendly waterless printing system with soy ink.

Printed in Japan M05E-0906B

**OLYMPUS**<sup>®</sup>

Your Vision, Our Future

SYSTEM MICROSCOPE

**CX41**

CX2 SERIES


**UIS2**  
World-leading optics

*Optics and performance outstanding in its class*



# Advanced optical and system performance with excellent cost-efficiency

The evolution of this bestselling microscope delivers new, cost-efficient improvements in both optics and system performance. Newly equipped with Olympus' leading-edge UIS2 optical system, it provides significantly enhanced image clarity in a variety of observation methods, from brightfield to reflected light fluorescence. For all inspection and training applications in the fields of biology and medicine, the CX41 sets the pace for its class in both basic and system performance.






## Outstanding flat images from PLCN objectives

The CX41 provides images of outstanding brightness and clarity in a variety of observation modes. As well as Olympus' renowned UIS2 optics infinity system, it employs the PLCN series of Plan Achromat objectives, which are made from carefully selected top quality glass and manufactured with the most rigorous precision. The result is a major improvement in image flatness, with the 10X and 40X objectives in particular providing images that are among the very best in this class of microscope. Transmitted light illumination is from a 6V, 30W high-intensity halogen light source.

### Flatness area comparison



# Excellent cost-performance in reflected light fluorescence and other observation methods





**Slide condenser / CX-SLC  
Brightfield condenser / CH3-CD**  
These Abbe type condensers allow brightfield observations from 4X to 100X. Accurate centering is provided by the attachment lens (CX-AL) and the iris diaphragm, to exclude unnecessary light and obtain bright Koehler illumination right across the magnification range. These highly economical condensers enable phase contrast and darkfield observations by simply adding basic accessories.



**Simple phase contrast attachment / CX-PH1, 2, 3**  
For phase contrast observations at 10X, 40X and 100X.




**Darkfield central stop / CH2-DS**  
For darkfield observations from 4X to 40X.  
\* Separate filter holder (CH2-FH) or attachment lens (CX-AL) required.

**Low magnification adapter /CX-LA**  
Use of the 2X low magnification objective allows macro observation.

## Reflected light fluorescence attachment / CX-RFA-2

Users can choose between blue or green excitation and transmitted light observations. UIS2 optics provide bright fluorescence images, with no intermediate magnifications when changing from transmitted light to fluorescence observation. Standard PLCN objectives can be used without replacement.



## Simple polarizing condenser / CH3-CDP

With the optional plate adapter U-TAD, polarizing observations from 4X to 100X using a tint plate can be performed. A U-GAN analyzer is provided for gout inspection. Polarizing objectives from 4X to 100X are available.

\* Separate polarizer U-POT and analyzer U-ANT required.



## Dry darkfield condenser / CX-DCD

This dry-type darkfield condenser gives a superior darkfield effect without the need for immersion in oil. Suitable for use at 10X and 40X magnifications.





# Dependable basic performance assures outstanding operational convenience



**Move the specimen with just one finger**  
Rubber grips are provided for the stage handles, allowing the specimen to be moved smoothly with just one finger. The slim body and conveniently positioned controls ensure that everything is within easy reach, so operators can maintain a natural posture.



**Tilting binocular tube**  
The tilting binocular tube lets each operator select the most suitable and comfortable eyepoint — a valuable contribution to reducing fatigue in extended observation sessions. High cost-efficiency type (U-CTBI) is also available.



U-TBI-3



U-CTBI

**Inward-facing quintuple revolving nosepiece**  
The quintuple revolving nosepiece enables a wide range of magnification observations and various combinations of objectives. Unrestricted access to the space in front of the objective allows specimens to be exchanged quickly and easily.



**Rackless stage with enhanced operability**  
To keep the work area clear, and to avoid interference with observation operations, the X-direction travel guide does not extend out from the side of the stage. The main and sub-scale displays are designed for easy read-out.



## Phase-contrast condenser / CX-PCD

The multi-purpose CX-PCD condenser allows observation of brightfield, phase-contrast and darkfield images without exchanging condensers. Phase-contrast observation from 10x to 100x and darkfield observation from 10x to 40x is allowed.



**Torque adjustable focusing knob**  
The torque of the coarse focusing knob can be adjusted, to suit different operators' needs and to make focusing smooth and easy while keeping the hands on the desk. A stage upper limit stopper is also provided.



**Easy transportation and installation**  
The CX41 is eminently portable, with convenient handgrips at the front and back of the frame and no inconvenient protrusion of the stage guide.



## Anti-fungus treatment

The treatment applied to the observation tubes, eyepieces and objectives, protects quality of optical parts even in high humidity regions.

# More accessories, more observation versatility



**Dual observation attachment/ U-DO3**  
Enables dual, simultaneous observation of a single specimen from the same direction with equal magnification and brightness for both operators. A pointer can be used to indicate specific sections of the specimen to simplify the training process and enhance discussion.

**Trinocular tube/ U-CTR30-2**  
Digital (DP20) or video cameras can be attached for on-site/ remote illustration of educational/ discussion meetings.



CX41+DP20

## Eyepoint adjuster/ U-EPA2

Allows the eyepoint position to be raised by 30mm. Up to two eyepoint adjusters can be attached between the frame arm and observation tube.



## Arrow pointer/ U-APT

Enables insertion of an LED arrow for display in a digital image.



## 2x magnification changer/ U-ECA

Magnification is doubled by engaging the auxiliary 2x lens.



## Drawing attachment/ U-DA

The drawing attachment projects an image of the pencil and drawing surface into the visual field.




## Trinocular intermediate attachment/ U-TRU

A binocular tube on its own allows digital imaging in combination with this trinocular intermediate attachment.



# CX41 SYSTEM DIAGRAM




\*1 Please consult your Olympus dealer for detail. \*2 U-TV1x cannot be attached for technical reason. Combine U-TV1x-2 for 1x adapter. \*3 10x eyepieces incorporated. F.N. 18

## Specifications

|                     |  |                                     |                                      |                           |
|---------------------|--|-------------------------------------|--------------------------------------|---------------------------|
| Item                | CX41   |                                     |                                      |                           |
| Optical system      | UIS2 (Universal Infinity-corrected) optical system   |                                     |                                      |                           |
| Illumination        | Built-in transmitted Koehler illuminator<br>6V30W halogen bulb<br>100-120V/220-240V~ 0.85/0.45A 50/60Hz  |                                     |                                      |                           |
| Focusing            | <ul style="list-style-type: none"> <li>• Stage height movement by roller guide (rack &amp; pinion)</li> <li>• Stroke per rotation: 36.8mm</li> <li>• Full stroke range: 25mm</li> <li>• Upper limit stopper</li> <li>• Tension adjustment on coarse focus adjustment knob</li> </ul> |                                     |                                      |                           |
| Revolving nosepiece | Fixed quintuple nosepiece with inward tilt   |                                     |                                      |                           |
| Observation tube    | Type   | U-CBI30-2, Binocular                | U-CTR30-2, Trinocular                | U-CTBI, Tilting binocular |
|                     | Field number   | 20                                  | 20                                   | 18                        |
|                     | Tube inclination   | 30°                                 | 30°                                  | 30°—60°                   |
|                     | Interpupillary distance adjustment range   | 48—75mm                             | 48—75mm                              | 48—75mm                   |
|                     | Light path selector  | None                                | None (Bi 50%, Video/Photo 50% fixed) | None                      |
| Stage               | Size   | 188(W)X134(D)mm                     |                                      |                           |
|                     | Movement range   | 76mm X-direction X 50mm Y-direction |                                      |                           |
|                     | Specimen holder  | Double slide holder                 |                                      |                           |
|                     | Rubber grip  | Equipped as standard                |                                      |                           |
| Condenser (CH3-CD)  | Type   | Abbe condenser                      |                                      |                           |
|                     | N.A.   | 1.25 with oil immersion             |                                      |                           |
|                     | Aperture iris diaphragm  | Built-in                            |                                      |                           |
| Dimensions & weight | 233(W) X 432(H) X 367.5(D)mm, approximately 8.5kg (approximately 18.7 lb.)   |                                     |                                      |                           |

## Dimensions

(Unit: mm)



## Objectives, Plan Achromat

| PLCN   | Numerical Aperture (N.A.) | Working Distance (W.D.) |
|--------|---------------------------|-------------------------|
| 4X     | 0.10                      | 18.5mm                  |
| 10X    | 0.25                      | 10.6mm                  |
| 20X    | 0.4                       | 1.2mm                   |
| 40X    | 0.65                      | 0.6mm                   |
| 60X    | 0.8                       | 0.2mm                   |
| 100XO  | 1.25                      | 0.13mm                  |
| 100XOI | 1.25-0.6                  | 0.13mm                  |

## Eyepiece

|                      | Field Number (F.N.) |
|----------------------|---------------------|
| WHB10X               | 20                  |
| CWH10X (for U-TBI-3) | 20                  |

CX41 is the environmental conscious product according to OLYMPUS's own standards.

Main features of OLYMPUS Eco-products are as follows.

- Lead-free and arsenic-free Eco-glass for optics, such as lenses and prisms.
- Exclusion of hexavalent chrome, mercury, lead and cadmium from metal materials and surface treatment of metal.
- Exclusion of lead solders.
- Adoption of cardboard for packing materials without styrene foam for promoting the recycling.
- \* A definition of exclusion depends on olympus standard. Some accessories are inapplicable.

Please visit our web site for further information:  
<http://www.olympus.co.jp/en/eco-products/>



OLYMPUS CORPORATION obtains ISO9001/ISO14001

ECO-PRODUCTS

Specifications are subject to change without any obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo, Japan  
**OLYMPUS LIFE AND MATERIAL SCIENCE EUROPA GMBH**  
Postfach 10 49 08, 20034, Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
3500 Corporate Parkway, Center Valley, Pennsylvania 18034-0610, U.S.A.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS AUSTRALIA PTY. LTD.**  
31 Gilby Road, Mt. Waverley, VIC 3149, Melbourne, Australia.

**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.  
**OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.**  
12-13F, NCI Tower, A12 Jianguomenwai Avenue, Chaoyang District, Beijing, 100022, China



This catalog is printed by environmentally-friendly waterless printing system with soy ink.

www.olympus.com

Printed in Japan M183E-0706B

# Specifications

| Item                                  | Specifications  | CH20i<br>Binocular/Trinocular  |
|---------------------------------------|---|--|
| <b>CH20i<br/>MICROSCOPE<br/>FRAME</b> | <ul style="list-style-type: none"> <li>Coaxial coarse/fine knobs: Tension adjustment on the right side</li> <li>Fine focus knob graduated rack and pinion</li> <li>Stage movement (XY direction) on fixed) Plane stage 120 x 132 mm</li> <li>Quadrant revolving nosepiece</li> <li>With right hand mechanical stage</li> <li>Abbe condenser N.A. 1.25 (oil immersion), with aperture iris diaphragm</li> <li>Blue filter</li> <li>Universal Power Supply (100V to 240V) for 6V 20W illuminator</li> <li>8cc immersion oil</li> <li>Dust cover</li> <li>Mirror unit (Plano-concave)</li> </ul> | <ul style="list-style-type: none"> <li>Binocular observation tube (inclination 45°, interpupillary distance, adjustment range 53-75 mm), diopter adjustment on the left</li> </ul> |
|                                       | <ul style="list-style-type: none"> <li>Trinocular observation tube (inclination 45°), diopter adjustment on the left</li> </ul>   | ○  |
| <b>Power Cord</b>                     |   | ○  |
| <b>Lamp</b>                           | 6V 20W halogen lamp ( x 2)  | ○  |
| <b>Objectives</b>                     | /NEA-Achromat 4X (anti-fungus)  | ○  |
|                                       | /NEA-Achromat 10X (anti-fungus)   | ○  |
|                                       | /NEA-Achromat 40X (anti-fungus) spring  | ○  |
|                                       | /NEA-Achromat 100X (anti-fungus) spring, oil  | ○  |
| <b>Eyepiece</b>                       | /CWHK10X (LB eyepiece 10X), F.N. 18mm, (anti-fungus) ( x 2)   | ○  |

Specifications are subject to change without any obligation on the part of the manufacturer

# OLYMPUS

BIOLOGICAL MICROSCOPE

## CH20i

In Japan, OLYMPUS prides itself as a company that manufactures quality optical products with a FOCUS ON LIFE. Olympus (India) extends that corporate philosophy as it brings the renowned Japanese technical expertise to its making of microscopes. It is progress through precision.



©2009 is manufactured under license from Olympus Corporation, Japan

## OLYMPUS

OLYMPUS (INDIA) PVT. LTD.


A-5, Michan Co-operative Industrial Estate, Mathura Road, New Delhi-110044  
Tel.: 30686766, 30686741, 30686744 Fax: 91-11-30686737 Email: oisales@idsimgate.com



ISO 9001:2000


Quality that's pure

**OLYMPUS**



CH20i

## Features



Sealed optical parts

**Anti Fungus Optics**  
From the eyepiece and objectives to the interior of the observation tube, comprehensive anti fungus treatments are applied to every area which effects the clarity of the observed images. This Tropicalized treatment ensures image excellence for long periods in conditions favouring to fungus growth.




Main operating controls within easy reach

**Ergonomic Design**  
All the main operating functions are within easy reach, thereby allowing extended observation without fatigue of eye or posture.

### Window in the arm section

The specimen can be observed with the naked eye from the back. Also, convenient to observe the specimen while using microscope through natural light by using reflector.



Optimization of illumination with aspheric lenses

### High Performance Aspheric Lenses

The abbe condenser & the light relay system are fitted with high performance aspheric lenses, which ensure uniform illumination, resulting in a bright and crisp image.

### Illumination System

Built-in illumination through SMPS circuit for constant voltage output to cover input voltage fluctuations from 100 volts to 240 volts. This prolongs bulb life and prevent frequent bulb change.



Easily interchangeable objectives

### Interchangeability of Objectives

The high precision objectives can be mounted on the nosepiece in any sequence without effecting their parfocality or centering.

## Technology From Japan


The CH-20i is an accumulation of advanced Japanese manufacturing technology that turns creative designs into quality products. Now, made in India under licence from Olympus Corporation, Japan, the CH-20i is the manifestation of the uncompromising quality standards and user friendly features which OLYMPUS products are known for.

The manufacturing process, under the strict supervision of a team of Japanese engineers, ensures compliance with demanding performance standards. OLYMPUS technology is constantly applied to every aspect of design, maximising the interchangeability of parts and minimising adjustment requirements during assembly. Consequently, the CH-20i is the only microscope in India conforming to internationally accepted standards.


CH20i

Quality that's pure

**OLYMPUS**



CH20i



4-axis Machining Centre



Clean Room Assembly area



Laser Interferometer for optics inspection



Lens Inspection

## Optional Accessories



**Digital Imaging Compliance**  
An optional adapter to mount digital cameras is provided to allow easy, cost-efficient digital imaging.




**Darkfield central stop/  
CH2-DS**  
Darkfield observation from 10X-40X is possible. This is placed underneath the condenser with attachment lens CHS-AL or filter holder CH2-FH.

## Attachments

**Phase-contrast Attachment/10X,40X**  
This is attached underneath the condenser. Phase-contrast objectives (10X & 40X) are provided.

**CCTV Attachment**  
Supplied with optical C-mount Video Adapter  
**Micro Image Projection System**



## Specifications

| Item                                  | Specifications  | CH20/<br>Binocular | Trinocular |
|---------------------------------------|---|--------------------|------------|
| <b>CH20/<br/>MICROSCOPE<br/>FRAME</b> | <ul style="list-style-type: none"> <li>Coaxial coarse/fine knobs: Tension adjustment on the right side</li> <li>Fine focus knob graduated</li> <li>Stage movement (XY direction) on rack and pinion</li> <li>Quadruple revolving nosepiece (fixed) Plane stage 120 x 132 mm</li> <li>With right hand mechanical stage</li> <li>Abbe condenser N.A. 1.25 (oil immersion), with aperture iris diaphragm</li> <li>Blue filter</li> <li>Universal Power Supply (100V to 240V) for 6V 20W illuminator</li> <li>8cc immersion oil</li> <li>Dust cover</li> <li>Mirror unit (Plano-concave)</li> </ul> | ○                  | ○          |
| <b>Power Cord</b>                     |   | ○                  | ○          |
| <b>Lamp</b>                           | 6V 20W halogen lamp (x 2)   | ○                  | ○          |
| <b>Objectives</b>                     | /NEA Achromat 4X (anti-fungus)  | ○                  | ○          |
|                                       | /NEA Achromat 10X (anti-fungus)   | ○                  | ○          |
|                                       | /NEA Achromat 40X (anti-fungus) spring  | ○                  | ○          |
| <b>Eyeiece</b>                        | /NEA Achromat 100X (anti-fungus) spring, oil  | ○                  | ○          |
|                                       | /CWHK10X (LB eyepiece 10X), F.N. 18mm, (anti-fungus) (x 2)  | ○                  | ○          |

Specifications are subject to change without any obligation on the part of the manufacturer

OLYMPUS is manufactured under license from Olympus Corporation, Japan

# OLYMPUS

OLYMPUS (INDIA) PVT. LTD.  
A-5, Mohan Co-operative Industrial Estate, Mathura Road, New Delhi-110044  
Tel.: 30886766, 30886741, 30886744 Fax: 91-11-30886737 Email: olsales@dssimage.com



ISO 9001:2000



**OLYMPUS**<sup>®</sup>


Your Vision, Our Future

Biological Microscope

**CX21i**

CX2 Series

**UIS2**  
World-leading optics



**Step up to  
Higher Performance  
with the Olympus  
UIS2 Infinity Optics**


The CX21i sets new standards for  
educational and laboratory applications

# Superb image clarity coupled with Oly

## Superior imaging

### Plan objectives are packaged as standard, providing image flatness that is among the best in this class

The CX21i is equipped with the same UIS2 optical system used in Olympus' top-level microscopes. In addition, Plan Achromat objectives are included as standard for the first time in this class of microscope, providing sharp images with high contrast right up to the edge of the field of view.




### Bright, uniform illumination

The CX21i's 6V20W high intensity halogen lamp delivers clear, stable illumination. The built-in aspheric collector lens provides images that are bright and uniformly illuminated over the entire field of view.

### Abbe condenser for optimized contrast


The optimal aperture stop position for each objective is clearly marked on the Abbe condenser, resulting in high-resolution, high-contrast images.



## Easy, safe operation


### Rackless stage for durability and ease of use

The wire-driven stage movement avoids any protrusion of racks — preventing the possibility of accidental damage or hand injury. Abrasion-resistant rackless design ensures continuous smooth movement.



### Component security — Student proof!

The eyepieces, objectives and condenser are all factory attached to the microscope body ensuring that no component is dropped or detached during transportation or regular use.



### Ultra-smooth, quadruple revolving nosepiece

High-precision machining ensures smooth operation and durability. Nosepiece includes a rubber ring for an easy ergonomic grip.

### Mechanical stage focus-lock prevents mishaps

The focus position can be locked making it easy to refocus when the specimen is changed. With the upper stage position locked, there is no chance of the objective accidentally touching the slide and damaging either the specimen or the objective itself.




# Olympus' renowned reliability and rigidity

## Use in any environment

### Anti-fungus treatment for component durability

An effective anti-fungus treatment is applied to the objectives, eyepieces and observation tube for consistent image clarity and long operating life even in hot humid work conditions.




### Compact design, easy to store and carry

The CX21i is only 391 mm high, compact enough to be stored in a typical classroom cabinet.



### Individual choice of adjustment

Binocular observation tube is inclined 30 degrees allowing for an ergonomic, comfortable posture during observation. Eyepieces can be easily set for multiple users with a scaled interpupillary distance adjustment and diopter setting. High eyepoint design eyepieces with F.N. 20 can be used comfortably without removing eyeglasses.



## Suitable for any user

### Tension adjustment for smoother focusing

The coaxial coarse/fine focusing knob is operable from either the left or right side. Coarse movement tension can be adjusted in accordance with personal preference.

### Ergonomic design for user convenience

The focusing knob, light control adjustment knob and stage handle are located close together, the user can work with hands on the desk and operate the CX21i with minimal movement.



## Specifications: CX21i Biological Microscope


|                      |  |
|----------------------|--|
| Microscope Frame     | CX21 FS1   |
| Optical system       | UIS2 (Universal Infinity System) optical system  |
| Illumination System  | Built-in transmitted illumination system<br>6 V 20 W halogen bulb<br>100-240 V 50/60 Hz universal voltage  |
| Focusing             | Stage height movement (coarse movement stroke 20 mm)<br>Fine focus graduation: 2.5 $\mu$ m   |
| Revolving Nosepiece  | Fixed quadruple nosepiece  |
| Stage                | Wire movement mechanical fixed stage: 120x132 mm<br>Traveling range: 76 mm (X) x 30 mm (Y)<br>Single specimen holder   |
| Observation Tube     | 30° inclined binocular tube<br>Interpupillary distance adjustment range 48–75 mm   |
| Condenser            | Abbe type with aperture iris diaphragm N.A.: 1.25  |
| Objective            | Plan Achromatic objectives (anti-fungus)<br>4x N.A.: 0.10 W.D.: 18.5 mm<br>10x N.A.: 0.25 W.D.: 10.6 mm<br>40x N.A.: 0.65 W.D.: 0.6 mm<br>100x N.A.: 1.25 W.D.: 0.13 mm (option)             |
| Eyepiece (10x)       | Field Number (F.N.): 20 (anti-fungus)  |
| Optional Accessories | Mirror unit, 15x eyepiece (F.N. 12, anti-fungus), micrometer, filter holder, darkfield stop, Micro image projection system, Phase-contrast Attachment/10X, 40X, Simple Polarizing attachment |

## Features:

- Ergonomic & compact design for user convenience
- UIS2 Plan infinity optics providing image flatness
- Seidentopf observation head for individual choice of adjustment
- Anti fungus treatment for component durability
- High eyepoint design eyepiece (F.N.20)
- High performance aspheric lenses in the abbe condenser & light relay system for bright & uniform illumination
- SMPS circuit for constant voltage output
- Rackless stage for durability and ease of use
- Mechanical stage focus-lock prevents mishaps
- Component security – Student proof

## Dimensions

(unit: mm)



- Over 200,000 highly satisfied CX21 customers world over
- CX21i is now exclusively for India at affordable cost

CX21i is the environmental conscious product according to OLYMPUS' own standards.

- Adoption of cardboard for packing materials without styrene foam for promoting the recycling.
- Lead-free and arsenic-free Eco-glass for optics, such as lenses and prisms.

• The CX21i meets CE standards for safety.

- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

CX21i is manufactured under license from Olympus Corporation, Japan

# OLYMPUS®

## OLYMPUS (INDIA) PVT. LTD.

Works : A-3, Sec-81, Phase-II, Noida-201305 U.P (India)

Regd Office:

A-5, Mohan Co-operative Indl. Estate, (Opp. Sarita Vihar), Mathura Road, New Delhi-110044

Ph.: 011 3088 6742 / 43 / 44 / 66, Fax: 011-2695 9382, email: oisales@dssimaging.com



**OLYMPUS**<sup>®</sup>

Your Vision, Our Future


BIOLOGICAL MICROSCOPE

**CX31**

CX2 SERIES

**UIS2**  
World-leading optics

*New standards of sharpness, clarity and flatness*




# Advanced optical performance with

The Olympus CX2 microscopes, which have gained an outstanding worldwide reputation in many medical and educational arenas, now evolve with new UIS2 infinity optics. The CX31 microscopes improve all-round performance and offer excellent cost-efficiency.

## Outstanding flat images from PLCN objectives

The PLCN series of UIS2 objectives improves flatness dramatically, producing sharp, clear images right up to the edge of the field of view. Ideal for the 10X and 40X objectives so frequently used for inspection work, the flatness ranks among the very best in this class of microscopes.



Flatness area comparison

## Bright, even observation images

The illumination system employs a high intensity 6V, 30W halogen lamp for bright images. The aperture iris diaphragm with built-in condenser and standard field stop combine to provide bright and even illumination at all levels of magnification.

## Rigid, durable construction for high performance and long service life

Construction quality is excellent throughout, with objectives, eyepiece, observation tube, revolving nosepiece, highly reliable rackless stage and other components all fixed firmly to the body — so there's nothing to come loose or fall off when the microscope is being transported. The highly reputable rackless stage is employed, and since the X-axis guide does not protrude, both transportation and operation are performed easily and safely.

## For digital imaging (optional)

An optional adapter (sold separately) to mount Olympus digital cameras is provided to allow easy, cost-efficient digital imaging.



## Anti-fungus treatment

The treatment applied to the observation tubes, eyepieces and objectives, protects quality of optical parts even in high humidity regions.



# h excellent cost-efficiency.

## UIS2 eyepiece

Provides wide field of view (F.N. 20) and allows easy observation with eyeglasses.

## Observation tube (Binocular/Tilting binocular)

Diopter adjustment is performed via the knurled collar on the left eyetube. Wide interpupillary distance adjustment (48-75mm) ensures the best and most comfortable observation conditions for every user. Tilting binocular is also available for your comfortable posture in operation.

## Quadruple inward-facing revolving nosepiece

Complies with a wide range of magnifications, allows unrestricted use of the space in front of the objectives and makes it easy to confirm observation magnifications.

## Coarse upper limit stopper

Locks the upper position of the stage, preventing contact between the objective and the specimen in high magnification observations.

## UIS2 objectives

4X, 10X, 40X and 100X PLCN objectives, providing world-class image flatness.

## Rackless stage

Rack-free stage with no side protrusions. Low-positioned control knob allows smooth and comfortable specimen movement. Scale gradations are in white lettering on a black background.

## Light intensity adjustment

Continuous light adjustment is possible.

## Abbe condenser

An Abbe condenser, with N.A. 1.25 and built-in aperture diaphragm, provides the appropriate diaphragm setting to suit different specimens and magnifications.

## Field stop

Frame integrated, with ability to accept a  $\varnothing 45\text{mm}$  filter.

## Coaxial coarse/fine focusing knob

The coaxial coarse/fine focusing knob allows each operator to adjust the torque for coarse focusing operations. Focusing is smooth and easy, with the user's hands placed on the desk.

## Stage handle with tactile grips


Tactile grips are fitted to the X and Y stage controls to provide a "light touch" stage movement.

## Highly rigid frame

Extra rigid frame withstands frequent use and repositioning.

## Hand grips for easy portability

Convenient hand grips at the front and back of the frame make it easy to carry the microscope.




## Specifications

|                     |   |   |                   |
|---------------------|---|---|-------------------|
| Item                | CX31  |   |                   |
| Optical system      | UIS2 (Universal Infinity-corrected) optical system  |   |                   |
| Illumination        | Built-in transmitted Koehler illuminator<br>6V30W halogen bulb<br>100-120V/220-240V~ 0.85/0.45A 50/60Hz   |   |                   |
| Focusing            | <ul style="list-style-type: none"> <li>•Stage height movement by roller guide (rack &amp; pinion)</li> <li>•Stroke per rotation: 36.8mm</li> <li>•Full stroke range: 25mm</li> <li>•Upper limit stopped by simplified pre-focusing dial</li> <li>•Tension adjustment on coarse focus adjustment knob</li> </ul> |   |                   |
| Revolving nosepiece | Fixed quadruple nosepiece with inward tilt  |   |                   |
| Observation tube    | Type  | Binocular                                     | Tilting binocular |
|                     | Field number  | 20  | 18                |
|                     | Tube inclination  | 30°   | 30°—60°           |
|                     | Interpupillary distance adjustment range  | 48—75mm                                       | 48—75mm           |
| Stage               | Size  | 188(W) X 134(D)mm                             |                   |
|                     | Movement range  | 76mm X-direction X 50mm Y-direction           |                   |
|                     | Specimen holder   | Double slide holder                           |                   |
|                     | Rubber grip   | Equipped as standard                          |                   |
| Condenser           | Type  | Abbe condenser, with built-in daylight filter |                   |
|                     | N.A.  | 1.25 with oil immersion                       |                   |
|                     | Aperture iris diaphragm   | Built-in                                      |                   |
| Dimensions & weight | 233(W) X 411(H) X 367.5(D)mm, approximately 8kg (approximately 17.6 lb.)  |   |                   |

## Dimensions

(Unit: mm)



## Objectives, Plan Achromat

| PLCN  | Numerical Aperture (N.A.) | Working Distance (W.D.) |
|-------|---------------------------|-------------------------|
| 4X    | 0.10                      | 18.5mm                  |
| 10X   | 0.25                      | 10.6mm                  |
| 40X   | 0.65                      | 0.6mm                   |
| 100XO | 1.25                      | 0.13mm                  |

## Eyepiece

|                         | Field Number (F.N.) |
|-------------------------|---------------------|
| 10X (Binocular)         | 20                  |
| 10X (Tilting binocular) | 18                  |




ECO-PRODUCTS

### Environmentally-friendly product

All optical components in CX31 microscope use lead-free eco-glass.

Specifications are subject to change without any obligation on the part of the manufacturer.



**OLYMPUS**

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome, Shinjuku-ku, Tokyo, Japan  
**OLYMPUS EUROPA GMBH**  
Postfach 10 49 08, 20034, Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
Two Corporate Center Drive, Melville, NY 11747-3157, U.S.A.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS UK LTD.**  
2-8 Honduras Street, London EC1Y 0TX, United Kingdom.

**OLYMPUS AUSTRALIA PTY. LTD.**  
31 Gilby Road, Mt. Waverley, VIC 3149, Melbourne, Australia.  
**OLYMPUS LATIN AMERICA, INC.**  
6100 Blue Lagoon Drive, Suite 390 Miami, FL 33126-2087, U.S.A.



This catalog is printed by environmentally-friendly waterless printing system with soy ink.

www.olympus.com

Printed in Japan M1546E-0105B



# OLYMPUS



Since the past 44 years Olympus (India) Pvt. Ltd. - in collaboration with Olympus Corporation, Tokyo, Japan has pursued the path of excellence in the field of optical engineering. To maintain international standards set for us by our Japanese collaborators, Olympus India has developed a completely integrated factory wherein optical lens processing, vacuum coating, finishing and other processes are carried out under the strict vigil of Japanese trained engineering personnel. This self-contained infrastructure alongwith a rigid quality control system ensures the production of high resolution optics and trouble-free instruments complying to the high OLYMPUS standards.

The OLYMPUS range of optical instruments comprise of:

- Biological Microscopes
- Research Microscopes

# OLYMPUS



## HSA

STUDENT MICROSCOPE

The HSA is based on a compact design and is extremely reliable in mechanical and optical performance. A wide variety of standard and optional accessories enhance the characteristics features of the HSA. The sphere of its application have been sufficiently extended to meet various requirements at high schools, colleges and universities.

|                      |  |
|----------------------|--|
| <b>Body</b>          | : Monocular die-cast body inclinable upto a horizontal position (90°)                            |
| <b>Magnification</b> | : 100 X – 600 X  |
| <b>Stage</b>         | : Plain square stage of 100 X 100 mm, with two clips, mechanical stage (Optional)                |
| <b>Nosepiece</b>     | : Triple hole nosepiece with positive click stops for precise and smooth alignment of objectives |
| <b>Focusing</b>      | : By coarse and fine focusing knobs  |
| <b>Eyeieces</b>      | : Huygenian 10X & 15X  |
| <b>Objectives</b>    | : Hard anti-reflection coated, colour coded Parfocal Achromat: 10X & 40X                         |
| <b>Condenser</b>     | : Built-in Sub-stage condenser with iris diaphragm   |
| <b>Illumination</b>  | : Plano-concave mirror in adjustable fork mount  |
| <b>Finish</b>        | : Olympus textured grey  |

Magnus

invi

Long Working Distance (LWD), Plan Infinity Optics

4x, 40x Objectives & Phase 10x, 20x as Standard

Pre-centred Phase Annulus

Trinocular Port in  
Standard Unit

FOV 22mm with 10x Eyepiece

Optional Items :

Phase 40x Objective

Attachable Mechanical Stage

Warm Plate

Terasaki, Petri Dish

& Slide Glass Holders



Digital Camera Attachment (Optional)

INVERTED MICROSCOPE FOR TISSUE CULTURE APPLICATIONS

## Specifications :

| MODEL INVI     | MODEL DESCRIPTIONS  |                                      |
|----------------|---|--------------------------------------|
| VIEWING HEAD   | Trinocular Head Inclined at 30 deg., Interpupillary distance 48 - 75mm      |                                      |
| EYEPIECE       | High-point, Extra Wide Field Eyepiece EW10X/22                              |                                      |
| OBJECTIVE      | LWD Plan Infinity Objective   | 4X/0.1 WD 18mm                       |
|                |   | 40X/0.6 WD 2.6mm (Cover Glass 1.2mm) |
|                | LWD Plan Infinity Phase Objective   | PH10X/0.25 WD 10mm                   |
|                |   | PH20X/0.4 WD 5.1mm                   |
| NOSEPIECE      | Quintuple Nosepiece   |                                      |
| CONDENSER      | ELWD Condenser NA 0.3, LWD 72mm, (without condenser 150mm)                  |                                      |
| PHASE ANNULUS  | 10X - 20X Phase Annulus Plate   |                                      |
| STAGE          | Plane Stage 160 x 250mm   |                                      |
|                | Glass Insert  |                                      |
|                | Auxiliary Stage 70 x 180mm  |                                      |
| FOCUSING       | Coaxial Coarse and Fine Adjustment  |                                      |
|                | Coaxial Stroke: 37.7mm per Rotation, Fine Stroke: 0.2mm per Rotation        |                                      |
| ILLUMINATION   | Halogen Lamp 6V30W  |                                      |
| FILTER         | Blue, Green and Frosted Glass, 45mm dia                                     |                                      |
| OPTIONAL ITEMS | PH40X/0.6 WD2.6mm (Cover Glass 1.2mm) Phase Contrast Objective              |                                      |
|                | Attachable Mechanical Stage, X-Y Co-axial Control, Moving Range 120 x 78 mm |                                      |
|                | Terasaki Holder, 38mm dia Petri Dish Holder, 54mm dia Slide Glass Holder    |                                      |
|                | Warm Plate  |                                      |
|                | Digital Camera / MIPS (Micro Image Projection System) Attachment            |                                      |
|                | Time Lapse Recording System   |                                      |



### MAGNUS ANALYTICS

A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886764, 30886739, Fax: 91-11-26959382

E-mail: [exports@magnusanalytics.com](mailto:exports@magnusanalytics.com) Website: [www.magnusanalytics.com](http://www.magnusanalytics.com)

All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment.

Magnus

invi

Long Working Distance (LWD), Plan Infinity Optics

4x, 40x Objectives & Phase 10x, 20x as Standard

Pre-centred Phase Annulus

Trinocular Port in  
Standard Unit

FOV 22mm with 10x Eyepiece

Optional Items :

Phase 40x Objective

Attachable Mechanical Stage

Warm Plate

Terasaki, Petri Dish

& Slide Glass Holders



Digital Camera Attachment (Optional)

INVERTED MICROSCOPE FOR TISSUE CULTURE APPLICATIONS

## Specifications :

| MODEL INVI     | MODEL DESCRIPTIONS  |                                      |
|----------------|---|--------------------------------------|
| VIEWING HEAD   | Trinocular Head Inclined at 30 deg., Interpupillary distance 48 - 75mm      |                                      |
| EYEPIECE       | High-point, Extra Wide Field Eyepiece EW10X/22                              |                                      |
| OBJECTIVE      | LWD Plan Infinity Objective   | 4X/0.1 WD 18mm                       |
|                |   | 40X/0.6 WD 2.6mm (Cover Glass 1.2mm) |
|                | LWD Plan Infinity Phase Objective   | PH10X/0.25 WD 10mm                   |
|                |   | PH20X/0.4 WD 5.1mm                   |
| NOSEPIECE      | Quintuple Nosepiece   |                                      |
| CONDENSER      | ELWD Condenser NA 0.3, LWD 72mm, (without condenser 150mm)                  |                                      |
| PHASE ANNULUS  | 10X - 20X Phase Annulus Plate   |                                      |
| STAGE          | Plane Stage 160 x 250mm   |                                      |
|                | Glass Insert  |                                      |
|                | Auxiliary Stage 70 x 180mm  |                                      |
| FOCUSING       | Coaxial Coarse and Fine Adjustment  |                                      |
|                | Coaxial Stroke: 37.7mm per Rotation, Fine Stroke: 0.2mm per Rotation        |                                      |
| ILLUMINATION   | Halogen Lamp 6V30W  |                                      |
| FILTER         | Blue, Green and Frosted Glass, 45mm dia                                     |                                      |
| OPTIONAL ITEMS | PH40X/0.6 WD2.6mm (Cover Glass 1.2mm) Phase Contrast Objective              |                                      |
|                | Attachable Mechanical Stage, X-Y Co-axial Control, Moving Range 120 x 78 mm |                                      |
|                | Terasaki Holder, 38mm dia Petri Dish Holder, 54mm dia Slide Glass Holder    |                                      |
|                | Warm Plate  |                                      |
|                | Digital Camera / MIPS (Micro Image Projection System) Attachment            |                                      |
|                | Time Lapse Recording System   |                                      |



### MAGNUS ANALYTICS

A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India  
 Tel.: 91-11-30886764, 30886739, Fax: 91-11-26959382  
 E-mail: [exports@magnusanalytics.com](mailto:exports@magnusanalytics.com) Website: [www.magnusanalytics.com](http://www.magnusanalytics.com)

All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer.  
 Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment.

### Magnus – Fraen Partnership

Since its inception in 2000, Fraen Srl has defined the concept of MAXIMIZING LIGHT! and is committed to developing high-quality, innovative and technologically advanced optical solutions to satisfy customer requirements.

Based on its rich knowledge of optical systems, Fraen Srl has developed automotive lighted instrumentation pointers, fiber optic couplers, light pipes, high-efficacy TIR collimators, etc and has become the industry benchmark for the top quality optical solutions for high-powered LEDs.

Magnus has a significant presence in the healthcare & education sectors in the Indian market and its microscopes have been the preferred choice for a number of World Bank funded projects. Since 1995, more than 15,000 Magnus microscopes have been supplied to important World Bank public health programs.

Guided by Japanese production methods and manufactured in a TUV-certified ISO 9001 : 2000 facility, Magnus microscopes are today recognized for their precision engineering and high-performance optics. The core team at Magnus has been specially trained in Japan & Germany by experts from leading optical engineering companies in the world.

The Magnus-Fraen partnership strives to combine the optical design and manufacturing strengths to produce high precision and cost competitive products to serve the education and healthcare markets all over the world.

### MicroLED Specifications

|                   |  |
|-------------------|--|
| Lifetime          | : >30,000 hrs  |
| LED Power         | : Typically 3W, depending on LED type  |
| Excitation        | : See table ▶  |
| Modules available | : Clamp-on modules available for the following microscopes:<br>• Magnus Icon • Magnus MLXi • Olympus CX21 • Olympus CX31 |
| Emission filters  | : 2 or 3 positions sliding filter carrier depending on the excitation spectra  |
| Mirror            | : Enhanced Al + SiO coating  |
| AC adaptor        | : Input Voltage 220/110V AC<br>Output Voltage - 7.5V DC / 12V DC<br>Power - Max. 15 - 18W                                |

| LED Cassettes | Excitation |
|---------------|------------|
| Royal Blue    | 450nm      |
| Blue          | 480nm      |
| Green         | 535nm      |

### Optionals

|                   |  |
|-------------------|--|
| Battery Pack      | : Rechargeable NiMH battery pack for Fluorescence module   |
| Intensity Control | : Fluorescence illumination intensity control through variable potentiometer   |
| Digital Cameras   | : Choice of Digital SLR camera or USB/Firewire camera.<br>(Requires Trinocular Head and adapter as per microscope model) |



# MicroLED

Fluorescence in a new light

## Magnus

## FRAEN

MAXIMIZING LIGHT

### MAGNUS ANALYTICS

A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India  
Tel.: 91-11-30886764, 30886739, Fax: 91-11-26959382  
E-mail: exports@magnusanalytics.com Website: www.magnusanalytics.com

## Magnus

## FRAEN

MAXIMIZING LIGHT

## THE NEW VISION OF FLUORESCENCE MICROSCOPY

### Overview

Magnus Analytics and Fraen Corporation, Italy bring you an integrated optical solution for fluorescence microscopy consisting of a unique, proprietary illumination system with high power solid-state (LED) sources to replace the mercury and xenon arc-lamps found in traditional epi fluorescence microscopy.

This approach allows significant increase of performance and light source lifetime, reduction of initial costs and operating costs, reduction of maintenance and less heat production.

The module is designed to attach to a number of standard bright field microscope and fluorescence microscopy can be done by simply inserting a mirror in the light path.

Bright field microscopy is not affected since the halogen white light function remains intact. Magnus MicroLED LED fluorescence modules are available for the Olympus CX series and the Magnus Icon as well as MLXi microscopes.

The standard fluorescence module is available with 480nm (blue) module. Other options include the 455nm (deep Blue) and the 535nm (green) LED cassettes.



MicroLED attachment on a Magnus MLXi microscope. The Digital SLR camera with appropriate adapter mounted on a Trinocular head provides high resolution images for archival and sharing. Options are available to attach other cameras for live image view on a laptop


### Tuberculosis Application

Transmitted light microscopy on sputum samples is the most widely used method to diagnose pulmonary tuberculosis. However, this method is complex and has low sensitivity compared to culture, while the more-sensitive fluorescence microscopy method is a far more reliable and effective diagnostic tool.

Till now, the fluorescence method has not found favour in public health programmes, due to the high initial investment required for a fluorescence microscope. But now, with MicroLED, fluorescence microscopy can be employed on small microscopes at an affordable price.

### Advantages of Fluorescence Microscopy over ZN staining

- Observation at low magnification provides high throughput – Due to the use of 40x magnification in fluorescence instead of 100x, the user can view a much larger field and thus make the process of diagnosis faster.
- The output of the Royal Blue LED (455nm) used for tuberculosis applications provides a perfect match for the excitation peak of the Auramine O Dye, resulting in a high contrast image with excellent Signal-to-Noise ratio




- With the use of a dry 40x objective, there is no need for using oil resulting in overall savings as well as lesser microscope maintenance issues associated with the use of immersion oil.
- The fluorescence method achieves higher sensitivity than the ZN method making detection of TB pathogens easier.


### Key Benefits

The MicroLED has been developed to provide equivalent performance and capability delivered in standard fluorescence microscopy equipment, but with a series of enhancements designed to make the technology accessible to more users, easier to operate and maintain, and significantly smaller to make it portable. Some of the key advantages offered by this technology are as under :

- The LED modules are light sources emitting an extremely efficient spectrum only in the desired bandwidth, thus ensuring a very good signal-to-noise ratio.
- Light source lifetime: typically 30,000 hrs, thus allowing many years of operation and cost savings.
- No warm-up time required for the light source.
- No need of any special alignment procedure.
- Variable light control allows adjustment of illumination intensity to reduce photobleaching
- Allows transmitted light observation without removing the fluorescence module
- Choice of Blue, Royal Blue and Green LED cassettes
- Battery pack option for field operation.



The mycobacteria appear as bright luminous rods on a dark background.



Single Colour Excitation Muscle, Alexa Fluor 546.

DISCOVER A NEW WORLD

**Magnus**  
MLM



**Anti-fungus treated optics:**  
Unique chemical treatment keeps fungus away

**Parfocal & Centered Objectives:**  
Minimises use of the fine focussing and stage-control knobs during objective change-over

**Mechanical Stage:**  
Sturdy, double-plate mechanical stage design

**Tamper Proof:**  
Eyepiece and condenser are pre-fixed to microscope body eliminating the possibility of loss & misplacement

MICROSCOPE FOR EDUCATION

## Specifications:

| ITEM                                 | SPECIFICATIONS   | MLM |
|--------------------------------------|--|-----|
| Body                                 | Aluminium die-cast monocular body  | ●   |
| Inclined Observation Head            | 45 degree monocular, rotatable through 360 degrees   | ●   |
| Eyepiece (widefield) for observation | WF 10x (F.N. 18) paired compensating eyepiece. Provides relief from eye fatigue and renders color-compensated images of utmost clarity. Compatible with an optionally available eyepiece micrometer.   | ●   |
| Nosepiece                            | Quadruple revolving nosepiece.   | ●   |
| Objectives                           | Achromat objectives 4X/N.A.0.10, 10X/N.A.0.25, 40X/N.A.0.65, (spring loaded) & 100X/N.A.1.25, (oil immersion)  | ●   |
| Mechanical stage                     | Left hand co-axial double plate mechanical stage, movement 60mm X 25 mm, size 120mm X 120mm  | ●   |
| Focussing system                     | Separate fine and coarse focussing adjustment knobs with tension control mechanism   | ●   |
| Condenser                            | Condenser N. A. 1.25 with iris diaphragm focusable by spiral movement  | ●   |
| Illumination base                    | Built-in illumination base with pre-centered 6V 15W halogen light source. Pre-centered bulb is coupled with an efficient collector lens system to provide optimum brightness along the optical path. A conveniently positioned rotatable knob enables variable light control | ●   |



OLYMPUS (INDIA) PVT. LTD.  
Works - A-3, Sector-61, Phase-II, Noida (UP)  
Corporate Office - S. Mahabaleshwar, Infratech Estate, Mehra Road, New Delhi-110044  
Tel: 91-11-30886766; 30886741; 30886744 Fax: 91-11-26959392  
E-mail: oia@olympusimg.com

**Magnus**



# Magnus MIPS

Computer Connectivity with USB2 Port

Connectivity to LCD Projection for Teaching and Training Application

Can be Mounted on a Binocular Microscope

Large Field of View

Plug and Play Design

Option of Built-in Micrometer Scale

Option of Software for Image Analysis

MIPS+MLX+PC

CE MICROSCOPE IMAGE PROJECTION SYSTEM

About Magnus:

World class manufacturing needs world class production infrastructure with a well trained and motivated work force - MAGNUS ensures it has one. Magnus microscopes are today recognized for their precision engineering and high-performance optics. Magnus has an uncompromising commitment towards quality. A group of highly qualified and experienced engineering personnel work closely with a highly motivated workforce in providing products and services of international quality. Magnus brings you convenience & performance in one single exciting package equipped with the Long life LED Light source. Also, Magnus & Fraen Srl Corporation, Italy bring you an integrated optical solution for Fluorescence Microscopy. Magnus has a very significant presence in the healthcare & education sectors.

## Magnus

OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in



C MR 1K 01 12

Accessories shown may not be part of standard equipment • Specifications are subject to change without notice • Optics are Anti-fungus Treated & with Multi Layer Coatings.

# Magnus invi

Long Working Distance (LWD), Plan Infinity Optics

4x, 40x Objectives & Phase 10x, 20x as Standard

Pre-centred Phase Annulus

Trinocular Port in Standard Unit

FOV 22mm with 10x Eyepiece

Optional Items :

Phase 40x Objective

Attachable Mechanical Stage

Warm Plate

Terasaki, Petri Dish

& Slide Glass Holders

INVI



Digital Camera Attachment (Optional)

INVERTED MICROSCOPE FOR TISSUE CULTURE APPLICATIONS

# Magnus MLX-M

Parfocal & Centered Optics

Minimises use of the fine focusing and stage-control knobs during objective change-over

Optics with Multi-layer Coating

Maximises transmission of light for crisp & bright image

Illumination

Efficient light collector-lens system for optimum brightness

Micron Sensitive Movements

Enables easy manipulation of specimen

MLX-M Monocular Version

Options Available

- With halogen or LED light source
- Freedom series with LED light & battery backup



Multilayer Vacuum Coating Plant

CE MICROSCOPES FOR EDUCATION & LABORATORIES

# Magnus HB



The HB & HSA based on a compact, precise and practical design and are extremely reliable in mechanical and optical performance. These compact and reliable microscopes incorporated some of the most effective of international optical techniques.

| Model         | HB Laboratory Microscope  | HSA Student Microscope  |
|---------------|---|---|
| Body          | Monocular die-cast body inclinable upto a horizontal position (90°)   |   |
| Magnification | 100X – 1000X (1500X)  | 100X – 600X   |
| Stage         | Built-in graduated Mechanical stage 120 X 125 mm with convenient adjustment for manipulation of slides  | Plain square stage of 100 X 100 mm, with two clips, mechanical stage (Optional) |
| Nosepiece     | Quadruple   | Triple hole   |
| Focusing      | By coarse and fine focusing knobs   |   |
| Eyepieces     | Widfield Eyepiece (NWF 10X)   | Huygenian 10X & 15X   |
| Objectives    | Achromatic 10X, 40X & 100X (spring loaded, oil immersion)   | Achromatic 10X & 40X  |
| Condenser     | Sub-stage abbe condenser of 1.25 N.A. focusable with rack and pinion; fitted with an iris diaphragm and filter holder (for blue filter) to facilitate optimum adjustment of light | Built-in Sub-stage condenser with iris diaphragm                                |
| Illumination  | Plano-concave mirror in adjustable fork mount   |   |
| Finish        | Black durable acid resistant epoxy coated finish  | Textured grey   |
| Optional      | Eyepieces widefield 15X & Huygenian 5X, 10X & 15X   |   |

Model HB & HSA, now manufactured under 'Magnus' brand with no change in performance, quality & model name

MICROSCOPES FOR LABORATORIES & STUDENTS

# Magnus HSA



# Magnus Accessories



Heating Stage  
Magnus Digital Warm Stage  
(Temp. from 27°C to 45°C) - for MLX / CH20i



Darkfield Attachments  
Magnus Immersion Darkfield Condenser with Objective 100x having iris diaphragm - for MLX  
Magnus Dry Dark field Condenser - for MLX



Phase Contrast Attachments  
Magnus Turret Type Phase Contrast Attachment Model PC-4 with Objectives 10x, 20x, 40x & 100x (oil) - for MLX. Also available : Model PC-2 with Long Barrel achromatic phase objectives 10X & 40X - for MLX



Gems - Inspection Accessories  
Magnus Darkfield Attachment (with specimen holder) for inspection of Gems & Micro Samples under transmitted light (To be used with Stereo microscopes MS / MSZ)



Magnus Microscope Image Projection System Model MIPS-USB  
Consisting of light weight compact imaging sensor with in-built widefield optical microscope adapter & with inbuilt USB connectivity (without any power supply, plug-n-play design) can be fitted in any Binocular / Trinocular microscope (suitable for NT based OS)



Magnus Digital Camera System Consisting of :  
Olympus Interchangeable Lens Type Live View Digital Camera (above 12 mega pixels) alongwith one zoom lens, rechargeable battery, carry bag, compatible 4GB memory card & Magnus optically corrected microscope adapter. Note: Presently, the camera model E-PL1 is being supplied. In case if it is discontinued, the camera available at the time of supply will be supplied.



Image Analysis Software  
Magnus-Pro Image Analysis Software for arranging and measuring captured images, creating databases and preparing reports

Discover a new world



# Magnus

International Quality Microscopes

# Magnus MicroLED

Technical Collaboration  
with Fraen SRL, Italy

Long Life LED Allows Many  
Years of Operation  
(Life Time 30000 hrs.)

No Warm Up Time Required



Modules Available for Olympus  
CX Series and Magnus  
MLX-i / ICON Microscopes

Royal Blue, Blue & Green Excitations

Battery Pack Option for Field Operation

Also allows transmitted light brightfield  
observation without removing the  
fluorescence module

No need of any special alignment procedure

Choice of Digital Camera (For BF Applications) or USB 5mp  
CMOS Color Camera System (For Fluorescence Applications)  
(Requires Trinocular head and adapter as per microscope model)

| LED cassettes | Excitation |
|---------------|------------|
| ROYAL BLUE    | 450 nm     |
| BLUE          | 480 nm     |
| GREEN         | 535 nm     |



LED FLUORESCENCE ATTACHMENTS  
FOR TB & AUTO-IMMUNE TESTING

# Magnus ICON

Plan Infinity Optics

6V 20W Koehler illumination

Inward Nosepiece

Centrable Condenser



Dual Slide Holder

Option of Quintuple nosepiece

Seidentopf Head  
with FN 20

ICON Binocular Version



Also available in  
Trinocular version ICON-Tr

Digital Camera Attachment  
(Optional)

MICROSCOPES FOR RESEARCH APPLICATIONS

# Magnus MLX*i*

Plan Infinity Corrected Optics

Seidentopf Head

Wide Field Eye Piece  
with FN 20



International Build Quality

Options Available

- With halogen or LED source
- Freedom series with LED light & battery backup

MLXi Binocular Version



Also available in  
Trinocular version MLXi-Tr

Digital Camera Attachment  
(Optional)

MICROSCOPES FOR RESEARCH APPLICATIONS

# Magnus MLX-B

Parfocal & Centered Optics  
Minimises use of the fine focusing  
and stage-control knobs during  
objective change-over

Optics with Multi-layer Coating  
Maximises transmission of light  
for crisp & bright image



Illumination  
Efficient light collector-lens  
system for optimum  
brightness

Micron Sensitive Movements  
Enables easy manipulation of specimen

Options Available

- With halogen or LED light source
- Freedom series with LED light & battery backup

MLX-B Binocular Version



Also available in  
Trinocular version MLX-Tr

Digital Camera Attachment  
(Optional)

MICROSCOPES FOR EDUCATION & LABORATORIES

# Magnus MSZ

Parfocality  
One Time Focus for all  
Zoom Magnifications

Ergonomic Design for  
Ease of Operation

Zoom Ratio 1:7



Options for Reflected and Transmitted  
illumination - Top halogen  
lamp 6v15w adjustable  
brightness bottom  
fluorescent lamp 5w

MSZ-Tr Trinocular Version



MSZ-BI Binocular Version

MICROSCOPES FOR BIOLOGICAL & INDUSTRIAL USE


# Magnus MS13/MS24

Superior Optics for Enhanced  
Image Flatness and Contrast

Ergonomic Design for  
Ease of Operation

Long Working Distance of 105 mm

Parfocal Objectives



Smooth Magnification  
Changeover between  
1X/3X or 2X/4X

Options for Reflected  
and Transmitted illumination -  
Top halogen lamp 6v10w adjustable  
brightness bottom fluorescent lamp 5w

Option for  
MIPS Attachments


MICROSCOPES FOR BIOLOGICAL & INDUSTRIAL USE

**Specifications: CX21i Biological Microscope**

|                      |  |
|----------------------|--|
| Microscope Frame     | CX21 FS1   |
| Optical system       | UIS2 (Universal Infinity System) optical system  |
| Illumination System  | Built-in transmitted illumination system<br>6 V 20 W halogen bulb<br>100-240 V 50/60 Hz universal voltage  |
| Focusing             | Stage height movement (coarse movement stroke 20 mm)<br>Fine focus graduation: 2.5 $\mu$ m   |
| Revolving Nosepiece  | Fixed quadruple nosepiece  |
| Stage                | Wire movement mechanical fixed stage: 120x132 mm<br>Traveling range: 76 mm (X) x 30 mm (Y)<br>Single specimen holder   |
| Observation Tube     | 30° inclined binocular tube<br>Interpupillary distance adjustment range 48–75 mm   |
| Condenser            | Abbe type with aperture iris diaphragm N.A.: 1.25  |
| Objective            | Plan Achromatic objectives (anti-fungus)<br>4x N.A.: 0.10 W.D.: 18.5 mm<br>10x N.A.: 0.25 W.D.: 10.6 mm<br>40x N.A.: 0.65 W.D.: 0.6 mm<br>100x N.A.: 1.25 W.D.: 0.13 mm (option)             |
| Eyepiece (10x)       | Field Number (F.N.): 20 (anti-fungus)  |
| Optional Accessories | Mirror unit, 15x eyepiece (F.N. 12, anti-fungus), micrometer, filter holder, darkfield stop, Micro image projection system, Phase-contrast Attachment/10X, 40X, Simple Polarizing attachment |

**Dimensions**

(unit: mm)



- Over 200,000 highly satisfied CX21 customers world over
- CX21i is now exclusively for India at affordable cost

CX21i is the environmental conscious product according to OLYMPUS' own standards.

- Adoption of cardboard for packing materials without styrene foam for promoting the recycling.
- Lead-free and arsenic-free Eco-glass for optics, such as lenses and prisms.

• The CX21i meets CE standards for safety.

- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

CX21i is manufactured under license from Olympus Corporation, Japan

**OLYMPUS****OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.**

Works : A-3, Sec-81, Phase-II, Noida-201305 U.P (India)

Regd. Office:

A-5, Mohan Co-operative Indl. Estate, (Opp. Sarita Vihar), Mathura Road, New Delhi-110044

Ph.: +91-11-30886743, 30886744 Fax: +91-11-26959382

Email: sales@olympusoptosystems.in Website: www.olympusoptosystems.in

**Features:**

- Ergonomic & compact design for user convenience
- UIS2 Plan infinity optics providing image flatness
- Seidentopf observation head for individual choice of adjustment
- Anti fungus treatment for component durability
- High eyepoint design eyepiece (F.N.20)
- High performance aspheric lenses in the abbe condenser & light relay system for bright & uniform illumination
- SMPS circuit for constant voltage output
- Rackless stage for durability and ease of use
- Mechanical stage focus-lock prevents mishaps
- Component security – Student proof

**OLYMPUS**

Your Vision, Our Future

Biological Microscope

**CX21i**

CX2 Series

**UIS2**  
World-leading optics

**Step up to  
Higher Performance  
with the Olympus  
UIS2 Infinity Optics**

**The CX21i sets new standards for  
educational and laboratory applications**

CA-1K-0812




# Superb image clarity coupled with Olympus' renowned reliability and rigidity

## Superior imaging

### Plan objectives are packaged as standard, providing image flatness that is among the best in this class

The CX21i is equipped with the same UIS2 optical system used in Olympus' top-level microscopes. In addition, Plan Achromat objectives are included as standard for the first time in this class of microscope, providing sharp images with high contrast right up to the edge of the field of view.



### Bright, uniform illumination

The CX21i's 6V20W high intensity halogen lamp delivers clear, stable illumination. The built-in aspheric collector lens provides images that are bright and uniformly illuminated over the entire field of view.

### Abbe condenser for optimized contrast


The optimal aperture stop position for each objective is clearly marked on the Abbe condenser, resulting in high-resolution, high-contrast images.



## Easy, safe operation

### Rackless stage for durability and ease of use

The wire-driven stage movement avoids any protrusion of racks — preventing the possibility of accidental damage or hand injury. Abrasion-resistant rackless design ensures continuous smooth movement.



### Component security — Student proof!

The eyepieces, objectives and condenser are all factory attached to the microscope body ensuring that no component is dropped or detached during transportation or regular use.



### Ultra-smooth, quadruple revolving nosepiece

High-precision machining ensures smooth operation and durability. Nosepiece includes a rubber ring for an easy ergonomic grip.

### Mechanical stage focus-lock prevents mishaps


The focus position can be locked making it easy to refocus when the specimen is changed. With the upper stage position locked, there is no chance of the objective accidentally touching the slide and damaging either the specimen or the objective itself.



## Use in any environment

### Anti-fungus treatment for component durability

An effective anti-fungus treatment is applied to the objectives, eyepieces and observation tube for consistent image clarity and long operating life even in hot humid work conditions.




### Compact design, easy to store and carry

The CX21i is only 391 mm high, compact enough to be stored in a typical classroom cabinet.



### Individual choice of adjustment

Binocular observation tube is inclined 30 degrees allowing for an ergonomic, comfortable posture during observation. Eyepieces can be easily set for multiple users with a scaled interpupillary distance adjustment and diopter setting. High eyepoint design eyepieces with F.N. 20 can be used comfortably without removing eyeglasses.




## Suitable for any user

### Tension adjustment for smoother focusing

The coaxial coarse/fine focusing knob is operable from either the left or right side. Coarse movement tension can be adjusted in accordance with personal preference.

### Ergonomic design for user convenience

The focusing knob, light control adjustment knob and stage handle are located close together, the user can work with hands on the desk and operate the CX21i with minimal movement.



Magnus

MLXi

Plan Infinity Corrected Optics

Seidentopf Head

Wide Field Eye Piece with FN 20

International Build Quality

Options Available

- With halogen or LED light source
- Freedom series with LED light & battery backup



MLXi Binocular Version



Also available in Trinocular version MLXi-Tr  
Digital Camera Attachment (Optional)

MICROSCOPES FOR RESEARCH APPLICATIONS

# Optional Accessories



LED based Fluorescence Attachment



MIPS

## Specifications :

| ITEM  | SPECIFICATIONS  | MLXi<br>Binocular<br>Version | MLXi-Tr<br>Trinocular<br>Version |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
|---|---|------------------------------|----------------------------------|------|----|------|-------|-----|------|------|---------------------|------|--------|---------------------------|------|---------|---|---|
| Body  | Aluminium die-cast body with all critical movements based on ball bearing & wire guides thereby ensuring smooth & precise manipulation  | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Mechanical Stage  | Co-axial low drive mechanical stage (125mm x 145mm) (+/-5mm) with traverse area of 50mm x 76mm (+/-5mm)   | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Focusing System   | Co-axial coarse & fine controls with a focus adjustment and fine adjustment knobs. Coarse Focus range 28mm. Fine focus range 0.2mm  | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Condenser holder  | Rack & pinion mounted condenser holder  | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Condenser   | Abbe condenser with aperture iris diaphragm (N.A. 1.25 ) focusable with rack & pinion through 20mm and a continuously variable iris diaphragm with a removable blue filter for daylight observation   | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Illumination base with option (a), (b), (c)   | (a) Built-illumination base with pre-centered 6V 20W halogen light source coupled with an efficient collector lens system. Power supply 230V AC 50Hz single phase.  | MLXi                         | MLXi-Tr                          |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
|   | (b) LED Light source High brightness, longlife (30,000hrs) 1w LED.  | MLXi LED                     | MLXi-Tr LED                      |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
|   | (c) LED light source (with battery back-up) High brightness, longlife (30,000hrs) 1w LED. Battery back-up in-built NiMH Rechargeable batteries provide 6 to 8 hrs back-up on full charge.   | MLXi Freedom series          | MLXi-Tr Freedom series           |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Nose Piece  | Quadruple revolving nosepiece based on precision ball-bearing mechanism with positive click stop Quintuple nosepiece option is also available.  | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Objectives  | <table border="1"> <thead> <tr> <th>Plan Achromat Objectives</th> <th>N.A.</th> <th>W.D.</th> </tr> </thead> <tbody> <tr> <td>4X</td> <td>0.10</td> <td>25 mm</td> </tr> <tr> <td>10X</td> <td>0.25</td> <td>5 mm</td> </tr> <tr> <td>40X (spring loaded)</td> <td>0.65</td> <td>0.5 mm</td> </tr> <tr> <td>100X (oil, spring loaded)</td> <td>1.25</td> <td>0.14 mm</td> </tr> </tbody> </table> | Plan Achromat Objectives     | N.A.                             | W.D. | 4X | 0.10 | 25 mm | 10X | 0.25 | 5 mm | 40X (spring loaded) | 0.65 | 0.5 mm | 100X (oil, spring loaded) | 1.25 | 0.14 mm | ● | ● |
|   | Plan Achromat Objectives  | N.A.                         | W.D.                             |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| 4X  | 0.10  | 25 mm                        |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| 10X   | 0.25  | 5 mm                         |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| 40X (spring loaded)   | 0.65  | 0.5 mm                       |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| 100X (oil, spring loaded)   | 1.25  | 0.14 mm                      |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Infinity corrected plan optics<br>Uniformly centered, Interchangeable & Parfocal<br>Anti fungus treated<br>Tropicalized anti fungus treatment ensures image excellence for long periods in conditions favoring to fungus growth |   |                              |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Inclined Observation Head<br><small>With a special anti-fungus treatment and an anti-reflection optical coating of the prism (to enhance the image brightness)</small>  | Binocular (30 degree inclined seidentopf), 360 degree rotatable, diopter adjustment   | ●                            |                                  |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
|   | Trinocular head 30 degree inclined seidentopf   |                              | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |
| Eyepiece (wide field) for observation   | WF 10x (FN 20mm) paired eyepiece. The unique optical design of the compensating eyepiece provides relief from eye fatigue and renders color-compensated wide-field images of utmost clarity. Compatible with optionally available eyepiece micrometer   | ●                            | ●                                |      |    |      |       |     |      |      |                     |      |        |                           |      |         |   |   |



OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in



All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment. • CE Certification is only for MLXi, MLXi-Tr

# Magnus

## MLX-B

### Parfocal & Centered Optics

Minimises use of the fine focusing and stage-control knobs during objective change-over

### Optics with Multi-layer Coating

Maximises transmission of light for crisp & bright image

### Illumination

Efficient light collector-lens system for optimum brightness

### Micron Sensitive Movements

Enables easy manipulation of specimen

### Options Available

- With halogen or LED light source
- Freedom series with LED light & battery backup




Also available in Trinocular version MLX-Tr  
Digital Camera Attachment (Optional)

MICROSCOPES FOR EDUCATION & LABORATORIES

# Specifications :

| ITEM                                 | SPECIFICATIONS   | MLX-B<br>Binocular<br>Version | MLX-Tr<br>Trinocular<br>Version |
|--------------------------------------|--|-------------------------------|---------------------------------|
| Body                                 | Aluminium die-cast body with all critical movements based on ball bearing & wire guides thereby ensuring smooth & precise manipulation   | ●                             | ●                               |
| Inclined Observation Head            | Binocular 45 degree inclined, rotatable through 360 degrees  | ●                             |                                 |
|                                      | Trinocular 45 degree inclined, rotatable through 360 degrees   |                               | ●                               |
| Eyepiece (widefield) for observation | WF 10x (FN 18) compensating paired eyepiece. Provides relief from eye fatigue and renders color-compensated images of utmost clarity. Compatible with an optionally available eyepiece micrometer    | ●                             | ●                               |
| Nosepiece                            | Quadruple revolving nosepiece based on precision ball-bearing mechanism with positive click stop   | ●                             | ●                               |
| Objectives                           | <b>Achromat objectives</b>   |                               |                                 |
|                                      |  | N.A.                          | W.D.                            |
|                                      | 4X   | 0.10                          | 29.0mm                          |
|                                      | 10X  | 0.25                          | 6.3mm                           |
|                                      | 40X (spring loaded)  | 0.65                          | 0.53mm                          |
| 100X (oil, spring loaded)            | 1.25   | 0.20mm                        |                                 |
| Mechanical stage                     | Stage size 125mm X 145mm with traverse area of 50mm X 76 mm with right hand low drive mechanical stage   | ●                             | ●                               |
| Focusing system                      | Co-axial coarse & fine controls with a focus adjustment range of 25 mm   | ●                             | ●                               |
| Condenser holder                     | Rack & pinion mounted condenser holder with height displacement upto 20mm  | ●                             | ●                               |
| Condenser                            | Abbe condenser with aperture iris diaphragm (N.A. 1.25 ) focusable with rack & pinion through 20 mm and a continuously variable iris diaphragm with a removable blue filter for daylight observation | ●                             | ●                               |
| Illumination base with option        | (a) Built-illumination base with pre-centered 6V 20W halogen light source coupled with an efficient collector lens system. Power supply 230V AC 50Hz single phase.                                   | MLX-B                         | MLX-Tr                          |
|                                      | (b) LED Light source High brightness, longlife (30,000hrs) 1w LED.   | MLX-B<br>LED                  | MLX-Tr<br>LED                   |
|                                      | (c) LED light source (with battery back-up) High brightness, longlife (30,000hrs) 1w LED. Battery back-up in-built NiMH Rechargeable batteries provide 6 to 8 hrs back-up on full charge.            | MLX-B<br>Freedom<br>series    | MLX-Tr<br>Freedom<br>series     |


## Optional Accessories • Simple Polarizing attachment • Eyepiece Widefield WF 15x • Image Analysis software



Micro Image Projection System  
MIPS



Digital Camera System



Phase Contrast Attachment PC-4



Darkfield attachment

- Imm. DF Condenser with OBJ 100x having Iris Diaphragm
- Dry. DF condenser

**Magnüs**

OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in



All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment. • CE Certification is only for MLX-B, MLX-Tr



# Magnus

## MLX-M

### Anti-Fungus treated Optics

Unique chemical treatment ensures image excellence for long periods in conditions favoring to fungus growth

### Parfocal & Centered Objectives

Minimises use of the fine focussing and stage-control knobs during objective change-over

### Optics with multi-layer coating

High-grade optics provides optimum brightness & contrast for long hours of comfortable viewing

### Micron Sensitive mechanical movements

Enables easy manipulation of specimen

### Illumination

Efficient light collector-lens system for optimum brightness

### Options Available

- With halogen or LED light source
- Freedom series with LED light & battery backup



Excellent optical clarity and compact size – which keeps all the key controls within easy reach – makes the Magnus MLX-M an ideal choice of microscope for education

MICROSCOPES FOR EDUCATION & LABORATORIES

## Specifications :

| ITEM  | SPECIFICATIONS   | MLX-M Monocular Version   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
|---|--|---|---------------------|------|------|----|------|--------|-----|------|-------|---------------------|------|--------|---------------------------|------|--------|---|
| Body  | Aluminum die-cast body with all critical movements based on ball-bearing & wire-guides thereby ensuring smooth & precise manipulation  | ●   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| Mechanical Stage                            | Co-axial low drive mechanical stage (125mm x 145mm) with traverse area of 50mm x 76mm  | ●   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| Focussing System                            | Co-axial coarse & fine focusing control with a focus adjustment and find adjustment knobs.   | ●   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| Condenser Holder                            | Rack & pinion mounted condenser holder.  | ●   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| Illumination base with option (a), (b), (c) | (a) Built-illumination base with pre-centered 6V 20W halogen light source coupled with an efficient collector lens system. Power supply 230V AC 50Hz single phase.   | MLX-M   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
|   | (b) LED Light source High brightness, longlife (30,000hrs) 1w LED.   | MLX-M LED   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
|   | (c) LED light source (with battery back-up) High brightness, longlife (30,000hrs) 1w LED. Battery back-up in-built NiMH Rechargeable batteries provide 6 to 8 hrs back-up on full charge.  | MLX-M Freedom series  |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| Nose Piece                                  | Quadruple nosepiece based on precision ball-bearing mechanism, enables smooth objective change.  | ●   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| CONDENSER                                   | Sub-stage Condenser With its compound lens system ensures that the traverse of light along the microscope's optical path is optimised. The built-in iris diaphragm enables maximum contrast of the specimen under observation.                                     | Abbe condenser with aperture iris diaphragm, N.A. 1.25, provided with a filter holder and blue filter.  | ●                   |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| OBJECTIVES                                  | Achromatic Objectives <ul style="list-style-type: none"> <li>• Anti fungus treated</li> <li>• Made from high quality Japanses optical glass</li> <li>• Precision engineered for parfocal &amp; centred viewing</li> </ul>  | <table border="1"> <thead> <tr> <th>Achromat objectives</th> <th>N.A.</th> <th>W.D.</th> </tr> </thead> <tbody> <tr> <td>4X</td> <td>0.10</td> <td>29.0mm</td> </tr> <tr> <td>10X</td> <td>0.25</td> <td>6.3mm</td> </tr> <tr> <td>40X (spring loaded)</td> <td>0.65</td> <td>0.53mm</td> </tr> <tr> <td>100X (oil, spring loaded)</td> <td>1.25</td> <td>0.20mm</td> </tr> </tbody> </table> | Achromat objectives | N.A. | W.D. | 4X | 0.10 | 29.0mm | 10X | 0.25 | 6.3mm | 40X (spring loaded) | 0.65 | 0.53mm | 100X (oil, spring loaded) | 1.25 | 0.20mm | ● |
| Achromat objectives                         | N.A.   | W.D.  |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| 4X  | 0.10   | 29.0mm  |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| 10X   | 0.25   | 6.3mm   |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| 40X (spring loaded)                         | 0.65   | 0.53mm  |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| 100X (oil, spring loaded)                   | 1.25   | 0.20mm  |                     |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| MONOCULAR INCLINED OBSERVATION HEAD         | <ul style="list-style-type: none"> <li>• 45 degree inclined</li> <li>• Rotatable through 360 degree</li> <li>a) With a special anti-fungus treatment</li> <li>b) With anti-reflection optical coatings of prisms to enhance the brightness of the image</li> </ul> | A unique design of observation head provides uniform illumination for ease of observation & eye comfort during extended usage.  | ●                   |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |
| WIDEFIELD EYEPIECE                          | WF 10x (F.No.18) compatible with an optionally available eyepiece micrometer   | Unique compensating eyepiece provides relief from eye fatigue & renders color-compensated images of utmost clarity. Color corrected in all magnifications including high power objective, resulting in better defined images.   | ●                   |      |      |    |      |        |     |      |       |                     |      |        |                           |      |        |   |

**Optional Accessories** Eyepiece WF 15x, Phase contrast, dark field and polarization accessories & imaging solutions with micro image projection system, Image analysis Software etc.

**Magnüs**

OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in

All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment. • CE Certification is only for MLX-M



Magnus

MSZ

Parfocality

One Time Focus for all  
Zoom Magnifications

Ergonomic Design for Ease of Operation

Zoom Ratio 1:7

Options for Reflected and Transmitted illumination -  
Top halogen lamp 6v15w adjustable brightness  
bottom fluorescent lamp 5w



MSZ-TR Trinocular Version



MSZ-BI Binocular Version

STEREO ZOOM MICROSCOPES FOR BIOLOGICAL & INDUSTRIAL USE

# Features

## High Performance:

The high performance MSZ Series stereo microscopes utilize a Greenough design and feature a 1 : 7 zoom ratio with a 100mm working distance and a large visual field. The parfocality allows an erect image of the observed specimens to stay in focus through repeated changes in the magnification. It is ergonomically designed for ease of operation.

## Durability:

In addition to performance and comfort, a highlight of the MSZ is its durability. Promoting one of the finest ball bearing guideways making the focusing mechanism extremely robust.

# Optional Accessories

## Accessory for Inspection of Gems & Micro Samples:

Darkfield Attachment (with specimen holder) for inspection of GEM & MICRO SAMPLES under transmitted light  
(To be used with Stereo Microscopes MSZ)

## Optional Lights:


- Fluorescent Ring Light
- Fibre Optics Illumination (Ring Light)
- Fibre Optics Illumination (Bifurcated)

## Auxiliary Objectives & Eyepieces:

- Eyepieces (Paired) Widefield 15x or 20x or 25x
- Auxiliary Lens 1.5x or 2x

# Standard Configurations

|                                    | <b>MODEL MSZ-BI</b>   | <b>MODEL MSZ-TR</b> |
|------------------------------------|---|---------------------|
| Optical body                       | Binocular body  | Trinocular body     |
| Zoom ratio                         | 1 : 7   |                     |
| Objective zoom range               | 0. 65X - 4.5X   |                     |
| Eyepiece                           | SWF10X/22mm (High eye-point Super Widefield Eyepieces)  |                     |
| Working distance                   | 100mm   |                     |
| Interpupillary distance adjustment | 55mm-75mm   |                     |
| Binocular head inclination of      | 45°   |                     |
| Diopter adjustment                 | +/- 5 diopter   |                     |
| Optical Body Rotation              | 360° with reflected and transmitted illumination Voltage input 220V/50HZ<br>Top halogen lamp 6V15W adjustable brightness bottom fluorescent lamp 5W |                     |



MIPS (Micro Image Projection System) :  
Optional attachments on MSZ & MSZ-TR



Gems - Inspection Accessories  
Darkfield attachment (with specimen holder) for inspection of  
Gems & Micro samples under transmitted light

# MSZ Series Optical Data

| Auxiliary Lens (optional) | WF 10X (standard) |            | WF 15X (optional) |            | WF 20X (optional) |            | WF 25X (optional) |            | Working Distance mm |
|---------------------------|-------------------|------------|-------------------|------------|-------------------|------------|-------------------|------------|---------------------|
|                           | Mag               | FOV mm     | Mag               | FOV mm     | Mag               | FOV mm     | Mag               | FOV mm     |                     |
| <b>1.5x</b>               | 9.8x - 67.5x      | 22.5 - 3.3 | 14.6x - 101.3x    | 16.4 - 2.3 | 19.5x - 135x      | 13.3 - 1.9 | 24.4x - 168.8x    | 12.3 - 1.7 | <b>45</b>           |
| <b>2.0x</b>               | 13x - 90x         | 16.9 - 2.5 | 19.5x - 135x      | 12.3 - 1.8 | 26x - 180x        | 10 - 1.4   | 32.5x - 225x      | 9.2 - 1.3  | <b>40</b>           |

The user may select different objective and eyepiece as per the requirement of magnification, working distance, visual field and depth of field.



OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in

All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment.



# Magnus MIPS

Computer Connectivity  
with USB2 Port

Connectivity to LCD Projection  
for Teaching and Training  
Application

Can be Mounted on a  
Binocular Microscope

Large Field of View

Plug and Play Design


Option of Built-in Micrometer Scale

Option of Software for Image Analysis




---

MICROSCOPE IMAGE PROJECTION SYSTEM



## The Perfect Educational Tool



Models available : • MIPS for direct connectivity to TV • MIPS-USB for USB 2 connectivity

**Magnüs**

OLYMPUS OPTO SYSTEMS INDIA PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida-201305 (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India

Tel.: 91-11-30886743, 30886744 Fax: 91-11-26959382

E-mail: sales@olympusoptosystems.in

All standard sets are supplied as per specifications which are subject to change without any obligation on the part of the manufacturer. Optics are anti-fungus treated & with multi-layer coatings. Accessories shown may not be part of standard equipment.



### Magnus – Fraen Partnership

Since its inception in 2000, Fraen Srl has defined the concept of MAXIMIZING LIGHT! and is committed to developing high-quality, innovative and technologically advanced optical solutions to satisfy customer requirements.

Based on its rich knowledge of optical systems, Fraen Srl has developed automotive lighted instrumentation pointers, fiber optic couplers, light pipes, high-efficacy TIR collimators, etc and has become the industry benchmark for the top quality optical solutions for high-powered LEDs.

Magnus has a significant presence in the healthcare & education sectors in the Indian market and its microscopes have been the preferred choice for a number of World Bank funded projects. Since 1995, more than 15,000 Magnus microscopes have been supplied to important World Bank public health programs.

Guided by Japanese production methods and manufactured in a TUV-certified ISO 9001 : 2008 facility, Magnus microscopes are today recognized for their precision engineering and high-performance optics. The core team at Magnus has been specially trained in Japan & Germany by experts from leading optical engineering companies in the world.

The Magnus-Fraen partnership strives to combine the optical design and manufacturing strengths to produce high precision and cost competitive products to serve the education and healthcare markets all over the world.

### MicroLED Specifications

|                   |  |
|-------------------|--|
| Lifetime          | : >30,000 hrs  |
| LED Power         | : Typically 3W, depending on LED type  |
| Excitation        | : See table ►  |
| Modules available | : Clamp-on modules available for the following microscopes:<br>• Magnus Icon • Magnus MLXi • Olympus CX21 • Olympus CX31 |
| Emission filters  | : 2 or 3 positions sliding filter carrier depending on the excitation spectra  |
| Mirror            | : Enhanced Al + SiO coating  |
| AC adaptor        | : Input Voltage 220/110V AC<br>Output Voltage - 7.5V DC / 12V DC<br>Power - Max. 15 - 18W                                |

| LED Cassettes | Excitation |
|---------------|------------|
| Royal Blue    | 450nm      |
| Blue          | 480nm      |
| Green         | 535nm      |

### Optionals

|                   |  |
|-------------------|--|
| Battery Pack      | : Rechargeable NiMH battery pack for Fluorescence module   |
| Intensity Control | : Fluorescence illumination intensity control through variable potentiometer   |
| Digital Cameras   | : Choice of Digital SLR camera or USB/Firewire camera.<br>(Requires Trinocular Head and adapter as per microscope model) |



**Magnus**

**FRAEN**  
MAXIMIZING LIGHT

### MAGNUS ANALYTICS

A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044 India  
Tel.: 91-11-30886764, 30826028, Fax: 91-11-26959382  
E-mail: sales@magnusanalytics.com Website: www.magnusanalytics.com

**Magnus**

**FRAEN**  
MAXIMIZING LIGHT

## THE NEW VISION OF FLUORESCENCE MICROSCOPY

### Overview

Magnus Analytics and Fraen Corporation, Italy bring you an integrated optical solution for fluorescence microscopy consisting of a unique, proprietary illumination system with high power solid-state (LED) sources to replace the mercury and xenon arc-lamps found in traditional epi fluorescence microscopy.

This approach allows significant increase of performance and light source lifetime, reduction of initial costs and operating costs, reduction of maintenance and less heat production.

The module is designed to attach to a number of standard bright field microscope and fluorescence microscopy can be done by simply inserting a mirror in the light path.

Bright field microscopy is not affected since the halogen white light function remains intact. Magnus MicroLED LED fluorescence modules are available for the Olympus CX series and the Magnus Icon as well as MLXi microscopes.

The standard fluorescence module is available with 480nm (blue) module. Other options include the 455nm (deep Blue) and the 535nm (green) LED cassettes.



MicroLED attachment on a Magnus MLXi microscope. The Digital SLR camera with appropriate adapter mounted on a Trinocular head provides high resolution images for archival and sharing. Options are available to attach other cameras for live image view on a laptop


### Tuberculosis Application

Transmitted light microscopy on sputum samples is the most widely used method to diagnose pulmonary tuberculosis. However, this method is complex and has low sensitivity compared to culture, while the more-sensitive fluorescence microscopy method is a far more reliable and effective diagnostic tool.


Till now, the fluorescence method has not found favour in public health programmes, due to the high initial investment required for a fluorescence microscope. But now, with MicroLED, fluorescence microscopy can be employed on small microscopes at an affordable price.

### Advantages of Fluorescence Microscopy over ZN staining


- Observation at low magnification provides high throughput – Due to the use of 40x magnification in fluorescence instead of 100x, the user can view a much larger field and thus make the process of diagnosis faster.
- The output of the Royal Blue LED (455nm) used for tuberculosis applications provides a perfect match for the excitation peak of the Auramine O Dye, resulting in a high contrast image with excellent Signal-to-Noise ratio



— Auromine-O Excitation  
— Auromine O Emission  
■ Royal Blue LED Spectrum



The mycobacteria appear as bright luminous rods on a dark background.



Single Colour Excitation Muscle, Alexa Fluor 546.

- With the use of a dry 40x objective, there is no need for using oil resulting in overall savings as well as lesser microscope maintenance issues associated with the use of immersion oil.
- The fluorescence method achieves higher sensitivity than the ZN method making detection of TB pathogens easier.

### Key Benefits

The MicroLED has been developed to provide equivalent performance and capability delivered in standard fluorescence microscopy equipment, but with a series of enhancements designed to make the technology accessible to more users, easier to operate and maintain, and significantly smaller to make it portable. Some of the key advantages offered by this technology are as under :

- The LED modules are light sources emitting an extremely efficient spectrum only in the desired bandwidth, thus ensuring a very good signal-to-noise ratio.
- Light source lifetime: typically 30,000 hrs, thus allowing many years of operation and cost savings.
- No warm-up time required for the light source.
- No need of any special alignment procedure.
- Variable light control allows adjustment of illumination intensity to reduce photobleaching
- Allows transmitted light observation without removing the fluorescence module
- Choice of Blue, Royal Blue and Green LED cassettes
- Battery pack option for field operation.



# OLYMPUS INDIA

*A Presentation*

WELCOME TO THE WORLD OF

OLYMPUS

&

Magnüs

*microscopes*

OLYMPUS & Magnüs

Series of microscopes are manufactured in a

✓TUV-Certified ISO 9001 : 2000

production facility by :

OLYMPUS INDIA PVT. LTD  
Noida

# OLYMPUS INDIA

Recognized for their

- ✓ precision engineering
- &
- ✓ high-performance optics

# OLYMPUS INDIA

## *Microscopes*

- ✓ A Preferred choice for World Bank Funded Projects
- ✓ More than 20,000 Microscopes supplied to Important World Bank public Health Programs

# OLYMPUS INDIA

- ✓ Influenced by a Japanese Legacy
- ✓ Uncompromising Commitment towards Quality

**Inauguration of new factory at NOIDA by  
Global Head, OLYMPUS Microscope, Japan**



**OLYMPUS - Noida factory**





# Horizontal 4-axis Machining Centre

[CNC operations with 5~10 microns accuracy]



# High Precision CNC Lathes



# High Speed Polishing

From OLYMPUS - Japan



Process Skills in Olympus India Manufacturing Unit



- Laser Interferometer

# Optics Assembly



# Optics Assembly 2



# Inclined Scope Assembly



# Olympus CH20i Assembly





# Opto-mechanical centring of Obj.



# CORPORATE OFFICE



# Marketing & Sales





**OLYMPUS**

OLYMPUS



**OLYMPUS**

**HSA**

STUDENT MICROSCOPE

Since the past 44 years Olympus (India) Pvt. Ltd. - in collaboration with Olympus Corporation, Tokyo, Japan has pursued the path of excellence in the field of optical engineering. To maintain international standards set for us by our Japanese collaborators, Olympus India has developed a completely integrated factory wherein optical lens processing, vacuum coating, finishing and other processes are carried out under the strict vigil of Japanese trained engineering personnel. This self-contained infrastructure alongwith a rigid quality control system ensures the production of high resolution optics and trouble-free instruments complying to the high OLYMPUS standards.

The OLYMPUS range of optical instruments comprise of:

- Biological Microscopes
- Research Microscopes

The HSA is based on a compact design and is extremely reliable in mechanical and optical performance. A wide variety of standard and optional accessories enhance the characteristics features of the HSA. The sphere of its application have been sufficiently extended to meet various requirements at high schools, colleges and universities.

|                      |  |
|----------------------|--|
| <b>Body</b>          | : Monocular die-cast body inclinable upto a horizontal position (90°)                            |
| <b>Magnification</b> | : 100X ~ 600X  |
| <b>Stage</b>         | : Plain square stage of 100 X 100 mm, with two clips, mechanical stage (Optional)                |
| <b>Nosepiece</b>     | : Triple hole nosepiece with positive click stops for precise and smooth alignment of objectives |
| <b>Focusing</b>      | : By coarse and fine focusing knobs  |
| <b>Eyepieces</b>     | : Huygenian 10X & 15X  |
| <b>Objectives</b>    | : Hard anti-reflection coated, colour coded Parfocal Achromat: 10X & 40X                         |
| <b>Condenser</b>     | : Built-in Sub-stage condenser with iris diaphragm   |
| <b>Illumination</b>  | : Plano-concave mirror in adjustable fork mount  |
| <b>Finish</b>        | : Olympus textured grey  |

# OLYMPUS



## HB

LABORATORY MICROSCOPE

This compact and reliable microscope incorporated some of the most effective of OLYMPUS optical techniques. The precise and practical design of this instrument has resulted in considerable extension of its applications in laboratories, universities and hospitals.

|                      |   |
|----------------------|---|
| <b>Body</b>          | : Monocular, die cast body, inclinable upto a horizontal position (90°)   |
| <b>Magnification</b> | : 100 X ~ 1000 X (1500 X)   |
| <b>Stage</b>         | : Built-in graduated Mechanical stage 120 X 125 mm with convenient adjustment for manipulation of slides  |
| <b>Nosepiece</b>     | : Quadruple nosepiece with positive click stops for precise and smooth alignment of objectives  |
| <b>Focusing</b>      | : By coarse and fine focusing knobs   |
| <b>Eyepiece</b>      | : Widefield Eyepiece (NWF 10X)  |
| <b>Objectives</b>    | : Hard anti-reflection coated, colour coded Parfocal Achromat: 10X, 40X & 100X (spring loaded, oil immersion)   |
| <b>Condenser</b>     | : Sub-stage Abbe condenser of 1.25 N.A. focusable with rack and pinion; fitted with an iris diaphragm and filter holder (for blue filter) to facilitate optimum adjustment of light |
| <b>Illumination</b>  | : Plano-concave mirror in adjustable fork mount   |
| <b>Finish</b>        | : Black durable acid resistant epoxy-coated finish  |
| <b>Optional</b>      | : Eyepieces widefield 15X & Huygenian 5X, 10X & 15X   |

# OLYMPUS



## GB

PATHOLOGICAL LABORATORY MICROSCOPE

This is a standard laboratory microscope which eliminates the operational complexities of conventional instruments, with features such as a graduated draw tube, and in-built graduated mechanical stage. The GB thus meets all the requirements of advanced microscopy.

|                      |   |
|----------------------|---|
| <b>Body</b>          | : Monocular die cast body, with graduated draw tube, inclinable upto a horizontal position (90°)  |
| <b>Magnification</b> | : 100 X ~ 1000 X (1500 X)   |
| <b>Stage</b>         | : Built-in graduated mechanical stage 120 X 125 mm with convenient adjustment for manipulation of slide   |
| <b>Nosepiece</b>     | : Quadruple nosepiece with positive click stops for precise and smooth alignment of objectives  |
| <b>Focusing</b>      | : By coarse and fine focusing knobs   |
| <b>Eyepiece</b>      | : Widefield Eyepiece (NWF 10X)  |
| <b>Objectives</b>    | : Hard anti-reflection coated, colour coded parfocal Achromat: 10X, 40X & 100X (spring loaded, oil immersion)   |
| <b>Condenser</b>     | : Sub-stage Abbe condenser of 1.25 N.A. focusable with rack and pinion; fitted with an iris diaphragm and filter holder (for blue filter) to facilities optimum adjustment of light |
| <b>Illumination</b>  | : Plano-concave mirror on adjustable fork mount   |
| <b>Finish</b>        | : Black durable acid resistant epoxy-coated finish  |
| <b>Optional</b>      | : Eyepieces Widefiled 15X, Huygenian 5X, 10X & 15X  |

OLYMPUS  
INDIA

Biological Microscope

CH20*i*



# Specifications

|                              |   | CH20i   |            |
|------------------------------|---|---|------------|
| Item                         | Specifications  | Binocular   | Trinocular |
| CH20i<br>MICROSCOPE<br>FRAME | <ul style="list-style-type: none"> <li>Coaxial coarse/fine knobs: Tension adjustment on the right side</li> <li>Fine focus knob graduated</li> <li>Stage movement (XY direction) on rack and pinion</li> <li>Quadruple revolving nosepiece (fixed) Plane stage 120 x 132 mm</li> <li>With right hand mechanical stage</li> <li>Abbe condensor N.A. 1.25 (oil immersion), with aperture iris diaphragm</li> <li>Blue filter</li> <li>Universal Power Supply (100V to 240V) for 6V 20W illuminator</li> <li>8cc immersion oil</li> <li>Dust cover</li> <li>Mirror unit (Plano-concave)</li> </ul> | <ul style="list-style-type: none"> <li>Binocular observation tube (inclination 45°, Interpupillary distance adjustment range 53-75 mm), diopter adjustment on the left</li> </ul> | ○          |
|                              | <ul style="list-style-type: none"> <li>Trinocular observation tube (inclination 45°), diopter adjustment on the left</li> </ul>   | ○   |            |
| Power Cord                   |   | ○   | ○          |
| Lamp                         | 6V 20W halogen lamp ( x 2)  | ○   | ○          |
| Objectives                   | NEA Achromat 4X (anti-fungus)   | ○   | ○          |
|                              | NEA Achromat 10X (anti-fungus)  | ○   | ○          |
|                              | NEA Achromat 40X (anti-fungus) spring   | ○   | ○          |
|                              | NEA Achromat 100X (anti-fungus) spring, oil   | ○   | ○          |
| Eye-piece                    | CWHK10X (LB eyepiece 10X), F.N. 18mm, (anti-fungus) ( x 2)  | ○   | ○          |

Specifications are subject to change without any obligation on the part of the manufacturer

CH20i is manufactured under license from Olympus Corporation, Japan

## OLYMPUS

OLYMPUS (INDIA) PVT. LTD.

A-5, Mohan Co-operative Industrial Estate, Mathura Road, New Delhi-110044  
Tel.: 30886766, 30886741, 30886744 Fax: 91-11-30886737 Email: oisales@dssimage.com



ISO 9001:2000

# OLYMPUS

BIOLOGICAL MICROSCOPE

## CH20i

In Japan, OLYMPUS prides itself as a company that manufactures quality optical products with a FOCUS ON LIFE. Olympus (India) extends that corporate philosophy as it brings the renowned Japanese technical expertise to its making of microscopes. It is progress through precision.




Quality that's pure

**OLYMPUS**



CH20i

## Features



Sealed optical parts

### Anti Fungus Optics

From the eyepiece and objectives to the interior of the observation tube, comprehensive anti fungus treatments are applied to every area which effects the clarity of the observed images. This Tropicalized treatment ensures image excellence for long periods in conditions favouring to fungus growth.




Main operating controls within easy reach

### Ergonomic Design

All the main operating functions are within easy reach, thereby allowing extended observation without fatigue of eye or posture.

### Window in the arm section

The specimen can be observed with the naked eye from the back. Also, convenient to observe the specimen while using microscope through natural light by using reflector.



Optimization of illumination with aspheric lenses

### High Performance Aspheric Lenses

The abbe condenser & the light relay system are fitted with high performance aspheric lenses, which ensure uniform illumination, resulting in a bright and crisp image.

### Illumination System

Built-in illumination through SMPS circuit for constant voltage output to cover input voltage fluctuations from 100 volts to 240 volts. This prolongs bulb life and prevent frequent bulb change.



Easily Interchangeable objectives

### Interchangeability of Objectives

The high precision objectives can be mounted on the nosepiece in any sequence without effecting their parfocality or centering.



## Technology From Japan

The CH-20i is an accumulation of advanced Japanese manufacturing technology that turns creative designs into quality products. Now, made in India under licence from Olympus Corporation, Japan, the CH-20i is the manifestation of the uncompromising quality standards and user friendly features which OLYMPUS products are known for.

The manufacturing process, under the strict supervision of a team of Japanese engineers, ensures compliance with demanding performance standards. OLYMPUS technology is constantly applied to every aspect of design, maximising the interchangeability of parts and minimising adjustment requirements during assembly. Consequently, the CH-20i is the only microscope in India conforming to internationally accepted standards.

CH20i



4-axis Machining Centre



Clean Room Assembly area



Laser Interferometer for optics inspection



Lens Inspection

Quality that's pure

**OLYMPUS**




CH20i

## Optional Accessories



### Digital Imaging Compliance

An optional adapter to mount digital cameras is provided to allow easy, cost-efficient digital imaging



### Darkfield central stop/ CH2-DS

Darkfield observation from 10X-40X is possible. This is placed underneath the condenser with attachment lens CH3-AL or filter holder CH2-FH

### Attachments


#### Phase-contrast Attachment/10X,40X

This is attached underneath the condenser. Phase-contrast objectives (10X & 40X) are provided

#### CCTV Attachment

Supplied with optical C mount Video Adapter

#### Micro Image Projection System



## Specifications

|                              |   | CH20i   |   |
|------------------------------|---|---|---|
| Item                         | Specifications  | Binocular   | Trinocular  |
| CH20i<br>MICROSCOPE<br>FRAME | <ul style="list-style-type: none"> <li>Coaxial coarse/fine knobs: Tension adjustment on the right side</li> <li>Fine focus knob graduated</li> <li>Stage movement (XY direction) on rack and pinion</li> <li>Quadruple revolving nosepiece (fixed) Plane stage 120 x 132 mm</li> <li>With right hand mechanical stage</li> <li>Abbe condensor N.A. 1.25 (oil immersion), with aperture iris diaphragm</li> <li>Blue filter</li> <li>Universal Power Supply (100V to 240V) for 6V 20W illuminator</li> <li>8cc immersion oil</li> <li>Dust cover</li> <li>Mirror unit (Plano-concave)</li> </ul> | <ul style="list-style-type: none"> <li>Binocular observation tube (inclination 45°, Interpupillary distance adjustment range 53-75 mm), diopter adjustment on the left</li> </ul> | <ul style="list-style-type: none"> <li>Trinocular observation tube (inclination 45°), diopter adjustment on the left</li> </ul> |
|                              |   | Power Cord  | ○   |
| Lamp                         | 6V 20W halogen lamp ( x 2)  | ○   | ○   |
| Objectives                   | iNEA Achromat 4X (anti-fungus)  | ○   | ○   |
|                              | iNEA Achromat 10X (anti-fungus)   | ○   | ○   |
|                              | iNEA Achromat 40X (anti-fungus) spring  | ○   | ○   |
| Eyepiece                     | iNEA Achromat 100X (anti-fungus) spring, oil  | ○   | ○   |
|                              | iCWHK10X (LB eyepiece 10X), F.N. 18mm, (anti-fungus) ( x 2)   | ○   | ○   |

Specifications are subject to change without any obligation on the part of the manufacturer

CH20i is manufactured under license from Olympus Corporation, Japan



ISO 9001:2000

# OLYMPUS

OLYMPUS (INDIA) PVT. LTD.

A-5, Mohan Co-operative Industrial Estate, Mathura Road, New Delhi-110044  
Tel.: 30886766, 30886741, 30886744 Fax: 91-11-30886737 Email: oisales@dssimage.com

# CH20i FEATURES

- **ERGONOMIC DESIGN**

- Aluminum die-cast body with new generation ergonomic design which **enables long hours of fatigue free work**

## **The slim body base design**

allows the user to operate the stage, coarse/fine focusing knob & illumination easily and comfortably while keeping the hands on the desk

# CH20i FEATURES

## A WINDOW IN THE ARM SECTION

**Allows the slide to be observed from the back**  
to see the objective which is currently engaged

The presence of window in the microscope frame plays an important role  
**for observation of the specimen from the back**  
while using external light to access the reflection mirror

# CH20i FEATURES

## HIGH EYEPOINT 10X WIDEFIELD EYEPIECES

- Easy to use by persons who wear glasses
- Eliminate Eye-Strain
- Enables fatigue free observation over long period of time
- Compensating Eyepieces  
Assure image sharpness in the visual field

## ANTI-FUNGUS TREATED - SEALED EYEPIECES

Hard anti-reflection coatings on eyepieces

Image excellence is long-lasting

## COLOR CORRECTIONS

At all magnifications

# CH20i FEATURES

## PARFOCALITY

**The precision engineered Parfocal objectives**

allow the specimen to remain in focus (within 1/4th turn of fine focus knobs) while changing objectives magnification via the revolving nosepiece.

## CENTERING

The precision engineered

**Objectives ensures that the specimen under observation remains in the centre of the field of view**

even while the objectives magnifications are changed via the revolving nose piece

## INTERCHANGEABILITY

**The objectives are interchangeable**

& can be mounted in any sequence or location on the quadruple revolving nose piece without effecting their Parfocality & Centering.

# CH20i FEATURES

## EXCELLENT DEFINITIONS

The higher image contrast is achieved by way of using special Objective design, Lens Curvature ratio setting, objective glass selection & lens construction followed by ultramodern multilayer coatings applied through a special lens coating machine to minimize lens surface reflections & eliminate flare and determine the best combination.

These are few ways Olympus achieves higher image contrast during Objective manufacture. This resulted

**minimum distance at which two minute dots in a specimen can be clearly distinguished**

## CO-AXIAL LOW DRIVE MECHANICAL STAGE

**HARD, ABRASION RESISTANT COATING ON STAGE**

**VERY SMOOTH & SENSITIVE MOVEMENT FOR EASE-OF OPERATION**

which enable accurate & precise manipulation of specimen slides

The Stage size 132 x 120mm, Movement 76 x 46mm

designed to accommodate the microscope slides which are 75x25 mm in size.

# CH20i FEATURES

- **OPTIMIZATION OF ILLUMINATION WITH ASPHERIC LENSE SYSTEM**

incorporated in Condenser & light relay system

**For brilliant and uniform light resulting in a  
Bright and crisp image**

The use of the aspheric lenses in the condenser & light relay system maximizes the collection of light emanating out of the bulb thereby condensing them as a strong light beam moving along the optical axis of the microscope which results in bright & even illumination of the observation field

- **Push fit filter holder**

**For safety & protection from breakage**

**TRULY DAY-LIGHT FILTER**

**for exceptional images contrast**



# CH20i FEATURES

- **BINOCULAR OBSERVATION TUBE**

Consists of a multilayer coated beam splitter prisms which ensures maximum transmittance / reflectance of light for negligible light loss and

**provides bright & Uniform illumination in Both the eyepieces for ease of observation & eye comfort during extended usage**

## **NO INTERNAL FLARE**

**Internal flare is completely eliminated**

thereby resulting in

**high contrast observation**

## **HIGH SENSITIVITY**

**Fine focus highly sensitive**

of sensitivity less than 2 microns

# CH20i FEATURES

- **EXCLUSIVE ANTI-FUNGUS TREATMENTS**

## **Maintain optimum performance**

Extensive anti-fungus application from the eyepiece and objectives to the interior of the observation tube, comprehensive anti fungus treatments are applied to every area which effects the clarity of the observation images

**Image excellence is extremely long-lasting even in conditions liable to fungus accumulation**

# CH20i FEATURES

## Overall Image Definition & Brilliance

Excellent at all magnifications throughout the field  
Seeing & Comparing is believing

## Collaboration

Manufactured under license agreement with  
a renowned International Co. of repute **Olympus  
Corporation, Tokyo, Japan**

## ISO 9000-2001

Manufactured in a ISO 9000-2001 certified  
company  
under strict vigil of Japanese trained quality control engineers

OLYMPUS  
INDIA

MLX*i*

Magnus  
MLX*i*

Infinity Corrected Optical System

Seidentopf Head

International Build Quality



MICROSCOPE WITH INFINITY CORRECTED OPTICS

*New introduction*

# Magnus

## MicroLED

Fluorescence in a New Light

[www.magnusanalytics.com](http://www.magnusanalytics.com)

**FRAEN**  
MAXIMIZING LIGHT

# Magnus MicroLED

New introduction

## MicroLED for Tuberculosis


Transmitted light microscopy on sputum samples is the most widely used method to diagnose pulmonary tuberculosis. However, this method is complex and has low sensitivity compared with culture, while the more-sensitive fluorescence microscopy method is a far more reliable and effective diagnostic tool.

Till now, the fluorescence method has not found favour in developing countries, due to the high initial investment required for a fluorescence microscope. But now, with MicroLED, fluorescence microscopy can be employed on small microscopes at an affordable price.


Furthermore it can be adapted on existing microscopes. Battery operation is also possible.

A unique, proprietary illumination system has been developed to utilize high-power solid-state (LED) light sources to replace the mercury arc-lamps used in traditional epi-fluorescence microscopy. The LED fluorescence module is designed to attach to a standard bright field microscope and is used in transmission mode. Bright field microscopy capability is not affected since the halogen white light function remains intact.

Spectral Response Curve of Auromine O Dye. LEDs emit an extremely efficient and narrow spectrum only in the desired bandwidth, thus producing excellent signal- to-noise ratio




▶ The mycobacteria appear as bright luminous rods on a dark background



MLXi with the MicroLED Attachment


MicroLED can be mounted on other microscopes like the Zeiss Axio Star\* as well as the Olympus CX31\* ▼



# Magnus MicroLED

New introduction

## MicroLED Overview



- ▶ Single Colour Excitation  
Muscle, Alexa Fluor 546

Magnus Analytics and Fraen Corporation, Italy bring you an integrated optical solution for fluorescence microscopy consisting of a unique, proprietary illumination system with high power solid-state (LED) sources to replace the mercury and xenon arc-lamps found in traditional epi-fluorescence microscopy.

This approach allows significant increase of performance and light source lifetime, reduction of initial costs and operating costs, reduction of maintenance and heat production.


The module is designed to attach to a standard bright field microscope and does not change the characteristics of the microscope in any way. The fluorescence light source is used in transmission mode and will not void any warranties. Bright field microscopy is not affected since the halogen white light function remains intact, which means that transmitted light observation is possible without major changes in the optical configuration.

Magnus MicroLED LED fluorescence is available as a complete ready-to-use instrument with the Magnus MLXi microscope or as an add-on kit for the following microscopes\* :

- Olympus
- Carl Zeiss
- Nikon
- Leica

The standard fluorescence module is available with 480nm (blue) module.

Other optional add-ons include 365nm (UV), 455nm (Deep Blue), 535nm (green), 590 nm (yellow) and 630nm (red) excitation wavelengths.



- ▶ Example of 3 colour Excitation  
BPAC Cells, DAPI / BODIPY, Mito Tracker Red

## MicroLED Benefits

- The MicroLED products were developed to provide equivalent performance and capability delivered in standard fluorescence microscopy equipment, but with a series of enhancements designed to make the technology accessible to more users, easier to operate and maintain, and significantly smaller to make it portable.
- The LED modules are light sources emitting an extremely efficient spectrum only in the desired bandwidth, thus ensuring a very good signal-to-noise ratio.
- No warm-up time required for the light source.
- Light source lifetime: typically 30000 hours, thus allowing many years of operation and cost savings.
- No need of any special alignment procedure.
- Variable light control for adjustment of illumination intensity allows regulation of each single color channel via adjustable electronic driver (reduce photobleaching).
- Allows transmitted light observation without removing the fluorescence module.
- One, two and three color excitation. Grab up to three colors in one time (no pixel shift).
- Smaller instrument footprint .
- Battery pack option for field operation.

# SYSTEMS SPECIFICATION

New introduction

## • Magnus MicroLED Cassettes

Customer determines appropriate excitation cassette(s) according to fluorescent dyes.

Lifetime : > 30,000 hours  
LED Power : Typically 3 W, depending on LED type  
Excitation  $\lambda$  : See table

| FRAEN AFTER cassettes | Excitation |
|-----------------------|------------|
| UV                    | 365 nm     |
| ROYAL BLUE            | 450 nm     |
| BLUE                  | 480 nm     |
| CYAN (on request)     | 505 nm     |
| GREEN                 | 535 nm     |
| YELLOW                | 590 nm     |
| RED                   | 630 nm     |

## • MicroLED Driver unit

Different driver units available for single, dual or triple color control.

AC adaptor : Input Voltage : 220/110V AC  
Output Voltage : 7.5VDC / 12V DC  
Power : Max. 15-18W

Features : Self detection of LED power  
Option : Battery pack for use on field

## • MicroLED Additional Optics

UV blocking filter : In filter carrier housing  
Emission filters : In 3 to 6 positions sliding filter carrier  
Mirror : Enhanced Al + SiO<sub>2</sub> coating

## • MicroLED Modules

Clamp-on adaptors for the following microscopes available:

- Olympus
- Carl Zeiss
- Nikon
- Leica

\*Micro LED modules for Olympus, Carl Zeiss, Nikon and Leica available for sale in India only.

# Magnus

**Magnus Analytics**

A-5 Mohan Coop. Industrial Estate Mathura Road New Delhi-110044 INDIA  
Voice : +91-11-30886739 Fax : +91-11-26959382 Mobile : +91-93505443204

E-mail : sales@magnusanalytics.com Website : www.magnusanalytics.com





# Magnus

## MLX-U unplugged

### Battery Backup:

A long-life LED coupled with a rechargeable battery allows 6 hrs usage without electrical power

### Anti-fungus treated optics:

Unique chemical treatment keeps fungus away

### Parfocal & Centered Objectives:

Minimises use of the fine focussing and stage-control knobs during objective change-over

### Optics with Multi-layer coating:


Maximises transmission of light for crisp, bright images



### Specifications:

|                                      |  |
|--------------------------------------|--|
| Body                                 | Aluminium die-cast body. All critical movements based on ball bearing & wire guides that ensure smooth & precise manipulation  |
| Inclined Observation Head            | Binocular 45 degree inclined, rotatable through 360 degrees  |
| Eyepiece (widefield) for observation | WF 10x (F.N. 18) compensating paired eyepieces that provides relief from eye fatigue and renders color-compensated images of utmost clarity. Compatible with an optionally available eyepiece micrometer |
| Nosepiece                            | Quadruple revolving nosepiece with click stop  |
| Objectives                           | Long barrel achromatic 4X, 10X, 40X (spring ) & 100X (oil immersion - spring)  |
| Mechanical stage                     | Stage size 125mm X 145mm with traverse area of 50mm X 76 mm with right hand low drive knobs  |
| Focusing system                      | Co-axial coarse & fine controls with a focus adjustment range of 25 mm and fine adjustment knobs   |
| Condenser holder                     | Rack & pinion mounted condenser holder with height displacement upto 20mm  |
| Condenser                            | Abbe condenser with aperture iris diaphragm (N.A. 1.25 ) focusable with rack & pinion through 20 mm and a continuously variable iris diaphragm   |
| Illumination base                    | Built-in illumination base with new 12V LED to provide adequate brightness at minimal power consumption. Supported by inbuilt battery backup during power failure  |

MICROSCOPES FOR FIELD



# Magnus

Manufactured by

OLYMPUS (INDIA) PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate,  
Mathura Road, New Delhi-110044

Tel.: 91-11-30886766, 30886741, 30886744 Fax: 91-11-26959382

E-mail: oisales@dssimage.com

OLYMPUS  
INDIA

MLX

Magnus :  
MLX

Anti-fungus treated optics:

Unique chemical treatment  
keeps fungus away

Parfocal & Centered Objectives:

Minimises use of the fine focussing  
and stage-control knobs during  
objective change-over

Optics with Multi-layer coating:

Maximises transmission of light  
for crisp, bright images

Illumination:

Condenser lens system  
for optimum brightness

Micron sensitive  
mechanical movements:

Enables easy manipulation of specimen



MICROSCOPES for Education & Laboratories

# Magnus MLX-B

## Key Features

- High-grade optics provides  
**Optimum brightness & contrast for long hours of comfortable viewing**
- **Optics are Uniformly Centered, Interchangeable And Parfocal**

Durable & rust-free construction with high quality paint finish ensures

**trouble free performance for years**

- Machining of components on high precision CNC machines to ensure better accuracy.
- **The ventilated light relay system resulting in the`minimum heating of base.**
- .

# Magnus MLX-B

## Key Features

- **Meets International safety standards of CE.**
- **Tropicalized Anti-Fungus treatment ensures image excellence** for long periods in conditions favoring to fungus growth.
- Coaxial focus mechanism provides precision at all magnifications.
- The superior design nosepiece system allows smooth rotation & easy access to specimens.
- A unique design of observation head provides **uniform illumination in both the eyepieces for ease of observation**  
&  
**eye comfort during extended usage.**

# Magnus MLX-B

## Key Features

- Unique compensating eyepiece

- Provides **Relief from eye fatigue** & renders color-compensated images of utmost clarity.

- Color corrected in all magnifications

- including high power objective, resulting in

- **better-defined images.**

OLYMPUS  
INDIA

MIPS

Magnus  
MIPS

Computer connectivity  
with USB Port

Connectivity to LCD  
projector for Teaching  
and Training applications

Can be mounted on a  
Binocular Microscope

Large field of view

Plug and play design

Option of built-in  
micrometer scale

Option of software for  
image analysis



MICROSCOPE IMAGE PROJECTION SYSTEM

## Key Features

## MIPS Micro Image Projection System

- **Universal design for adaptation of any microscope**  
on Binocular eyepiece or Trinocular port
- Control Communication for complete handling through PC.  
**Usable through Laptop as well as desktops.**
- Interline Transfer CCD Sensor  
**provides true colour image**
- **Connectivity through USB 2.0**  
provides 25 frames/sec.
- **Large view area - Selectable through software**  
(from 240 x 240 up to 720 x 576 pixels)
- Sensitivity < 1.0 Lux for  
**excellent low light performance**
- Signal to noise ratio > 48 dB (AGC off) for  
**high quality image**
- **Back light compensation for balancing the light on the specimen**
- CCD Iris Control and Shutter speeds from 1/50 to 1/100000 for  
**Automatic Manual adjustment of light intensity**
- **Auto White Balance**  
and Manual White Balance modes for accommodating different light sources like halogen, tungsten, etc
- **Power supply through USB bus**  
for convenient single cable operation

OLYMPUS  
INDIA

MSZ

Magnus  
MSZ

**Anti-fungus treated optics:**  
Unique chemical treatment  
keeps fungus away

**Parfocality:**  
One time focus for all zoom  
magnifications

Ergonomic design for ease of operation

Zoom Ratio 1:7



STEREOZOOM MICROSCOPES



# Magnus MSZ



**Anti-fungus treated optics:**  
Unique chemical treatment  
keeps fungus away

**Parfocality:**  
One time focus for  
all zoom magnifications

Ergonomic design for ease of operation

Zoom Ratio 1:7

STEREOZOOM MICROSCOPES

## Features

### High Performance:

The high performance MSZ Series stereo microscopes utilize a Greenough design and feature a 1:7 zoom ratio with a 100mm working distance and a large visual field. The parfocality allows an erect image of the observed specimens to stay in focus through repeated changes in the magnification. It is ergonomically designed for ease of operation.

### Durability:

In addition to performance and comfort, a highlight of the MSZ is its durability. Promoting one of the finest ball bearing guideways making the focusing mechanism extremely robust.

## Standard Configurations

|                                    | MODEL MSZ-BI  | MODEL MSZ-TR    |
|------------------------------------|---|-----------------|
| Zoom ratio                         | Binocular body  | Trinocular body |
| Objective zoom range               | 1:7   | 0.65X - 4.5X    |
| Eye-piece                          | SWF10X/22mm (High eye-point Super Widefield Eyepieces)  |                 |
| Working distance                   | 100mm   |                 |
| Interpupillary distance adjustment | 55mm-75mm   |                 |
| Binocular head inclination of      | 45°   |                 |
| Diopter adjustment                 | +/- 5 diopter   |                 |
| Optical Body Rotation              | 360° with reflected and transmitted illumination Voltage input 220V/50HZ<br>Top halogen lamp 6V15W adjustable brightness bottom fluorescent lamp 5W |                 |

## MSZ Series Optical Data

| Auxiliary Lens (optional) | WF 10X (standard) |            | WF 15X (optional) |            | WF 20X (optional) |            | WF 20X (optional) |            | Working Distance mm |
|---------------------------|-------------------|------------|-------------------|------------|-------------------|------------|-------------------|------------|---------------------|
|                           | Mag               | FOV mm     | Mag               | FOV mm     | Mag               | FOV mm     | Mag               | FOV mm     |                     |
| 1.5x                      | 9.8x - 67.5x      | 22.5 - 3.3 | 14.6x - 101.3x    | 16.4 - 2.3 | 19.5x - 135x      | 13.3 - 1.9 | 24.4x - 168.8x    | 12.3 - 1.7 | 45                  |
| 2.0x                      | 13x - 90x         | 16.9 - 2.5 | 19.5x - 135x      | 12.3 - 1.8 | 26x - 180x        | 10 - 1.4   | 32.5x - 225x      | 9.2 - 1.3  | 40                  |

The user may select different objective and eyepiece as per the requirement of magnification, working distance, visual field and depth of field.

## Optional Accessories

### Accessory for Inspection of Gems & Micro Samples

Darkfield Attachment (with specimen holder) for inspection of GEMS & MICRO SAMPLES under transmitted light  
(To be used with Stereo Microscopes MSZ)

## Optional Lights

Fluorescent Ring Light  
Fibre Optics Illumination (Ring Light)  
Fibre Optics Illumination (Bifurcated)

## Auxiliary Objectives & Eyepieces

Eyepieces (Paired) Widefield 15x or 20x or 25x  
Auxiliary Lens 1.5x or 2x

**Magnus**

OLYMPUS (INDIA) PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044

Tel.: 91-11-30886766, 30886741, 30886744 Fax: 91-11-26959382

E-mail: oisales@dssimage.com



## APPLICATIONS FOR STEREO MICROSCOPES

| <b>BIOLOGICAL APPLICATIONS</b>          | <b>Segment</b>                                     |
|---|--|
| Dissecting                              | Schools & Colleges                                 |
| Fossils                                 | Geology / Museums / Paleobotany                    |
| Insects                                 | Malaria deptt. / Schools / College                 |
| Plant and animal specimens              | Agricultural Universities                          |
| Botanical and anatomical studies        | Colleges / Medical Colleges                        |
| Medical                                 | Hair transplant and dental clinics                 |
| Forensic studies                        | Forensics  |
| Tissue                                  | Medical  |
| IVF                                     | Medical  |
| <b>INDUSTRIAL APPLICATIONS</b>          |  |
| Viewing collections/inspection of coins | Mints / Museums                                    |
| stony or metallic object                | Geology  |
| Yarns / Fibres                          | Textile Industries                                 |
| for studying minerals                   | Geology  |
| stamps                                  | Forensics  |
| Gems & jewelry                          | Gemologists / Jewellers & Gem Stones Manufacturers |
| Antique artifacts                       | Museums  |
| Electronic circuitry                    | Electronic Industries                              |
| Miniatures                              | Watch Industries                                   |
| Mold identification                     | Molding Industries                                 |
| Paint Manufacturing (R&D / QC)          | Paint Industries                                   |
| Syringe (QC)                            | Syringe Manufacturers                              |
| Inspection (QC)                         | Intracocular lenses                                |
| Inspection (QC)                         | Picture Tube (TV / Monitors) Manufacturers         |
| Pen Tips (QC)                           | Pen Industries                                     |

# Magnus MS13/MS24



Superior Optics for enhanced image flatness and contrast


Ergonomic design for ease of operation

Long working distance of 105 mm

Parfocal objectives

Smooth magnification changeover between 1X/3X or 2X/4X

STEREO MICROSCOPES



Micro Image Projection System

## Specifications :

|                     |   | MS 13                                  |     | MS 24 |     |
|---------------------|---|--|-----|-------|-----|
| Microscope Body     | Objective   | 1X/3X                                  |     | 2X/4X |     |
|                     | Working Distance  | 105mm                                  |     |       |     |
| Eye Piece           | Diopter Adjustment  | +/- (left side)                        |     |       |     |
|                     | High Eyepoint Super Widefield eyepiece SWH 10X with Field Number 23 |  |     |       |     |
| Field of View       | Objective   | 1X                                     | 3X  | 2X    | 4X  |
|                     | Visual field (mm)   | 23                                     | 7.7 | 11.5  | 5.7 |
| Light Source        | Reflected   | 6V 10W lamp with adjustable brightness |     |       |     |
|                     | Transmitted   | 5W fluorescent lamp                    |     |       |     |
| Optical Accessories |   | Fluorescent ring light                 |     |       |     |
|                     |   | Eyepiece Micrometer                    |     |       |     |
|                     |   | Eyepiece 15X and 20X                   |     |       |     |
|                     |   | MIPS (Micro Image Projection System)   |     |       |     |

**Magnus**

OLYMPUS (INDIA) PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044

Tel.: 91-11-30886766, 30886741, 30886744 Fax: 91-11-26959382

E-mail: oisales@dssimage.com



# Magnus MLM

DISCOVER A NEW WORLD



**Anti-fungus treated optics:**  
Unique chemical treatment keeps fungus away

**Parfocal & Centered Objectives:**  
Minimises use of the fine focussing and stage-control knobs during objective change-over

**Mechanical Stage:**  
Sturdy, double-plate mechanical stage design

**Tamper Proof:**  
Eyepiece and condenser are pre-fixed to microscope body eliminating the possibility of loss & misplacement

MICROSCOPE FOR EDUCATION

## Specifications :

| ITEM                                 | SPECIFICATIONS   | MLM |
|--------------------------------------|--|-----|
| Body                                 | Aluminium die-cast monocular body  | ●   |
| Inclined Observation Head            | 45 degree monocular, rotatable through 360 degrees   | ●   |
| Eyepiece (widefield) for observation | WF 10x (F.N. 18) paired compensating eyepiece. Provides relief from eye fatigue and renders color-compensated images of utmost clarity. Compatible with an optionally available eyepiece micrometer  | ●   |
| Nosepiece                            | Quadruple revolving nosepiece  | ●   |
| Objectives                           | Achromat objectives 4X/N.A.0.10, 10X/N.A.0.25, 40X/N.A.0.65, (spring loaded) & 100X/N.A.1.25, (oil immersion)  | ●   |
| Mechanical stage                     | Left hand co-axial double plate mechanical stage, movement 60mm X 25 mm, size 120mm X 120mm  | ●   |
| Focusing system                      | Separate fine and coarse focusing adjustment knobs with tension control mechanism  | ●   |
| Condenser                            | Condenser N. A. 1.25 with iris diaphragm focusable by spiral movement  | ●   |
| Illumination base                    | Built-in illumination base with pre-centered 6V 15W halogen light source. Pre-centered bulb is coupled with an efficient collector lens system to provide optimum brightness along the optical path. A conveniently positioned rotatable knob enables variable light control | ●   |

Magnus

OLYMPUS (INDIA) PVT. LTD.

Works : A-3, Sector-81, Phase-II, Noida (U.P.)

Corporate Office : A-5, Mohan Co-operative Indl. Estate, Mathura Road, New Delhi-110044

Tel.: 91-11-30886766, 30886741, 30886744 Fax: 91-11-26959382

E-mail: oisales@dssimage.com



# ML-M

## Best option for HB/GB users

### **Affordable Price Difference with Add-on features**

- **Inclined head**
- **double layer mechanical stage**
- **bright illumination**
- **long barrel optics**
- **quadruple nosepiece**
- **4 objectives**
- **ergonomically design**
- **temper proof optics**

# ML-M

## **KEY FEATURES :**

### **Anti-Fungus Optics**

ensures image excellence for long periods in conditions favoring to fungus growth

**Optics are uniformly centered, interchangeable & Parfocal  
Optimum brightness & contrast for long hours  
of Comfortable viewing.**

**Sub stage Condenser with Spiral Movement for ease  
of operation**

**Double Layer Mechanical Stage for comfortable  
movement of slide**

**Temper proof optics for security. Easy handling  
& Ergonomically design**

# Accessories



Heating Stage  
Magnus Digital Warm Stage  
(Temp. From 27°C to 45°C) - for MLX/CH20i



Darkfield Attachments  
Magnus Immersion Darkfield Condenser with Objective 100x having iris diaphragm - for MLX  
Magnus Dry Darkfield Condenser - for MLX



Phase Contrast Attachments  
Magnus Turret Type Phase Contrast Attachment  
Model PC-4 with Objectives 10x, 20x, 40x & 100x (oil) - for MLX



Gems - Inspection Accessories  
Magnus Darkfield Attachment (with specimen holder) for inspection of Gems & Micro Samples under transmitted light (To be used with Stereo microscopes MS / MSZ)



Magnus Microscope Image Projection System Model MIPS-USB  
Consisting of light weight compact imaging sensor with in-built widefield optical microscope adapter & with inbuilt USB connectivity (without any power supply, plug-n-play design) can be fitted in any Binocular/Trinocular microscope. (suitable for NT based OS)



Magnus Digital Camera System  
Magnus Digital Camera system complete with any renowned available brand digital camera (above 7 Million Pixels) with rechargeable batteries with charger, & with Magnus Optical Microscope adapter (Olympus digital camera is discontinued & not available at the moment)



Image Analysis Software  
Magnus-Pro Image Analysis Software for arranging and measuring captured images, creating databases and preparing reports

world class Microscopes

# Magnus

the Magnus range ...



GEMOLOGY APPLICATIONS - ACCESSORIES



# Digital Camera Adaptors



ZERTIFIKAT ♦ CERTIFICATE ♦ CERTIFICADO ♦ CERTIFICAT

# CERTIFICATE



MANAGEMENT SERVICE

The Certification Body of  
**TÜV AMERICA, INC.**  
Danvers, Massachusetts USA

hereby certifies that

**OLYMPUS (INDIA) PVT. LTD.**  
A-5, MOHAN CO-OPERATIVE INDUSTRIAL ESTATE,  
MATHURA ROAD, NEW DELHI - 110 044, INDIA

has implemented a Quality Management System  
in accordance with:

## ISO 9001:2000

The scope of this Quality Management System includes:

### Manufacturing, Sales & Service of Microscopes

*Further Clarifications regarding the scope of this certificate and the applicability of  
ISO 9001:2000 requirements may be obtained by consulting the organization.*

This certificate is valid until: **December 31, 2006**

Certificate Registration No: **951 04 2437**

Issued: **February 2, 2004**



  
\_\_\_\_\_  
Gary V. Minks  
Director, TÜV America Certification Body  
Member of TÜV SÜD Deutschland Group





India

TÜV SÜddeutschland India • C-27, 11th Floor, Qutab Institutional Area, New Delhi - 110016 • Tel: +91-11-51688423-27 • Fax: +91-11-51688421

Ref: TUV/DEL/CE/ 033

Date:15 Jun 2005

To,

Olympus India Pvt Ltd  
A/4, Mohan Co-operative Industrial Estate  
Mathura Road, New Delhi-44, India

Attn: Mr H.S.Ahuja

Dear Sir,

Please refer our quotation Ref:P/279/1026/Q/06062004 Rev 1 Date: 28 Feb 2004 for CE marking of your MLX Series of Microscopes

We do hereby confirm that testing of MLX-DX model of microscope as per the following directives and standard was successfully completed on 3<sup>rd</sup> May 2005

| Directives      | Standards  |
|-----------------|--|
| 73/23/EEC (LVD) | EN 61010-1-90+A1 :92+A2:95                         |
| 89/336/EC (EMC) | CISPER-22, IEC 61000-4-2,4,5,11<br>IEC 61000-3-2,3 |

You may now affix CE symbol as shown below to above mentioned product after giving Declaration of Conformity and preparation of TCF

CE

Please note that the thickness of the letters should be at least 5mm

Thank You

Saraj Patel  
Dy. Gen Manager  
TUV Sudeutschland India Pvt. Ltd.,



TUV SÜddeutschland India  
TUV SUD Group

321 Solitare Corporate Park  
Chakala, Andheri (East),  
Mumbai - 400053

Tel: +91-22-56923415  
Fax: +91-22-56923418  
www.tuvindia.com

PAN No: AABCT0716G  
TAN No: MUMT09385F  
Service Tax No: M-IV/S1/T1/06

## Branches & Service Centres

### **Mumbai Office:**

BSEL Tech Park, B-611, 6th Floor,  
Opp. Vashi Station, Vashi,  
Navi Mumbai-400 705  
Tel: 022-39180452, 39180447  
Fax: 022-27812906  
E-mail: mumbai@dssimage.com

### **Bangalore Office:**

245, Raheja Arcade, 2nd Floor,  
1/1 Koramangla Industrial Layout,  
Bangalore-560 095  
Tel: 080-30278066, 30278071  
Fax: 080-25534792  
E-mail: natesh@dssimage.com

### **Chennai Office:**

Flat No. 5, 11rd Floor, 19,  
Conran Smith Road,  
Gopalapuram, Chennai-600 086  
Tel: 044-28350707  
Fax: 044-28352458  
E-mail: chennai@dssimage.com

### **Secunderabad Office:**

Flat No. 202, H. No. 6-1-190,  
Gharonda Kanaklaxmi Apts.,  
Padma Rao Nagar,  
Secunderabad-500 025  
Tel: 040-27506884  
Fax: 040-27508205  
E-mail: rajareddy.oi@dssimage.com

### **Kolkata Office:**

14, Bipin Behari Ganguli Street,  
(Bow Bazar Street), Kolkata-700 012  
Tel: 033-22360032, 22362410  
Fax: 033-22367787  
E-mail: kolkata@dssimage.com

### **Pune Office:**

Room No. 11, Rajnigandha Apartments,  
Revenue Colony, Opp. Bhushan Hotel,  
J.M. Road, Shivaji Nagar, Pune-411 005  
Tel: 020-32306391  
Fax: 020-25510233  
E-mail: pune@dssimage.com

### **Chandigarh Office:**

SCO-103, 11nd Floor,  
Sector-47C, Chandigarh  
Tel: 0172-5049894  
E-mail: jasmeen@dssimage.com

# OLYMPUS

# Magnüs



ISO 9001:2000  
Certified Company

## OLYMPUS INDIA PVT. LTD.

**Corporate Office & Mailing Address : A-5, Mohan Co-operative Industrial Estate, New Delhi-110044**

Ph. : 91-11-30886766, 30886741, 30886744 Fax : 91-11-26959382

email : oisales@dssimage.com

*Import & Trade Division : P-39, (first floor), South Ext. Part-II, New Delhi-110049*

*Works : A-3, Sector-81, Noida 201305 (U.P)*

**OLYMPUS**

Your Vision, Our Future

*thank you*



DON'T REPLACE YOUR EYES.  
REPLACE YOUR MICROSCOPES.

